



SURTEX
INSTRUMENTS

Crafted to Excellence





SURTEX
INSTRUMENTS

Crafted to Excellence

INTRODUCTION

• INTRODUCTION TO THE CATALOGUE

Thank you for your interest in Surtex Oral, Maxillofacial and ENT Instruments. We remain committed to providing the highest quality products and customer service.

In this special catalogue you will find a comprehensive range of Oral, Maxillofacial and ENT Instruments, which have been developed following detailed consultation with leading Oral, Maxillofacial and ENT surgeons. These instruments in conjunction with Surtex General Instruments provides complete solution for oral, maxillofacial, otology, rhinology, rhinoplasty, tonsillectomy and tracheotomy surgeries.

We appreciate your patronage and look forward to a successful relationship with you.

• ABOUT SURTEX INSTRUMENTS

Surtex Instruments is a privately owned medical technology company dedicated to innovation and excellence in crafting of surgical instruments. The company started as a specialist in reconstructive surgery instruments and based on its success and innovative drive gradually expanded to cover a vast range of surgical instruments.

We now manufacture over 8,000 surgical instruments and offer innovative solutions in trauma, gynecology, cardiovascular, microsurgery and dental surgery. All our products are available to view online where you can download our entire range of catalogues or request a printed version.

WHY IS SURTEX FOCUSED ON INSTRUMENTS QUALITY?

Surgeons are highly talented and skilled individuals. They operate under significant time constraints and are required to make split second decisions while performing complex operations. At Surtex we understand that in such a scenario the surgeon requires tools that are not just fit for the purpose but are reliable, dependable and crafted by a company that has a track record for relentless pursuit for excellence in the making of surgical instruments.

• AFFORDABLE QUALITY

Surgical instruments represent a significant proportion of hospital's total assets. Hospitals and health centers are often required to minimize costs of surgical instruments without impacting quality or patient's health. Decision makers tend to respond via group purchasing contracts & multiple quotes to shop around for the lowest price. However this approach has significant drawbacks as numerous instrument companies have focused towards competing on price alone with a disregard towards instruments quality while claiming their instruments to be guaranteed and free of defect.

The outcome for the buyer being a perceived notion of a quality instrument, which is actually little more than a knock-off and delivers only poor value for money.

At **Surtex Instruments** our success has been based on helping customers navigate successfully through the instrument procurement process by identifying real value for money through fair comparison. We believe that quality should be affordable and you should not pay too much for over priced brands and neither should you pay too much for a low quality instrument.

• SURTEX QUALITY STANDARDS

In European Union medical devices are subject to a series of uniform laws and the EC directive places rigid demands on surgical instruments. Surtex Instruments is proud to claim that we have always surpassed the legal requirements to ensure that each instrument reaching the surgeons hand is indeed a testament to surgical craftsmanship.

We are **BS EN ISO 13485:2012/ISO 9001:2008 (SGS UK)** certified and **CE** approved to ensure compliance with the standards defined in line with medical devices directive **93/42/EEC**.



• LIFE TIME WARRANTY

We encourage our customers to calculate lifetime costs prior to making purchase decisions and investing in surgical instruments.

Surtex Instruments are guaranteed for life* against manufacturing defects of material and workmanship. This means we will repair or replace free of charge any surgical instrument that does not meet its functional requirements when used for its intended purpose and maintained properly.



*Conditions apply.

[OUR PRODUCT BRANDS]



A new dimension in cutting and dissecting performance. Our class leading integral self-sharpening technology standard on all Surtex scissors enabling precise cutting with minimum force even after pro-longed use.



These Surtex premium scissors feature finely ground micro serrations along one blade of the scissor preventing tissue slippage. The ultra sharp edge of the other blade cuts through the tissue precisely and atraumatically.



Instruments with this symbol benefit from tungsten carbide (TC) inserts expertly bonded on the cutting edge of the instrument. This material is harder than stainless steel and enables the instruments to last up to 6 times longer than standard instruments.



A true craftsmanship masterpiece and surgeons favourite. These are the best scissors in our range and combine the features of SuperEdge™ & UltraCut™ into a perfect pair of scissors that is precise, durable and reliable.



The only choice when the perfect grip matters. Utilising Creep Feed technology individual teeth are grinded into hard wear resistant tungsten carbide (TC). The result is super sharp teeth with well-defined point and inserts with impeccable gripping properties and instruments that last up to 6 times longer.



These instruments have jaws, which are coated with diamond dust resulting in improved gripping enabling superior manipulation of sutures, needles and tissues.



Our specially designed high strength fibre grip handles resulting in lightweight instruments. These ergonomic handles provide the surgeon with excellent feel, better grip, reduced shock, and less hand fatigue.



Surtex ENT diagnostic and anaesthesia instruments featuring direct bulb illumination of xenon grade bulbs.



This range features distal fibre optic to carry light from xenon bulb allowing cold illumination on operating site in our ENT diagnostics and anaesthesia instruments.



Premium grade ENT diagnostics and anaesthesia instruments combining features of distal fibre optic technology and ultra bright LED illumination ensuring the best in clear vision of operating site.

**CURRENT REFERENCE STANDARDS
FOR SURTEX PRODUCTS**

Standard / Regulation	Details of Standard / Regulation	Date of Issue
REGULATION (EU) 2017/745	REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 5 April 2017 on medical devices, amending Directive 2001/83/EC, Regulation (EC) No 178/2002 and Regulation (EC) No 1223/2009 and repealing Council Directives 90/385/EEC and 93/42/EEC	05-05-2017
ISO 9001	Quality management systems Requirements	2015-09
EN ISO 13485:2016+A11	Medical devices; Quality management systems - Requirements for regulatory purposes	2021-09
EN ISO 14971	Medical devices - Application of risk management to medical devices	2019-12
ISO 7153-1	Surgical instruments - Materials Metals	2016-10
BS EN 10088-1	Stainless steels List of stainless steels	2014-10
ISO 15223-1	Medical devices Symbols to be used with information to be supplied by the manufacturer Part 1: General requirements	2021-07
BS EN ISO 10993-1	Biological evaluation of medical devices Evaluation and testing within a risk management process	2020-12
ISO 10993-5	Biological evaluation of medical devices Part 5: Tests for in vitro cytotoxicity	2009-06
BS EN ISO 10993-10	Biological evaluation of medical devices Tests for skin sensitization	2023-02
BS EN ISO 10993-18:2020+A1	Biological evaluation of medical devices Chemical characterization of medical device materials within a risk management process	2023-07
BS EN ISO 13402	Surgical and dental hand instruments. Determination of resistance against autoclaving, corrosion and thermal exposure	2001-03
ISO 17664-1	Processing of health care products Information to be provided by the medical device manufacturer for the processing of medical devices Part 1: Critical and semi-critical medical devices	2021-07
ISO 17665-1	Sterilization of health care products Moist heat Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices	2006-08
ISO 20417	Medical devices. Information to be supplied by the manufacturer	2021-12
IEC 62366-1	Medical devices - Part 1: Application of usability engineering to medical devices	2015-02
ISO 7151	Surgical instruments Non-cutting, articulated instruments - General requirements and test methods	1988-12
BS EN ISO 9173-1	Dentistry. Extraction forceps General requirements	2016-12
ISO 7741	Instruments for surgery Scissors and shears - General requirements and test methods	1986-02

TABLE OF CONTENTS

■ Oral Maxillo-Facial Surgery	1 - 38
Otology Instruments	39 - 92
Rhinology Instruments	93 - 165
Tonsillectomy, Laryngo-Bronchoscopy	166 - 189
Tracheostomy Instruments	190 - 195



Crafted to Excellence

OM

ORAL MAXILLO-FACIAL SURGERY

Mund-, Kiefer- und Gesichtschirurgie

Cirugía oral y maxilofacial

Chirurgie buccale et maxillo-faciale

Chirurgia maxillo-facciale

Finger Protectors, Mouth Gag
 Finger-Protektoren, Mund-Gags
 Protectores de dedo, Mordaza De la Boca
 Protections de doigt, Bâillon De Bouche
 Protezioni Finger, Bocca Bavaglio



OM-011-01
Adams Napkin Holder
 Consisting of Chain and
 Two Crocodile Forceps



OM-030-13
Heister Mouth Gag
 Overall Length 13 cm (5")



Langenbeck Finger Protector

- OM-015-01** Ø 18.0 mm
OM-015-02 Ø 21.0 mm
OM-015-03 Ø 22.0 mm



OM-018-10
Mauder Mouth Gag
 Wood



OM-032-11
Obwegeser Mouth Gag
 Overall Length 11 cm (4 1/4")



Mouth Gag
Mund-Gags
Mordaza De la Boca
Bâillon De Bouche
Bocca Bavaglio

OM

ORAL MAXILLO-FACIAL SURGERY



OM-035-11
11 cm (4 1/4")

OM-035-14
13 cm (5")

Molt Mouth Gag



OM-039-12
Denhart Mouth Gag
16.0 mm
Overall Length 13 cm (5")

Mouth Gag
 Mund-Gags
 Mordaza De la Boca
 Bâillon De Bouche
 Bocca Bavaglio



Roser-Koenig Mouth Gag

OM-041-16
 16 cm (6 1/4")

OM-041-20
 20 cm (8")



Fergusson-Ackland Mouth Gag

OM-043-18
 14.5 cm (5 3/4")

Mouth Gag
 Mund-Gags
 Mordaza De la Boca
 Bâillon De Bouche
 Bocca Bavaglio

OM


Jennings Mouth Gag

OM-047-09

10.5 cm (4")

OM-047-11

12.5 cm (5")

OM-047-13

14.5 cm (5 ¾")

OM-047-15

16.5 cm (6 ½")



Whitehead Mouth Gag

Without Tongue Depressor Ratchet

OM-049-11

11 cm (4 ¼")

OM-049-13

13 cm (5")

OM-049-15

15 cm (6")



OM-051-13

Whitehead Mouth Gag

With Ratchet for Tongue Depressors

Overall Length 13 cm (5")



OM-058-00

Negus Chest Support For Mouth Gag

Overall Length 20 cm (8")

Mouth Gag
 Mund-Gags
 Mordaza De la Boca
 Bâillon De Bouche
 Bocca Bavaglio



OM-060-00
Kilner-Doughty Mouth Gag
 Complete
 With 3 Tongue Depressors



OM-061-00
Kilner-Doughty Mouth Gag
 Frame Only
 Overall Length 14.5 cm (5 ¾")



Kilner-Doughty Tongue Depressor
 Fenestrated Blade for Suction
 Overall Length 13.5 cm (5 ¼")

- OM-063-01**
 70 x 29 mm
- OM-063-02**
 80 x 34 mm
- OM-063-03**
 85 x 38 mm



Doughty Tongue Depressor
 Overall Length 13.5 cm (5 ¼")

- OM-065-01**
 28 x 66 mm
- OM-065-02**
 32 x 75 mm
- OM-065-03**
 34 x 83 mm
- OM-065-04**
 37 x 93 mm
- OM-065-05**
 38 x 100 mm

Mouth Gag
 Mund-Gags
 Mordaza De la Boca
 Bâillon De Bouche
 Bocca Bavaglio

OM

ORAL MAXILLO-FACIAL SURGERY


OM-070-00
Davis-Boyle Mouth Gag Complete

With 5 Tongue Depressors

Ref:- OM-073-01 to OM-073-05


OM-071-00


Adult

OM-071-01

Child

Davis-Boyle Mouth Gag

Frame Only


Davis-Boyle Tongue Depressor

With Fixed Suction Tube

Overall Length 13.5 cm (5 1/4")

OM-073-01

110 x 26 mm

OM-073-02

95 x 25 mm

OM-073-03

80 x 24 mm

OM-073-04

65 x 23 mm

OM-073-05
Davis-Boyle Tongue Depressor

Without Suction Tube, 30 x 20 mm

Overall Length 13 cm (5")

Mouth Gag
Mund-Gags
Mordaza De la Boca
Bâillon De Bouche
Bocca Bavaglio



OM-080-00
Complete

OM-081-00
Frame Only

McIvor Mouth Gag
With 3 Tongue Depressors
Ref: OM-083-01 to OM-083-03
Overall Length 15 cm (6")



OM-083-01
100 x 26 mm

OM-083-02
85 x 25 mm

OM-083-03
75 x 23 mm

McIvor Tongue Depressor Blade
Overall Length 13.5 cm (5 1/4")



OM-090-00
Complete

OM-091-00
Frame Only

Dingmann Mouth Gag
With 3 Tongue Depressors
Ref: OM-093-01 to OM-093-03
Overall Length 14.5 cm (5 3/4")



OM-093-01
65 x 26 mm

OM-093-02
70 x 29 mm

OM-093-03
80 x 32 mm

Dingmann Tongue Depressor Blade
Overall Length 13.5 cm (5 1/4")

Mouth Gag
 Mund-Gags
 Mordaza De la Boca
 Bâillon De Bouche
 Bocca Bavaglio

OM

ORAL MAXILLO-FACIAL SURGERY


Negus Tongue Depressor

Overall Length 13.5 cm (5 1/4")

OM-097-01


70 x 20 mm

OM-097-02

80 x 30 mm

OM-097-03

95 x 35 mm


Russel-Davis Tongue Depressor

Overall Length 13.5 cm (5 1/4")

OM-099-01

18 x 55 mm

OM-099-05

32 x 75 mm

OM-099-02

22 x 65 mm

OM-099-06

37 x 80 mm

OM-099-03

25 x 75 mm

OM-099-07

40 x 90 mm


OM-099-04

28 x 65 mm

OM-099-08

45 x 100 mm

Lip and Cheek Retractors
 Lippen-und Wangenhalter
 Labios y Mejillas Retractores
 Lèvres et jugaux
 Labbra e guance Divaricatori



Draffin Bipod Complete

Standard, With 2 Draffin Bipod Rod Ref:- OM-100-48

OM-100-00
 Complete Set

OM-100-01
 Maguaran, Plate For Holding Draffin Bipod
 40 cm (15 ¾")

OM-100-48
 Draffin Bipod Rod, With 4 Rings for Varying Height
 48 cm (19")



McKesson Mouth Prop

OM-110-01
 Small

OM-110-02
 Medium

OM-110-03
 Large

Lip and Cheek Retractors
 Lippen-und Wangenhalter
 Labios y Mejillas Retractores
 Lèvres et jugaux
 Labbra e guance Divaricatori

OM

ORAL MAXILLO-FACIAL SURGERY



OM-121-11
Martin-Simplex Lip and Cheek Retractor
 Overall Length 10 cm (4")



OM-123-12
Obwegeser Lip and Cheek Retractor
 Overall Length 12 cm (4 ¾")



OM-127-13
Bishop Lip and Cheek Retractor
 27.0 mm
 Overall Length 13.5 cm (5 ¼")



OM-131-15
Cawood-Minnesota Cheek Retractor
 20 mm
 Overall Length 15 cm (6")





OM-135-14
Minnesota Cheek Retractor
 Overall Length 14.5 cm (5 ¾")



OM-141-15
Hartmann Tongue Depressor
 72 x 23 mm
 Overall Length 15 cm (6")

Tongue Depressors
 Zungesenker
 Bajalenguas
 Abaisse-langue
 Abbassalingua



OM-143-15
Tobold Tongue Depressor
 88 x 25 mm
 Overall Length 15 cm (6")



OM-145-14
Andrew Tongue Depressor
 Overall Length 14.5 cm (5 ¾")


Tongue Depressors
 Zungesenker
 Bajalenguas
 Abaisse-langue
 Abbassalingua

OM

ORAL MAXILLO-FACIAL SURGERY



Wieder Tongue Depressor


OM-147-01
 36 mm
 14 cm (5 1/2")

OM-147-02
 43 mm
 15 cm (6")

OM-151-14
Standard Tongue Depressor
 Overall Length 14.5 cm (5 3/4")

OM-153-15
Koenig Tongue Depressor
 Overall Length 15 cm (6")

OM-155-17
Mayo Tongue Depressor
 Overall Length 17 cm (6 3/4")

OM-159-19
Bruenings Tongue Depressor
 19 mm - 15 mm
 Overall Length 19 cm (7 1/2")

Tongue Holding Forceps
 Zungehaltezangen
 Pinza tira-lengua
 Pinces tire-langue
 Pinze per lingua



OM-170-17
Collin Tongue Holding Forcep
 24 x 27 mm
 Overall Length 17 cm (6 ¾")



OM-173-17
Young Tongue Holding Forcep
 Overall Length 15 cm (6")

Oral Retractors
 Orale Wundhaken
 Separadores orales
 Rétracteurs orales
 Divaricatori orali

OM

ORAL MAXILLO-FACIAL SURGERY


Henahan Retractor

Overall Length 18 cm (7")


OM-201-01
 Fig. 1

OM-201-02
 Fig. 2

OM-201-03
 Fig. 3

OM-251-24
Sigmoid Notch Extra Oral Retractor

Concave Blade Full Bent Tip, 9.5 mm x 16 mm

Overall Length 22 cm (8 3/4")

Oral Retractors
 Orale Wundhaken
 Separadores orales
 Rétracteurs orales
 Divaricatori orali



Soft Tissue Retractors

Weichteilhaken

Separadores de tejidos blandos

Rétracteurs des tissus mous

Divaricatori tessuti molli

OM
ORAL MAXILLO-FACIAL SURGERY


Obwegeser Soft Tissue Retractor
Bent Upwards
Overall Length 24 cm (9 1/2")



OM-271-16	OM-271-55
5 x 16 mm	12 x 55 mm
OM-271-25	OM-271-70
7 x 25 mm	14 x 70 mm
OM-271-35	OM-271-80
10 x 35 mm	16 x 80 mm
OM-271-42	
10 x 42 mm	




Obwegeser Soft Tissue Retractor
Bent Downwards
Overall Length 24 cm (9 1/2")




OM-273-16	OM-273-55
5 x 16 mm	12 x 55 mm
OM-273-25	OM-273-70
7 x 25 mm	14 x 70 mm
OM-273-35	OM-273-80
10 x 35 mm	16 x 80 mm
OM-273-42	
10 x 42 mm	

Soft Tissue Retractors
 Weichteilhaken
 Separadores de tejidos blandos
 Rétecteurs des tissus mous
 Divaricatori tessuti molli



Obwegeser Soft Tissue Retractor
 Bent Upwards
 Overall Length 21.5 cm (8 1/2")



Obwegeser Soft Tissue Retractor
 Overall Length 21 cm (8 1/4")

OM-279-30
 12 x 30 mm

OM-279-45
 12 x 45 mm

OM-279-60
 12 x 60 mm

Oral Retractors
Orale Wundhaken
Separadores orales
Rétlecteurs orales
Divaricatori orali

OM


ORAL MAXILLO-FACIAL SURGERY



Obwegeser Ramus Retractor
Overall Length 24.5 cm (9 3/4")



OM-281-22
15 x 70 mm




OM-283-22
22 x 70 mm


Oral Retractors
 Orale Wundhaken
 Separadores orales
 Rétracteurs orales
 Divaricatori orali



Obwegeser Ramus Retractor (Nasal Spine Retractor)
 Overall Length 22 cm (8 3/4")



OM-287-22
 10 x 45 mm



OM-289-22
 10 x 72 mm


Oral Retractors
 Orale Wundhaken
 Separadores orales
 Rétracteurs orales
 Divaricatori orali

OM


ORAL MAXILLO-FACIAL SURGERY



Steinhauser Zygomatic Retractor
 Overall Length 24 cm (9 1/2")



OM-291-01
 Right



OM-291-02
 Left



OM-301-16
Obwegeser Chin Retractor
 Overall Length 16.5 cm (6 1/2")

Oral Retractors
 Orale Wundhaken
 Separadores orales
 Rétracteurs orales
 Divaricatori orali



Obwegeser Channel Retractor
 Overall Length 17 cm (6 $\frac{3}{4}$ ')

OM-305-08
 8.0 mm


OM-305-10
 10.0 mm



Hargis Medial Ramus Retractor
 Overall Length 19.5 cm (7 $\frac{3}{4}$ ')

OM-315-01
 Left

OM-315-02
 Right



Oral Retractors, Periosteal Elevator
 Orale Wundhaken, Raspatorium
 Separadores orales, Elevador perióstico
 Rétracteurs orales, Ascenseur périostée
 Divaricatori orali, Scollaperiostio

OM
ORAL MAXILLO-FACIAL SURGERY

OM-319-18

Hargis Mandibular Body Retractor
 Overall Length 18 cm (7")


OM-331-19

Hargis Anterior Border Stripper
 Overall Length 19 cm (7 1/2")


OM-333-06

Right
 6.5 mm

OM-335-06

Left
 6.5 mm

OM-333-09

Right
 9.0 mm

OM-335-09

Left
 9.0 mm



Periosteal Elevator
Raspatorium
Elevador perióstico
Ascenseur périostée
Scollaperiostio



OM-337-18
Hargis Perform Rim Retractor
Overall Length 18.5 cm (7 $\frac{1}{4}$ ")



OM-339-22
Hargis Pterygoid Osteotome
Overall Length 22 cm (8 $\frac{3}{4}$ ")



OM-341-19
Hargis Bone Holding Forcep
With Chain
Overall Length 19 cm (7 $\frac{3}{4}$ ")

Periosteal Elevator
 Raspatorium
 Elevador perióstico
 Ascenseur périostée
 Scollaperiostio

OM

ORAL MAXILLO-FACIAL SURGERY



Obwegeser Periosteal Elevator
 Overall Length 17.5 cm (7")



OM-361-04
 4.0 mm



OM-361-06
 6.0 mm



OM-361-07
 7.0 mm




Obwegeser Elevator
 Overall Length 20.5 cm (8")



OM-381-21
 Fig. 1




OM-383-21
 Fig. 2



OM-361-09
 9.0 mm

OM-361-11
 11.0 mm





Ramus Pusher, Steinhauser Reposition Instrument
Knopfsonden, Hohlsonden
Sondas abotonadas, Sondas acanaladas
Sondes boutonnées, Sondes cannelées
Sonde a bottone, Sonde cannellate

OM



OM-415-24
Terry Ramus Pusher
Overall Length 24 cm (9 1/2")



OM-421-18
Steinhaus Reposition Instrument
Overall Length 18 cm (7")



Maxillary Awls
Unterkieferahlen
Leznas Maxilares
Alènes maxillaires
Punteruoli mascellari

OM
ORAL MAXILLO-FACIAL SURGERY


OM-501-15
15 cm (6")



OM-501-17
17.5 cm (7")



OM-503-24
Obwegeser Zygomatic Arch Awl
Overall Length 24 cm (9 ½")



OM-507-18
Obwegeser Maxillary Awl
Overall Length 18 cm (7")

Obwegeser Maxillary Awl

Wire Loop Forcep, Ramus Separator
 Knopfsonden, Hohlsonden
 Sondas abotonadas, Sondas acanaladas
 Sondes boutonnées, Sondes cannelées
 Sonde a bottone, Sonde cannellate



OM-540-15
Obwegeser Wire Loop Forcep
 Overall Length 15 cm (6")



OM-550-20
Smith Ramus Separator
 Overall Length 23.5 cm (9 1/4")

Disimpaction Forcep
 Repositionszangen
 Pinzas para reposición
 Pinces de reposicion
 pinze per riduzione

OM

ORAL MAXILLO-FACIAL SURGERY



OM-554-01
 Right


OM-554-02
 Left



OM-570-01
 Left



OM-570-02
 Right



Smith Sagittal Split Separator
 Overall Length 18.5 cm (7 1/4")

Rowe Disimpaction Forcep
 Overall Length 22.5 cm (8 3/4")

Zygoma Elevators
 Zygoma-Elevatorien
 Elevadores, cigomáticos
 Élevateurs zygomatiques
 Elevatori per zigomo



OM-601-21
Gillies Zygoma Elevator
 Overall Length 20 cm (8")

OM-602-21
Dingmann Zygoma Elevator
 Ø 5.0 mm
 Overall Length 21 cm (8 1/4")

OM-605-21
Rowe Zygoma Elevator
 Overall Length 21.5 cm (8 1/2")

OM-607-19
Boise Nasal Fracture Elevator
 Overall Length 19.5 cm (7 3/4")


Osteotomes
Osteotome
Osteótomas
Ostéotomes
Osteotomi

OM
ORAL MAXILLO-FACIAL SURGERY


Dingmann Bone Hook
Overall Length 20 cm (8")

OM-613-10
10.0 mm

OM-615-20
20.0 mm



OM-620-06
Straight
6.0 mm

OM-620-08
Straight
8.0 mm

Bauer Osteotome
Overall Length 20.5 cm (8")



OM-621-06
Curved
6.0 mm

OM-621-08
Curved
8.0 mm

Bauer Osteotome
Overall Length 20.5 cm (8")



OM-623-08
Epker-Bauer Osteotome
Curved, 8.0 mm
Overall Length 17.5 cm (7")



Dautrey-Munro Osteotome
Overall Length 16 cm (6 1/4")

OM-626-04 4.0 mm	OM-626-06 6.0 mm	OM-626-08 8.0 mm
OM-626-10 10.0 mm	OM-626-12 12.0 mm	OM-626-14 14.0 mm
OM-626-16 16.0 mm		

Osteotomes
Osteotome
Osteótomas
Ostéotomes
Osteotomi

OM
ORAL MAXILLO-FACIAL SURGERY


Dunn-Dautrey Osteotome

OM-628-04
Straight
4.0 mm
16.5 cm (6 1/2")

OM-628-06
Straight
6.0 mm
16 cm (6 1/4")

OM-628-08
Straight
8.0 mm
16 cm (6 1/4")

OM-629-06
Curved
6.0 mm
16 cm (6 1/4")

OM-631-06
Strongly Curved
6.0 mm
16 cm (6 1/4")



Dunn-Dautrey Elevator
Overall Length 15.5 cm (6")

OM-633-48
6 x 48 mm

OM-633-56
6 x 56 mm



OM-637-12
Dolwick-Type TMJ Disc Clamp
Overall Length 12 cm (4 $\frac{3}{4}$ ")

OM-639-14
Dolwick-Type TMJ Condylar Retractor
Overall Length 14.5 cm (5 $\frac{1}{2}$ ")

OM-641-17
Dolwick-Type TMJ Elevator
Overall Length 16 cm (6 $\frac{3}{4}$ ")


Osteotomes
Osteotome
Osteótomas
Ostéotomes
Osteotomi

OM
ORAL MAXILLO-FACIAL SURGERY


OM-700-17
Turvey Palatal Expander
Overall Length 17 cm (6 ¾")



OM-730-19
Terry Septal Osteotome
Overall Length 19 cm (7 ¾")



Obwegeser Septonasal Osteotome
Graduated
Overall Length 18.5 cm (7 ¼")

OM-735-04
4.0 mm



OM-735-06
6.0 mm



OM-735-08
8.0 mm





Osteotomes
Osteotome
Osteótomas
Ostéotomes
Osteotomi

OM

ORAL MAXILLO-FACIAL SURGERY



Obwegeser Osteotome
For Alveolar Process
Overall Length 15 cm (6")



SURTEX
INSTRUMENTS

Crafted to Excellence

OT

OTOLOGY INSTRUMENTS

Otologie

Otología

Otologie

Otologia

Ear Specula

Ohrspkula

Espéculos auriculares

Spéculums auriculaires

Specoli auricolari


OT-010-01
Hartmann Ear Speculum

Set of 3 Ref: OT-012-00 to OT-012-02


OT-010-02
Hartmann Ear Speculum

Set of 3 Ref: OT-012-02 to OT-012-04


OT-010-03
Hartmann Ear Speculum

Set of 4 Ref: OT-012-01 to OT-012-04


OT-010-04
Hartmann Ear Speculum

Set of 5 Ref: OT-012-00 to OT-012-04


OT-012-00

 Fig. 0
 \varnothing 3.0 mm

OT-012-01

 Fig. 1
 \varnothing 4.0 mm

OT-012-02

 Fig. 2
 \varnothing 5.0 mm

OT-012-03

 Fig. 3
 \varnothing 6.0 mm

OT-012-04

 Fig. 4
 \varnothing 7.0 mm

Hartmann Ear Speculum

Overall Length 3.6 cm (1 ½")

Ear Specula
Ohrspekula
Espéculos auriculares
Spéculums auriculaires
Specoli auricolari

OT
OTOLOGY INSTRUMENTS


OT-014-00
Fig. 0
2.5 x 3.0 mm

OT-014-01
Fig. 1
4.5 x 5.5 mm

OT-014-02
Fig. 2
5.5 x 6.0 mm

OT-014-03
Fig. 3
6.5 x 7.0 mm

OT-014-04
Fig. 4
7.5 x 8.0 mm

Hartmann Ear Speculum (Oval)
Overall Length 3.6 cm (1 ½")



OT-020-01
Toynbee Ear Speculum
Set of 3 Ref: OT-022-02 to OT-022-04



OT-020-02
Toynbee Ear Speculum
Set of 4 Ref: OT-022-01 to OT-022-04



OT-022-01
Fig. 1
Ø 4.0 mm

OT-022-02
Fig. 2
Ø 5.0 mm

OT-022-03
Fig. 3
Ø 6.0 mm

OT-022-04
Fig. 4
Ø 7.0 mm

Toynbee Ear Speculum
Overall Length 3.6 cm (1 ½")



OT-030-01
Hartmann Ear Speculum (Black)
Set of 3 Ref: OT-032-00 to OT-032-02



OT-030-02
Hartmann Ear Speculum (Black)
Set of 3 Ref: OT-032-02 to OT-032-04



Ear Specula
Ohrspekula
Espéculos auriculares
Spécums auriculaires
Specoli auricolari



OT-030-03

Hartmann Ear Speculum (Black)

Set of 4 Ref: OT-032-01 to OT-032-04



OT-032-00

Fig. 0
Ø 3.0 mm

OT-032-01

Fig. 1
Ø 4.0 mm

OT-032-02

Fig. 2
Ø 5.5 mm

OT-032-03

Fig. 3
Ø 6.0 mm

OT-032-04

Fig. 4
Ø 7.0 mm

Hartmann Ear Speculum (Black)
Overall Length 3.6 cm (1 1/2")



OT-040-01

Zoellner Ear Speculum (Black)

Set of 6 Ref: OT-042-01 to OT-042-06



OT-042-01

Fig. 1
5.5 x 6.5 mm

OT-042-02

Fig. 2
6.5 x 7.5 mm

OT-042-03

Fig. 3
7.5 x 8.5 mm

OT-042-04

Fig. 4
8.5 x 9.5 mm

OT-042-05

Fig. 5
9.5 x 10.5 mm

OT-042-06

Fig. 6
10.5 x 11.5 mm

Zoellner Ear Speculum (Black)
Overall Length 3.8 cm (1 1/2")

Ear Specula
Ohrspekula
Ohrspeküle
Espéculos auriculares
Spéculums auriculaires
Specoli auricolari

OT
OTOLOGY INSTRUMENTS

OT-050-01
Farror Ear Speculum (Black)

Set of 7 Ref: OT-052-01 to OT-052-07


Farror Ear Speculum (Black)
Overall Length 3.6 cm (1 1/2")

OT-052-01

Fig. 1
Ø 5.5 mm

OT-052-02

Fig. 2
Ø 6.5 mm

OT-052-03

Fig. 3
Ø 7.5 mm

OT-052-04

Fig. 4
Ø 8.5 mm

OT-052-05

Fig. 5
Ø 9.5 mm

OT-052-06

Fig. 6
Ø 10.5 mm

OT-052-07

Fig. 7
Ø 11.5 mm

Farror Ear Speculum
Oval/Oblique End

OT-054-00

Fig. 00
3 x 4 mm

OT-054-01

Fig. 0
4 x 5 mm

OT-054-02

Fig. 1
5 x 6 mm

OT-054-03

Fig. 2
6 x 7 mm

OT-054-04

Fig. 3
7 x 8 mm

OT-054-05

Fig. 4
8 x 10 mm

OT-054-06

Fig. 5
9 x 11 mm

OT-054-07

Fig. 6
10 x 12 mm

OT-054-08

Fig. 7
11 x 13 mm

Ear Specula
Ohrspekula
Espéculos auriculares
Spécums auriculaires
Specoli auricolari



OT-060-12
Rack for Storage of Ear Speculum
For 12 Specula



- OT-064-02**
Ø 3.0 mm
OT-064-04
Ø 4.0 mm
OT-064-06
Ø 5.0 mm

Richnau-Holmgren Ear Speculum



OT-070-00
Siegle Otoscope Complete
With 3 Specula Ref: OT-072-01 to OT-072-03



Siegle Ear Speculum

OT-072-01 Fig. 1 Ø 6.5 mm	OT-072-02 Fig. 2 Ø 7.0 mm	OT-072-03 Fig. 3 Ø 8.0 mm
--	--	--

Tuning Forks
Stimmgabeln
Diapasones
Diapasons
Diapason

OT

OTOTOLOGY INSTRUMENTS


OT-100-10
Hartmann Tuning Fork

Set of 5

 Ref: OT-102-03 to OT-102-05 and
OT-104-04 to OT-104-05

OT-100-20
Hartmann Tuning Fork

Set of 6

 Ref: OT-102-03 to OT-102-06 and
OT-104-04 to OT-104-05

OT-100-30
Hartmann Tuning Fork

Set of 6

 Ref: OT-102-03 to OT-102-06 and
OT-104-04 to OT-104-05

OT-102-01

 Frequency
C 128

OT-102-02

 Frequency
C¹ 256

OT-102-03

 Frequency
C² 512

OT-102-04

 Frequency
C³ 1024

OT-102-05

 Frequency
C⁴ 2048

OT-102-06

 Frequency
C⁵ 4096

OT-102-60
Hartmann Tuning Fork

Set of 6

Ref: OT-102-01 to OT-102-06

Hartmann Tuning Fork



Tuning Forks
Stimmgabeln
Diapasones
Dipasóns
Diapason



Hartmann Tuning Fork (French)



OT-110-00
Tuning Fork
Set of 5
Ref: OT-112-02 to OT-112-06

Tuning Forks
 Stimmgabeln
 Diapasones
 Diapasons
 Diapason

OT

OTOTOLOGY INSTRUMENTS


OT-112-02

 Frequency
 C 128

OT-112-03

 Frequency
 C¹ 256

OT-112-04

 Frequency
 C² 512

OT-112-05

 Frequency
 C³ 1024

OT-112-06

 Frequency
 C⁴ 2048

Tuning Fork

OT-122-01

 Frequency
 C 128

OT-122-02

 Frequency
 C¹ 256

OT-122-03

 Frequency
 C² 512

OT-122-04

 Frequency
 C³ 1024


OT-122-05

 Frequency
 C⁴ 2048

OT-122-06

 Frequency
 C⁵ 4096

Lucae Tuning Fork



Tuning Forks
Stimmgabeln
Diapasones
Diapasons
Diapason

OT
OTOLOGY INSTRUMENTS

OT-136-00

Neuro Tuning Fork
With Hammer, Pinwheel, Brush
and Needle


OT-138-00

Koenig Tuning Fork
Moveable Up to g
Frequency A 109


OT-140-00

Koenig Tuning Fork
Frequency a¹ 440

Ear Catheters, Polypus Snare
 Ohrkatheter, Polypenschnüre
 Sondas auriculares, Polipótomas
 Sondes auriculaires, Serre-noeuds
 Sonde auricolari, Serranodi



OT-150-01

Fig. 1

OT-150-02

Fig. 2

OT-150-03

Fig. 3

Hartmann Ear Catheter
 Overall Length 15 cm (6")



OT-161-18
Krause Polypus Snare
 Open Tip

Overall Length 18.5 cm (7 1/4")



OT-163-18
Krause-Voss Polypus Snare
 Close Tip
 Overall Length 18.5 cm (7 1/4")

Polypus Snare
Polypenschnüre
Polipótomas
Serre-noeuds
Serranodi

OT

OTOLOGY INSTRUMENTS

**OT-165-16**

Lange-Wilde Polypus Snare
Overall Length 16 cm (6 1/4")

OT-170-30

Ø 0.3 mm

OT-170-40

Ø 0.4 mm

Snare Wire

Overall Length 100 cm (39 1/2")



OT-201-11
Troeltsch (Wilde) Ear Forcep
Overall Length 11 cm (4 $\frac{3}{4}$ ")

OT-201-13
Troeltsch (Wilde) Ear Forcep
Overall Length 12.5 cm (5")

OT-203-11
Troeltsch (Wilde) Ear Forcep
Overall Length 11 cm (4 $\frac{3}{4}$ ")

Ear Forceps
 Ohrrinzelten
 Pinzas auriculares
 Pince auriculaires
 Pinze auricolari

OT

OTOTOLOGY INSTRUMENTS



OT-205-11
Troeltsch (Wilde) Ear Forcep
 1 x 2 Teeth
 Overall Length 11 cm (4 3/4")



OT-207-18
Adson Ear Forcep
 Bayonet Shaped
 Overall Length 17.5 cm (7")



OT-211-14
Gruenwald (Lucae) Ear Forcep
 Bayonet Shaped
 Overall Length 14.5 cm (5 3/4")



OT-213-16
Jansen Ear Forcep
Bayonet Shaped
Overall Length 16 cm (6 $\frac{1}{4}$ ")

OT-213-21
Jansen Ear Forcep
Bayonet Shaped
Overall Length 21.5 cm (8 $\frac{1}{2}$ ")

OT-215-14
Gruenwald (Lucae) Ear Forcep
1 x 2 Teeth, Bayonet Shaped
Overall Length 14.5 cm (5 $\frac{3}{4}$)

Ear Polypus Forceps
 Ohrrölpenzangen
 Pinzas para pólipos auriculares
 Pince pour polypes auriculaires
 Pinze per polipi auricolari

OT

OTOLGY INSTRUMENTS




OT-221-10
Littauer Ear Polypus Forcep
 Overall Length 10 cm (4")



OT-223-12
Hartmann Ear Polypus Forcep
 Fine Model
 Overall Length 12 cm (4 3/4")



OT-225-14
Hartmann Ear Polypus Forcep
 Regular Model
 Overall Length 14 cm (5 1/2")



Tilley Ear Polypus Forcep


OT-229-14
 14 cm (5 1/2")

OT-229-16
 16 cm (6 1/4")

Foreign Body Forcep, Alligator Forcep
 Knopfsonden, Hohlsonden
 Sondas abotonadas, Sondas acanaladas
 Sondes boutonnées, Sondes cannelées
 Sonde a bouton, Sonde cannellate



OT-235-11
Buck Foreign Body Forcep
 Overall Length 11.5 cm (4 1/2")



OT-237-10
Quire Foreign Body Forcep
 Overall Length 10.5 cm (4")



OT-251-08
Hartmann-Wullstein Alligator Forcep
 Overall Length 8.5 cm (3 1/4")



OT-253-08
Hartmann-Wullstein Alligator Forcep
 Cup Size Ø 1.0 mm
 Overall Length 8.5 cm (3 1/4")

Alligator Forcep
 Hohlsonden
 Sondas acanaladas
 Sondes cannelées
 Sonde cannellate

OT

OTOTOLOGY INSTRUMENTS




Hartmann Alligator Forcep
 Serrated Jaws

OT-255-08
 8.5 cm (3 1/4")

OT-255-14
 14 cm (5 1/2")

OT-255-16
 16 cm (6 1/4")




Hartmann Alligator Forcep
 1 x 2 Teeth, Serrated Jaws

OT-257-08
 8.5 cm (3 1/4")

OT-257-14
 14 cm (5 1/2")

Ear Polypus Forceps
 Ohrpolypenzangen
 Pinzas para pólipos auriculares
 Pinces pour polypes auriculaires
 Pinze per polipi auricolari



Ear Polypus Forceps
 Ohrpolypenzangen
 Pinzas para pólipos auriculares
 Pince pour polypes auriculaires
 Pinze per polipi auricolari

OT


OTOLGY INSTRUMENTS



OT-265-08
Struempel Alligator Forcep
 Overall Length 8.5 cm (3 1/4")



OT-267-08
Hartmann-Weingartner Alligator Forcep
 Overall Length 8.5 cm (3 1/4")



OT-269-08
Duckbill Alligator Forcep
 Overall Length 8.5 cm (3 1/4")



OT-271-08
Hartmann-Hoffmar Alligator Punch Forcep
 Hartmann-Hoffmar Alligator Punch Forcep
 Overall Length 8.5 cm (3 1/4")

Ear Polypus Forceps
 Ohrpolypenzangen
 Pinzas para pólipos auriculares
 Pinces pour polypes auriculaires
 Pinze per polipi auricolari



OT-273-08
Hartmann Alligator Forcep
 2 x 7 mm
 Oval Cup Shaped Jaws
 Overall Length 8.5 cm (3 1/4")




McGee Micro Alligator Forcep
 Overall Length 8 cm (3")




Mini-McGee Micro Alligator Forcep
 Overall Length 8 cm (3")



OT-281-08
 Smooth, Straight
 4.0 x 0.8 mm




OT-283-08
 Serrated, Straight
 4.0 x 0.8 mm



OT-285-08
 Serrated
 Bent Upwards
 4.0 x 0.8 mm



OT-287-08
 Serrated, Right
 4.0 x 0.8 mm



OT-289-08
 Serrated, Left
 4.0 x 0.8 mm



OT-291-08
 Smooth, Straight
 3.5 x 0.6 mm



OT-293-08
 Serrated, Straight
 3.5 x 0.6 mm

Ear Polypus Forceps
 Ohrrölpenzangen
 Pinzas para pólipos auriculares
 Pince pour polypes auriculaires
 Pinze per polipi auricolari

OT

OTOLGY INSTRUMENTS



Micro Alligator Forcep
 Overall Length 8 cm (3")



Micro Alligator Forcep
 Cup Shaped
 Overall Length 8 cm (3")


OT-301-08

 Smooth, Straight
 4.0 x 0.6 mm

OT-303-08

 Serrated, Straight
 4.0 x 0.6 mm

OT-305-08

 Serrated, Bent
 Upwards
 4.0 x 0.6 mm

OT-307-08

 Serrated, Right
 4.0 x 0.6 mm

OT-309-08

 Serrated, Left
 4.0 x 0.6 mm

OT-311-08

 Straight
 Cup Size
 1.0 x 0.9 mm
 Jaw Size
 4.0 mm

OT-313-08

 Bent Upwards
 Cup Size
 1.0 x 0.9 mm
 Jaw Size
 4.0 mm

OT-315-08

 Right
 Cup Size
 1.0 x 0.9 mm
 Jaw Size
 4.0 mm

OT-317-08

 Left
 Cup Size
 1.0 x 0.9 mm
 Jaw Size
 4.0 mm

OT-321-08

 Bent Upwards
 Cup Size
 0.6 x 0.5 mm
 Jaw Size
 3.5 mm

OT-323-08

 Bent Upwards
 Cup Size
 0.6 x 0.5 mm
 Jaw Size
 3.5 mm

OT-325-08

 Right
 Cup Size
 0.6 x 0.5 mm
 Jaw Size
 3.5 mm

OT-327-08

 Left
 Cup Size
 0.6 x 0.5 mm
 Jaw Size
 3.5 mm

Micro Scissors
Mikroscheren
Micro-Tijeras
Micro-Ciseaux
Micro-Forbici



Bellucci Micro Scissor
Overall Length 8 cm (3")



OT-331-08
Straight
Very Fine
4.0 x 0.8 mm

OT-335-08
Right
Very Fine
4.0 x 0.8 mm

OT-333-08
Bent Upwards
Very Fine
4.0 x 0.8 mm

OT-337-08
Left
Very Fine
4.0 x 0.8 mm

Wullstein Micro Scissor
Overall Length 8 cm (3")



OT-341-08
Straight
4.0 x 1.8 mm

OT-345-08
Right
4.0 x 1.8 mm

OT-343-08
Bent Upwards
4.0 x 1.8 mm

OT-347-08
Left
4.0 x 1.8 mm


Micro Scissors, Micro Ear Forceps
 Mikroscheren, Mikro-Ohrzängchen
 Micro-Tijeras, Micro pinzas auriculares
 Micro-Ciseaux, Micro Pince auriculaires
 Micro-Forbici, Micro Pinze auricolari

OT

OTOLGY INSTRUMENTS




Shea Micro Scissor
 Overall Length 8 cm (3'')



OT-351-08

Straight
 7.0 x 1.8 mm



OT-355-08

Right
 7.0 x 1.8 mm



OT-357-08

Left
 7.0 x 1.8 mm



McGee Wire Bending Forcep
 Bent Downwards
 Overall Length 8 cm (3'')



OT-361-01

3.5 x 0.8 mm



OT-361-02

6.2 x 0.8 mm



**House-Dieter Malleus
 Nipper**
 Overall Length 8 cm (3'')



OT-365-01

Up Cutting
 1.3 mm



OT-365-02

Down Cutting
 1.3 mm



OT-365-03

Right Cutting
 1.3 mm



OT-365-04

Left Cutting
 1.3 mm



OT-367-08
Fuller Introducing Forcep
 For Teflon Prosthesis
 Overall Length 8 cm (3'')

Micro Ear Knives

Mikro-Ohrmesser

Micro cuchillos auriculares

Micro Couteaux auriculaires

Micro Coltelli


OT-401-01

Schuknecht Rotary Knife
Overall Length 15 cm (6")

OT-401-02

OT-403-15
Schuknecht Sickle Knife
Overall Length 14.5 cm (5 ¾")

OT-405-15
Wullstein Sickle Knife
Overall Length 15.5 cm (6")

OT-407-15
Tabb Mucosa Knife
Overall Length 15.5 cm (6")


Micro Ear Knives
 Mikro-Ohrmesser
 Micro cuchillos auriculares
 Micro Couteaux auriculaires
 Micro Coltelli

OT

OTOTOLOGY INSTRUMENTS



OT-409-15
Plester Flap Knife
 Overall Length 15 cm (6")



OT-413-16
Fisch Tenotome (Sickle Shaped)
 Extra Delicate
 Overall Length 16 cm (6 1/4")



OT-415-15
Plester Sickle Knife
 Overall Length 15 cm (6")



OT-417-15
Sexton Myringotomy Knife
 Overall Length 15 cm (6")



Micro Ear Knives

Mikro-Ohrmesser

Micro cuchillos auriculares

Micro Couteaux auriculaires

Micro Coltelli


OT-419-16

Sexton Tympanum Knife

Overall Length 16 cm (6 1/4")

OT-420-03

Rosen Circular Cut Knife

Fig. 3, Ø 2.6 mm

Overall Length 15 cm (6")


Rosen Circular Cut Knife

Overall Length 15 cm (6")

OT-421-01

Ø 3.0 mm

Fig. 1

OT-421-02

Ø 2.8 mm

Fig. 2

OT-421-03

Ø 2.6 mm

Fig. 3

OT-421-04

Ø 2.0 mm

Fig. 4

OT-421-05

Ø 1.6 mm

Fig. 5

OT-421-06

Ø 1.4 mm

Fig. 6

OT-421-07

Ø 1.0 mm

Fig. 7

Micro Ear Hooks
 Mikro-Ohrhähkchen
 Micro ganchitos auriculares
 Micro Crochets auriculaires
 Micro Uncini auricolari

OT

OTOTOLOGY INSTRUMENTS


OT-422-03

Ø 0.3 mm

OT-422-04

Ø 0.4 mm

OT-422-05

Ø 0.5 mm

OT-422-06

Ø 0.6 mm

OT-422-07

Ø 0.7 mm

OT-422-08

Ø 0.8 mm


OT-440-16
Plester Wire Guide and Applicator

For Prothesis Flat End With V-Notch

Overall Length 15.5 cm (6")

Fisch Footplate Perforator

Trocar Shaped


Overall Length 16 cm (6 1/4")



Micro Ear Hooks, Micro Ear Needles
 Mikro-Ohrhaken, Mikro-Ohrnadeln
 Micro ganchitos auriculares, Micro agujas auriculares
 Micro Crochets auriculaires, Micro Aiguilles auriculaires
 Micro Uncini auricolari, Micro Aghi auricolari

OT

OTOLGY INSTRUMENTS



Barbara Micro Ear Needle
 Angled 90°
 Overall Length 16 cm (6 1/4")

Barbara Micro Ear Needle
 Curved Down
 Overall Length 16 cm (6 1/4")

Micro Ear Needle
 Angled 90°
 Overall Length 15.5 cm (6")



Micro Ear Needles

Mikro-Ohrnadeln

Micro agujas auriculares

Micro Aiguilles auriculaires

Micro Aghi auricolari

Micro Ear Needle
Angled 90°
Overall Length 15.5 cm (6")



OT-469-03
0.3 mm
OT-469-06
0.6 mm
OT-469-10
1.0 mm

McGee Micro Ear Needle
For Footplate
Overall Length 15.5 cm (6")



OT-473-05
0.5 mm
OT-473-10
1.0 mm

McGee Micro Ear Needle
For Footplate
Overall Length 15.5 cm (6")



OT-475-05
0.5 mm
OT-475-10
1.0 mm

Micro Ear Needles
 Mikro-Ohrnadeln
 Micro agujas auriculares
 Micro Aiguilles auriculaires
 Micro Aghi auricolari

OT

OTOTOLOGY INSTRUMENTS



Barbara Micro Ear Needle
 Overall Length 15.5 cm (6")

McGee Micro Ear Needle
 Angled 25°, Upwards
 Overall Length 15.5 cm (6")

McGee Micro Ear Needle
 Angled 25°, Downwards
 Overall Length 15.5 cm (6")



Micro Ear Needles
Mikro-Ohrnadeln
Micro agujas auriculares
Micro Aiguilles auriculaires
Micro Aghi auricolari



Micro Ear Needle
Angled 45°, Upwards
Overall Length 15.5 cm (6")

OT-485-03
0.3 mm
OT-485-06
0.6 mm
OT-485-10
1.0 mm



Micro Ear Needle
Angled 45°, Downwards
Overall Length 15.5 cm (6")

OT-487-03
0.3 mm
OT-487-06
0.6 mm
OT-487-10
1.0 mm



Micro Ear Needle
Angled 25°
Overall Length 15.5 cm (6")

OT-491-03
0.3 mm
OT-491-06
0.6 mm

Micro Ear Needles, Micro Ear Hooks
 Mikro-Ohrnadeln, Mikro-Ohrhähchen
 Micro agujas auriculares, Micro ganchitos auriculares
 Micro Aiguilles auriculaires, Micro Crochets auriculaires
 Micro Aghi auricolari, Micro Uncini auricolari

OT

OTOTOLOGY INSTRUMENTS


OT-491-10

Shambaugh Micro Ear Needle
 Angled 25°, 1.0 mm
 Overall Length 15.5 cm (6")


OT-493-03

0.3 mm

OT-493-06

0.6 mm


OT-493-10

Shambaugh Micro Ear Needle
 Angled 45°, 1.0 mm
 Overall Length 15.5 cm (6")


OT-497-20

Shea Exploring Hook
 Angled 45°, Blunt, 2.0 mm
 Overall Length 15.5 cm (6")



Micro Curettes, Micro Dissectors, Micro Measuring Rods
 Mikro-Küretten, Mikro-Dissektoren, Mikro-Messsonden
 Micro cucharillas, Micro disecadores, Micro Sondas para medir
 Micro Curettes, Micro Dissecteurs, Micro Sondes pour mesurer
 Micro cucchiali, Micro Dissettori, Micro Sonde per misurare

OT

OTOTOLOGY INSTRUMENTS



OT-509-04
Rosen Micro Curette
 Overall Length 15.5 cm (6")

Rosen Micro Dissector
 Overall Length 15.5 cm (6")

House Micro Measuring Rod
 Overall Length 15.5 cm (6")



Micro Curettes
 Mikro-Küretten
 Micro cucharillas
 Micro Curettes
 Micro cucchiai

OT


OTOLGY INSTRUMENTS



Dingmann Otobrader


 With Two Rows of 10 Fine Teeth
 Overall Length 15 cm (6")


House Micro Ear Curette

Overall Length 14.5 cm (5 ¾")


OT-520-01

 Straight
 1.0 mm, 1.2 mm

OT-520-02

 Straight
 1.5 mm, 1.8 mm

OT-520-03

 Angled
 1.0 mm, 1.2 mm

OT-520-04


 Angled
 1.5 mm, 1.8 mm




Micro Curettes
Mikro-Küretten
Micro cucharillas
Micro Curettes
Micro cucchiai



House Micro Ear Curette
Overall Length 14.5 cm (5 $\frac{3}{4}$ ")



- OT-521-01**
2.8 mm x 3.2 mm & 2.6 x 3.5 mm
- OT-521-02**
2.2 mm x 3.8 mm & 2.8 x 4.8 mm



- OT-521-03**
0.6 mm x 0.8 mm & 0.8 x 1.0 mm



OT-525-14
Shea Micro Ear Curette
Overall Length 14.5 cm (5 $\frac{3}{4}$ ")

Micro Ear Needles
 Mikro-Ohrnadeln
 Micro agujas auriculares
 Micro Aiguilles auriculaires
 Micro Aghi auricolari


OT

OTOTOLOGY INSTRUMENTS




Zoellner Micro Ear Needle
 Overall Length 13.5 cm (5 1/2")


OT-535-01
 Straight



OT-535-02
 Angled 45°
 Upwards
 0.5 mm



OT-535-03
 Angled 45°
 Upwards
 1.0 mm



OT-535-04
 Angled 45°
 Downwards
 0.5 mm




OT-535-05
 Angled 45°
 Downwards
 1.0 mm



OT-535-06
 Angled 90°
 Upwards
 0.5 mm



OT-535-07
 Angled 90°
 Downwards
 0.5 mm



OT-535-08
 For
 Footplate
 1.0 mm





Micro Ear Needles
 Mikro-Ohrnadeln
 Micro agujas auriculares
 Micro Aiguilles auriculaires
 Micro Aghi auricolari

OT

OTOTOLOGY INSTRUMENTS



Zoellner Tympanum Needle
 Overall Length 13.5 cm (5 1/2")

Zoellner Respiratory
 Overall Length 13.5 cm (5 1/2")



Ear Curettes
 Ohrkürettchen
 Curetas auriculares
 Curettes auriculaires
 Curette auricolari

OT

OTOLGY INSTRUMENTS



Buck Ear Curette
 Straight, Blunt
 Overall Length 17 cm (6 3/4")



OT-570-01
 Fig. 000

OT-570-02
 Fig. 00

OT-570-10
 Fig. 0

OT-570-11
 Fig. 1

OT-570-12
 Fig. 2

OT-570-13
 Fig. 3

OT-570-14
 Fig. 4



Buck Ear Curette
 Angled, Blunt
 Overall Length 17 cm (6 3/4")



OT-571-01
 Fig. 000

OT-571-02
 Fig. 00

OT-571-10
 Fig. 0



OT-571-11
 Fig. 1

OT-571-12
 Fig. 2

OT-571-13
 Fig. 3



OT-571-14
 Fig. 4

Ear Curettes

 Ohrkürettchen
 Curetas auriculares
 Curettes auriculaires
 Curette auricolari

Buck Ear Curette

Straight, Sharp

Overall Length 17 cm (6 3/4")

OT-572-01
 Fig. 000

OT-572-02
 Fig. 00

OT-572-10
 Fig. 0

OT-572-11
 Fig. 1

OT-572-12
 Fig. 2

OT-572-13
 Fig. 3

OT-572-14
 Fig. 4

Buck Ear Curette

Angled, Sharp

Overall Length 17 cm (6 3/4")

OT-573-01
 Fig. 000

OT-573-02
 Fig. 00

OT-573-10
 Fig. 0

OT-573-11
 Fig. 1

OT-573-12
 Fig. 2

OT-573-13
 Fig. 3

OT-573-14
 Fig. 4

Foreign Body Instruments
 Fremdkörperinstrumente
 Instrumentos para cuerpos extraños
 Instruments pour corps étrangers
 Strumenti per corpi estranei

OT

OTOTOLOGY INSTRUMENTS



OT-576-15
Jansen Ear Scoop
 Sharp
 Overall Length 17 cm (6 3/4")



OT-580-18
Barth Ear Scoop
 Sharp
 Overall Length 18 cm (7")



OT-583-13
Gross Ear Scoop / Hook
 Blunt
 Overall Length 13 cm (5")



OT-600-10
Lempert Ear Scoop
 Sharp
 Overall Length 21.5 cm (8 1/2")

O O
OT-600-50 **OT-600-40**
 Fig. 0 Fig. 00

O O
OT-600-30 **OT-600-20**
 Fig. 000 Fig. 0000

O
OT-600-10
 Fig. 00000

Foreign Body Instruments
 Fremdkörperinstrumente
 Instrumentos para cuerpos extraños
 Instruments pour corps étrangers
 Strumenti per corpi estranei



FiberGrip™ Lempert Ear Scoop
 Sharp
 Overall Length 21.5 cm (8 1/2")



- | | |
|------------------|------------------|
| ○ | ○ |
| OT-610-10 | OT-610-20 |
| Fig. 0 | Fig. 00 |
| | |
| ○ | ○ |
| OT-610-30 | OT-610-40 |
| Fig. 000 | Fig. 0000 |
| | |
| ○ | |
| OT-610-50 | |
| Fig. 00000 | |



OT-650-01
Weber-Loch Ear Loop
 Overall Length 16.5 cm (6 1/4")



OT-660-01
 2.0 mm **OT-660-02**
 3.0 mm **OT-660-03**
 4.0 mm



Billeau Ear Loop
 Overall Length 16 cm (6 1/4")

Ear Hooklet
 Ohrhakchen
 Ganchito auricular
 Levier pour corps étrangers
 Uncino auricolare

OT

OTOTOLOGY INSTRUMENTS




OT-671-16
Hartmann Ear Hook
 Sharp
 Overall Length 16 cm (6 1/4")



OT-673-01
 Small



OT-673-02
 Medium



OT-673-03
 Large



OT-675-01
 4.0 mm
 Fig. 1



OT-675-02
 3.5 mm
 Fig. 2



OT-675-03
 3.0 mm
 Fig. 3



OT-675-04
 2.5 mm
 Fig. 4



OT-675-05
 2.0 mm
 Fig. 5



Ear Hooklet
Ohrhäkchen
Ganchito auricular
Levier pour corps étrangers
Uncino auricolare

OT-677-01
Long

OT-677-02
Short

Day Ear Hook
Overall Length 17 cm (6 3/4")



OT-679-16
Seeker
With Ball Probe, Angled 45°
Overall Length 15.5 cm (6")



OT-681-19
Maxillary Ostium Sinus Seeker
With Ball Probe, Double Ended
Ø 1.2 and 2.0 mm
Overall Length 19 cm (7 1/2")

Ear Hooklet
 Ohrhakchen
 Ganchito auricular
 Levier pour corps étrangers
 Uncino auricolare

OT

OTOTOLOGY INSTRUMENTS


OT-683-01
Sinus Seeker

Straight, Malleable Double Ended Fig. 1
 One Tip with Ball End Other Tip without
 Ball End & Blunt
 Overall Length 24 cm (9 1/2")

OT-683-02

Fig. 2, Both Sides Curved 90°
 One Tip Straight Other Tip Reverse Angled

OT-683-03

Fig. 3, Both Sides Curved 77°
 Other Side Curved 90° Both Tips Forward Angled

OT-683-04

Fig. 4, Both Sides Curved 90°
 Tips Angled Right and Left

OT-683-05

Fig. 5, Both Sides Curved 77°
 tips Angled Right and Left

OT-683-06

Fig. 6, Both Sides Curved 77°
 One Tip Straight, Other Tip Reverse Angled

Kuhn Frontal Sinus Seeker

Double Ended,
 Overall Length 22 cm (8")

Tympanum Needles
 Trommelfellnadeln
 Agujas para timpano
 Instruments pour corps étrangers
 Strumenti per corpi estranei



OT-701-16
Politzer Furuncle Knife
 Overall Length 14.5 cm (5 ¾")



OT-704-16
Politzer Tympanum Needle
 Overall Length 17.5 cm (7")



Politzer Tympanum Needle

OT-707-16
 15.5 cm (6")
 Vertical

OT-709-16
 16 cm (6 ¼")
 Horizontal



Tympanum Needles
 Trommelfellnadeln
 Agujas para tímpano
 Instruments pour corps étrangers
 Strumenti per corpi estranei

OT

OTOTOLOGY INSTRUMENTS



Lucae Tympanum Needle
 Overall Length 17.5 cm (7")

OT-711-16
Moldenhauer Foreign Body Lever
 Overall Length 15.5 cm (6")

OT-723-16
Zaufal Foreign Body Lever
 Overall Length 15 cm (6")


Mastoid Rongeurs
 Ohr-Knochenstanzen
 Pinzas cortantes para huesos
 Ostéotomes emporte-pièces
 Pinze osteotomi



Lempert Periosteal Elevator

OT-727-20
 6.0 mm
 20.5 cm (8")

OT-729-18
 6.0 mm
 18 cm (7")



Lempert Periosteal Elevator
 Overall Length 17 cm (6 3/4")




OT-731-02
 2.0 mm

OT-731-03
 3.0 mm

OT-731-04
 4.0 mm

OT-731-05
 5.0 mm



Citelli Mastoid Rongeur
 Overall Length 17 cm (6 3/4")

□
OT-751-02
 2.0 mm

□
OT-751-03
 3.0 mm

Mastoid Rongeurs, Ear Syringes
 Ohr-Knochenstanzen, Ohrspritzen
 Pinzas cortantes para huesos, Jeringas auriculares
 Ostéotomes emporte-pièces, Seringues auriculaires
 Pinze osteotomi, Siringhe auricolari

OT

OTOTOLOGY INSTRUMENTS


Beyer Mastoid Rongeur

 ♂
OT-755-19
 2.0 mm
 19 cm (7 1/2")

 ♂
OT-759-14
 1.0 mm
 14 cm (5 1/2")

Ear Syringe Complete
 With Protection Disk and 2 Tips

OT-800-00
 50 ml

OT-802-00
 100 ml



SURTEX
INSTRUMENTS

Crafted to Excellence

RH

RHINOLOGY INSTRUMENTS

Rhinologie

Rinología

Rhinologie

Rinologia

Nasal Specula
 Nasenspekula
 Espéculos nasales
 Spéculums nasaux
 Specoli nasalí



RH-011-10
Thudichum Nasal Speculum
 Set of 3 Ref: RH-013-01 to RH-013-03



RH-011-20
Thudichum Nasal Speculum
 Set of 4 Ref: RH-013-01 to RH-013-04



RH-011-30
Thudichum Nasal Speculum
 Set of 3 Ref: RH-013-05 to RH-013-07



RH-011-40
Thudichum Nasal Speculum
 Set of 7 Ref: RH-013-01 to RH-013-07

Nasal Specula
Nasenspekula
Espéculos nasales
Spécums nasaux
Specoli nasal

RH
RHINOLOGY INSTRUMENTS

RH-013-01

12 mm

Fig. 1

RH-013-03

16.5 mm

Fig. 3

RH-013-05

20.5 mm

Fig. 5

RH-013-07

23 mm

Fig. 7

RH-013-02

14 mm

Fig. 2

RH-013-04

18.5 mm

Fig. 4

RH-013-06

22 mm

Fig. 6


RH-017-01

Fig. 1

RH-017-02

Fig. 2

RH-017-03

Fig. 3

Thudichum Nasal Speculum
Voltolini Nasal Speculum

RH-021-10
Mod. Vienna Nasal Speculum

Set of 3 Ref: RH-023-01 to RH-023-03


Mod. Vienna Nasal Speculum
RH-023-01

27 mm

Fig. 1

RH-023-02

31 mm

Fig. 2

RH-023-03

35 mm

Fig. 3

Nasal Specula
Nasenspekula
Espéculos nasales
Spécums nasaux
Specoli nasalí



RH-031-10

Hartmann-Halle Nasal Speculum

Set of 3 Ref: RH-033-01 to RH-033-03




Hartmann-Halle Nasal Speculum

Overall Length 15.5 cm (6")




RH-033-01

27 mm
Fig. 1



RH-033-02

31 mm
Fig. 2



RH-033-03

35 mm
Fig. 3

Nasal Specula
 Nasenspekula
 Espéculos nasales
 Spéculums nasaux
 Specoli nasal

RH

RHINOLOGY INSTRUMENTS


RH-035-13

Martin Nasal Speculum
 17 mm
 Overall Length 13 cm (5")

RH-037-13

Martin Nasal Speculum
 17.0 mm, Bent Sideways
 Overall Length 12.5 cm (5")

RH-041-13

Tieck-Halle Nasal Speculum
 16.0 mm
 Overall Length 13 cm (5")



Nasal Specula
Nasenspekula
Espéculos nasales
Spécums nasaux
Specoli nasalí



RH-051-10
Killian Nasal Speculum
Set of 4 Ref; RH-053-01 to RH-053-04




Killian Nasal Speculum
Overall Length 13 cm (5")




RH-053-01
Fig. 1, 35 mm



RH-053-02
Fig. 2, 50 mm



RH-053-03
Fig. 3, 75 mm



RH-053-04
Fig. 4, 85 mm

Nasal Specula
 Nasenspekula
 Espéculos nasales
 Spéculums nasaux
 Specoli nasal

RH
RHINOLOGY INSTRUMENTS

RH-061-10

Cottle Nasal Speculum
 Set of 4 Ref; RH-063-01 to
 RH-063-04

Cottle Nasal Speculum

Overall Length 13.5 cm (5 1/4")


RH-063-01

Fig. 1, 35 mm

RH-063-02

Fig. 2, 50 mm


RH-063-03

Fig. 3, 75 mm

RH-063-04

Fig. 4, 85 mm
 14.5 cm (5 3/4")

Nasal Specula
 Nasenspekula
 Espéculos nasales
 Spéculums nasaux
 Specoli nasalí



RH-065-10
Cottle Nasal Speculum
 With Screw Fixation
 Set of 4 Ref, RH-067-01
 to RH-067-04



Cottle Nasal Speculum
 With Screw Fixation
 Overall Length 13.5 cm (5 1/4")

RH-067-01
 Fig. 1, 35 mm

RH-067-02
 Fig. 2, 50 mm

RH-067-03
 Fig. 3, 75 mm

RH-067-04
 Fig. 4, 85 mm

Nasal Specula, Polypus Snares
 Nasenspekula, Polypenschnüre
 Espéculos nasales, Polipótomas
 Spécums nasaux, Serre-noeuds pour polypes
 Specoli nasali, Serranodi per polipi

RH
RHINOLOGY INSTRUMENTS


Cottle Nasal Speculum
Black



RH-073-45

45 mm
 13.5 cm (5 1/4")



RH-073-75

75 mm
 14 cm (5 1/2")



RH-091-26

Krause Polypus Snare

Ø 3.0 mm
 Overall Length 26 cm (10 1/4")

Polypus Snare
 Polypschnüre
 Polipótomas
 Serre-noeuds pour polypes
 Serranodi per polipi



RH-093-26
Krause-Voss Polypus Snare
 \varnothing 3.0 mm
 Overall Length 26 cm (10 $\frac{1}{4}$ ")



RH-095-24
Lange-Wilde Polypus Snare
 Overall Length 24 cm (9 $\frac{1}{2}$ ")



RH-097-12
Polypus Wire
 Bundle with 12 Pieces, \varnothing 0.3
 Overall Length 10.5 cm (4")

Nasal Tampon Forceps
 Nasentamponpinzetten
 Pinzas para taponamiento
 Pince à pansements
 Pinze per medicatura


RH
RHINOLOGY INSTRUMENTS


RH-101-15
Troeltsch Nasal Tampon Forcep
 Overall Length 14.5 cm (5 3/4")

RH-103-18
Troeltsch Nasal Tampon Forcep
 Overall Length 18 cm (7")

RH-105-21
Gruenwald (Jansen) Nasal Tampon Forcep
 Overall Length 21.5 cm (8 1/2")

Nasal Tampon Forceps
 Nasentamponpinzetten
 Pinzas para taponamiento
 Pince à pansements
 Pinze per medicatura



RH-107-16
Nasal Tampon Forcep
 Smooth Jaws
 Overall Length 16 cm (6 $\frac{1}{4}$ ")

RH-109-16
Nasal Tampon Forcep
 1 x 2 Teeth
 Overall Length 16 cm (6 $\frac{1}{4}$ ")

RH-111-18
Nasal Tampon Forcep
 Smooth Jaws
 Overall Length 18.5 cm (7 $\frac{1}{4}$ ")

RH-113-20
Nasal Tampon Forcep
 1 x 2 Teeth
 Overall Length 20 cm (8")

Nasal Tampon Forceps
 Nasentamponpinzetten
 Pinzas para taponamiento
 Pinces à pansements
 Pinze per medicature



RH-121-18
Westmacott Nasal Polypus Forcep
 Shaft Length 18.5 cm (7 1/4")

RH-123-20
Lubet-Barbon Nasal Polypus Forcep
 Shaft Length 20 cm (8")



RH-127-16
Hartmann Nasal Polypus Forcep
 Shaft Length 15 cm (6")

RH-129-18
Heymann-Knight Nasal Polypus Forcep
 Shaft Length 18 cm (8")

RH

RHINOLOGY INSTRUMENTS

Nasal Tampon Forceps
 Nasentamponpinzetten
 Pinzas para taponamiento
 Pince à pansements
 Pinze per medicatura



RH-131-20
Hartmann Nasal Polypus Forcep
 Shaft Length 18.5 cm (7 1/4")

RH-133-18
Littauer Nasal Polypus Forcep
 Shaft Length 18 cm (7")

RH-135-16
Killian Nasal Polypus Forcep
 Shaft Length 16 cm (6 1/4")

Septum Forceps
 Septumzangen
 Pinzas para tabique nasal
 Pince pour la cloison nasale
 Pinze per il setto nasale

RH
RHINOLOGY INSTRUMENTS

RH-140-01
Zurich Nasal Scissor

 Straight, Large Model, Cutting Length 18 mm
 Shaft Length 14 cm (4 1/2")

RH-161-01

Fig. 1

RH-161-02

Fig. 2


RH-161-03

Fig. 3

Heymann Septum Forcep
 Shaft Length 18 cm (7")

RH-165-17
Knight Septum Forcep
 Shaft Length 17 cm (6 3/4")

RH-167-16
Craig Septum Forcep
 Shaft Length 16 cm (6 1/4")



Ethmoid Bone Forceps
 Siebbeinzangen
 Pinzas para etmoides
 Pince éthmoïdales
 Pinze etmoidali

RH
RHINOLOGY INSTRUMENTS


B鲁恩尼氏鼻中隔成形鉗
 Shaft Length 16.5 cm (6 1/2")

RH-175-01

Fig. 1

RH-175-02

Fig. 2

RH-175-03

Fig. 3

RH-175-04

Fig. 4



海威斯耳咽管夹持鉗
 Shaft Length 12 cm (4 3/4")



RH-177-01

Jaw Curved Downwards
 Fixed Jaw Curved 90° Moveable Jaw
 Opening upto 120° Backward Opening



RH-177-02

Jaw Curved Upwards
 Fixed Jaw Curved 90° Moveable Jaw
 Opening upto 120° Backward Opening



RH-177-03

Jaw Curved Right
 Fixed Jaw Curved 90° Moveable Jaw
 Opening upto 120° Backward Opening



RH-177-04

Jaw Curved Left
 Fixed Jaw Curved 90° Moveable Jaw
 Opening upto 120° Backward Opening




RH-177-05

With Extra Long Curved For Alveolar
 Recess Fixed Jaw Curved 115° Moveable
 Jaw Opening upto 140° Backward Opening




RH-177-06

With Fixed Jaw Curved 140°
 Moveable Jaw Backward Opening
 155°



Nasal Cutting Forceps
 Scharfe Nasenzangen
 Pinzas nasales cortantes
 Pince nasale coupante
 Pinze taglienti per il naso

RH
RHINOLOGY INSTRUMENTS


Weil-Blakesley Nasal Cutting Forcep

Shaft Length 12 cm (4 ¾")



RH-201-01

Fig. 1
Straight
3.0 mm



RH-201-02

Fig. 2
Straight
3.5 mm



RH-201-03

Fig. 3
Straight
4.5 mm



RH-201-04

Fig. 4
Straight
5.0 mm



RH-201-05

Fig. 5
Straight
5.5 mm



RH-203-01

Fig. 1
Angled 45°
3.0 mm



RH-203-02

Fig. 2
Angled 45°
3.5 mm



RH-203-03

Fig. 3
Angled 45°
4.5 mm



RH-203-04

Fig. 4
Angled 45°
5.0 mm



RH-203-05

Fig. 5
Angled 45°
5.5 mm



RH-205-01

Fig. 1
Angled 90°
3.0 mm



RH-205-02

Fig. 2
Angled 90°
3.5 mm



RH-205-03

Fig. 3
Angled 90°
4.5 mm



RH-205-04

Fig. 4
Angled 90°
5.0 mm





Nasal Cutting Forceps
Scharfe Nasenzangen
Pinzas nasales cortantes
Pinces nasales coupantes
Pinze taglienti per il naso



Tilley-Henckel Nasal Cutting Forcep
Shaft Length 12 cm (4 3/4")

RH-207-01

Straight
4.5 mm
Fig. 5

RH-207-02

Straight
5.5 mm
Fig. 5

RH-207-03

Straight
6.5 mm
Fig. 5

Struycken Nasal Cutting Forcep
Shaft Length 12 cm (4 3/4")

RH-231-01
0.9 mm
Fig. 1

RH-231-02
1.2 mm
Fig. 2



Nasal Cutting Forcep
Shaft Length 12 cm (4 3/4")



RH-233-01
Through Cut, Straight Slender Model
with Flat Jaws Bite Size 3 x 4 mm



RH-233-02
Through Cut with Basket Straight Slender
Model with Flat Jaws Bite Size 3 x 4 mm



RH-233-03
Through Cut Straight Fine Model
with Flat Jaws Bite Size 1.25 x 5 mm



RH-233-04
Through Cut Straight Small Robust Model
with Flat Jaws Bite Size 2.5 x 6.5 mm


Nasal Cutting Forceps
Scharfe Nasenzangen
Pinzas nasales cortantes
Pince nasale coupante
Pinze taglienti per il naso

RH

RHINOLOGY INSTRUMENTS



RH-234-00
Stammberger Nasal Cutting Forcep
Retrograde 360° Rotatable
Bite Size 2.0 mm
Shaft Length 12 cm (4 3/4")



Stammberger Nasal Cutting Forcep
Bite Size 2.0 mm
Shaft Length 12 cm (4 3/4")

RH-235-01
Retrograde Up Cutting

RH-235-02
Retrograde Right Cutting

RH-235-03
Retrograde Left Cutting



Nasal Cutting Forceps
Scharfe Nasenzangen
Pinzas nasales cortantes
Pinces nasales coupantes
Pinze taglienti per il naso



Stammberger Circular Cutting Punch
Shaft Length 18 cm (7")



Stammberger Circular Cutting Punch
Shaft Length 18 cm (7")

Nasal Cutting Forceps
 Scharfe Nasenzangen
 Pinzas nasales cortantes
 Pince nasale coupante
 Pinze taglienti per il naso

RH
RHINOLOGY INSTRUMENTS


Hartmann Nasal Cutting Forcep (Conchotome)
 Shaft Length 12 cm (4 3/4")


RH-241-01
 Ø 5.0 mm
 Fig. 1


RH-241-02
 Ø 7.0 mm
 Fig. 2


RH-241-03
 Ø 9.0 mm
 Fig. 3


RH-241-04
 Ø 11.0 mm
 Fig. 4



RH-245-12
Schmeden Nasal Cutting Forcep
 8.0 mm
 Shaft Length 12 cm (4 3/4")

Nasal Cutting Forceps
 Scharfe Nasenzangen
 Pinzas nasales cortantes
 Pinces nasales coupantes
 Pinze taglienti per il naso



Blakesley Nasal Cutting Forcep
 Shaft Length 12 cm (4 3/4")



Blakesley Nasal Cutting Forcep
 Shaft Length 12 cm (4 3/4")



RH-251-25
 Through Cut Straight
 Bite Size 2.5 mm



RH-251-35
 Through Cut Straight
 Bite Size 3.5 mm



RH-251-42
 Through Cut Straight
 Bite Size 4.2 mm



RH-253-35
 Through Cut Curved
 Right Bite Size 3.5 mm



RH-251-48
 Through Cut Straight
 Bite Size 4.8 mm



RH-255-35
 Through Cut Curved
 Left Bite Size 3.5 mm



RH-257-25
 Through Cut
 Angled 45° Upwards
 Bite Size 2.5 mm



RH-257-35
 Through Cut
 Angled 45° Upwards
 Bite Size 3.5 mm



RH-257-42
 Through Cut
 Angled 45° Upwards
 Bite Size 4.2 mm



RH-257-48
 Through Cut
 Angled 45° Upwards
 Bite Size 4.8 mm

9
 220
 210
 8
 200
 190
 7
 180
 6
 170
 5
 160
 4
 150
 3
 140
 2
 130
 1
 120
 0
 110
 100
 90
 80
 70
 60
 50
 40
 30
 20
 10
 MILLIMETERS
 INCHES

Nasal Scissors
 Nasenscheren
 Tijeras nasales
 Ciseaux nasaux
 Forbici nasal

RH

RHINOLOGY INSTRUMENTS


RH-261-00

 Vertical Straight
 Cup Size 2.0 mm

RH-261-01

 Vertical Angled 45°
 Cup Size 2.0 mm

RH-261-02

 Vertical Angled 55°
 Cup Size 2.0 mm

RH-261-03

 Vertical Angled 70°
 Cup Size 2.0 mm

RH-261-04

 Vertical Angled 90°
 Cup Size 2.0 mm

RH-261-05

 Vertical Angled 110°
 Cup Size 2.0 mm

Giraffe Cup Shaped Biopsy Punch

Shaft Length 18 cm (7")


RH-261-11

Straight

RH-261-13

Curved Right

RH-261-14

Curved Left


Nasal Sinus Scissor

Shaft Length 12 cm (4 3/4")



Nasal Cutting Forceps
Scharfe Nasenzangen
Pinzas nasales cortantes
Pinces nasales coupantes
Pinze taglienti per il naso

RH-263-00

Vertical Straight
Cup Size 3.0 mm

RH-263-01

Vertical Angled 45°
Cup Size 3.0 mm

RH-263-02

Vertical Angled 55°
Cup Size 3.0 mm

RH-263-04

Vertical Angled 90°
Cup Size 3.0 mm

RH-261-03

Vertical Angled 70°
Cup Size 3.0 mm

RH-263-05

Vertical Angled 110°
Cup Size 3.0 mm

Giraffe Cup Shaped Biopsy Punch

Shaft Length 18 cm (7")



Nasal Cutting Forceps
 Scharfe Nasenzangen
 Pinzas nasales cortantes
 Pinces nasales coupantes
 Pinze taglienti per il naso

RH
RHINOLOGY INSTRUMENTS

RH-265-00

 Vertical Straight
 Cup Size 4.0 mm

RH-265-01

 Vertical Angled 45°
 Cup Size 4.0 mm

RH-265-02

 Vertical Angled 55°
 Cup Size 4.0 mm

RH-265-03

 Vertical Angled 70°
 Cup Size 4.0 mm

RH-265-04

 Vertical Angled 90°
 Cup Size 4.0 mm

RH-265-05

 Vertical Angled 110°
 Cup Size 4.0 mm

Giraffe Cup Shaped Biopsy Punch

Shaft Length 18 cm (7")


RH-267-00

 Horizontal Straight
 Cup Size 2.0 mm

RH-267-01

 Horizontal Angled 45°
 Cup Size 2.0 mm

RH-267-02

 Horizontal Angled 55°
 Cup Size 2.0 mm

RH-267-03

 Horizontal Angled 70°
 Cup Size 2.0 mm

RH-267-04

 Horizontal Angled 90°
 Cup Size 2.0 mm

RH-267-05

 Horizontal Angled 110°
 Cup Size 2.0 mm

Giraffe Cup Shaped Biopsy Punch

Shaft Length 18 cm (7")



Nasal Cutting Forceps
 Scharfe Nasenzangen
 Pinzas nasales cortantes
 Pinces nasales coupantes
 Pinze taglienti per il naso



RH-269-00
 Horizontal Straight
 Cup Size 3.0 mm

RH-269-01
 Horizontal Angled 45°
 Cup Size 3.0 mm


RH-269-02
 Horizontal Angled 55°
 Cup Size 3.0 mm

RH-269-03
 Horizontal Angled 70°
 Cup Size 3.0 mm

RH-269-04
 Horizontal Angled 90°
 Cup Size 3.0 mm

RH-269-05
 Horizontal Angled 110°
 Cup Size 3.0 mm

Giraffe Cup Shaped Biopsy Punch
 Shaft Length 18 cm (7")



RH-271-00
 Horizontal Straight
 Cup Size 4.0 mm

RH-271-01
 Horizontal Angled 45°
 Cup Size 4.0 mm

RH-271-02
 Horizontal Angled 55°
 Cup Size 4.0 mm

RH-271-03
 Horizontal Angled 70°
 Cup Size 4.0 mm

RH-271-04
 Horizontal Angled 90°
 Cup Size 4.0 mm

RH-271-05
 Horizontal Angled 110°
 Cup Size 4.0 mm

Giraffe Cup Shaped Biopsy Punch
 Shaft Length 18 cm (7")


Septum Forceps, Septum Morselizer
 Septumzangen, Septumquetschzange
 Pinzas para tabique nasal, Pinza para aplastamiento
 Pince pour la cloison nasale, Pince-gouge pour la cloison nasale
 Pinze per il setto nasale, Pinza contundente per il setto nasale

RH
RHINOLOGY INSTRUMENTS


RH-271-15
Cottle Masing Needle Holder
 Smooth Jaws
 Shaft Length 15 cm (6")



Jansen-Middleton Septum Forcep
 Overall Length 19 cm (7 1/2")



RH-283-19
Through
Cutting Jaws

RH-291-21
UltraGrip™ TC Rubin Septum Forcep
 Overall Length 21 cm (8 1/4")



RH-293-21
UltraGrip™ TC Lower Lateral Morselizer
 Overall Length 19 cm (7 1/2")

ULTRA GRIP™
ULTRA GRIP™

Sphenoid Bone Punches
 Keilbeinstanzen
 Pinzas sacabocadas para esfenoides
 Pince emporte-pièces pour sinus sphénoidal
 Pinze taglienti per seno sfenoidale



Hajek Sphenoid Punch
Overall Length 14 cm (5 1/2")



RH-301-01

Up Cutting
4.0 mm x 4.0 mm



RH-301-02

Down Cutting
4.0 mm x 4.0 mm



Hajek-Kofler Sphenoid Punch
Overall Length 14 cm (5 1/2")



RH-303-01

Up Cutting
3.5 mm x 3.5 mm



RH-303-02

Down Cutting
3.5 mm x 3.5 mm

Sphenoid Bone Punches

Keilbeinstanzen

Pinzas sacabocadas para esfenoides

Pince emporte-pièces pour sinus sphénoidal

Pinze taglienti per seno sfenoidale



Kerrison Sphenoid Punch
Up Cutting

RH-305-03

Ø 3.0 mm x 3.0 mm
9 cm (3 ½")

RH-305-04

Ø 4.0 mm x 4.0 mm
9 cm (3 ½")

RH-305-05

Ø 5.0 mm x 5.0 mm
9 cm (3 ½")

RH-305-06

Ø 6.0 mm x 6.0 mm
9 cm (3 ½")


RH-307-20

Hajek-Claus Antrum Punch
Bite Size 5.0 mm
Overall Length 20 cm (8")


Antrum Punch

Bite Size 4.0 mm x 5.0 mm
Overall Length 14 cm (5 ½")


RH-309-01

Angled 45° Backward Cutting


RH-309-02

Angled 45° Forward Cutting

RH
RHINOLOGY INSTRUMENTS

Rhinoplastic Knives
 Rhinoplastik-Messer
 Cuchillos rinoplásticos
 Couteaux pour plastie du nez
 Bisturi per plastica nasale



RH-350-15
 Straight
 Fig. 1




RH-351-15
 Curved
 Fig. 2

Joseph Rhinoplasty Knife
 Overall Length 15 cm (6")



RH-355-15
Fomon Rhinoplasty Knife
 Overall Length 16.5 cm (6 1/2")



RH-361-16
Converse Rhinoplasty Knife
 Overall Length 16 cm (6 1/4")



Rhinoplastic Knives,
 Rhinoplastik-Messer,
 Cuchillos rinoplásticos,
 Couteaux pour plastie du nez,
 Bisturi per plastica nasale,

RH
RHINOLOGY INSTRUMENTS

RH-365-15

Joseph Rhinoplasty Knife
 Overall Length 15 cm (6")


RH-367-15

Cottle Rhinoplasty Knife
 Overall Length 14 cm (5 1/2")


RH-369-15

Freer Mucosa Knife
 Overall Length 16.5 cm (6 1/2")


RH-371-01
 Right

RH-371-02
 Left

Freer Septum Knife

Overall Length 18 cm (7")

Septum Knives, Swivel Knives
 Septummesser, Schwingmesser
 Cuchillos para tabique nasal
 Bistouris pour la cloison nasale
 Bisturi per il setto nasale



Freer Septum Knife

RH-373-01
18 cm (7")

RH-373-02
19 cm (7 1/2")



Ballenger Swivel Knife


RH-380-03
3.0 mm

RH-380-04
4.0 mm

RH-380-05
5.0 mm



Ballenger Swivel Knife
 Bayonet Shaped
 Overall Length 21 cm (8 1/4")




RH-381-03
3.0 mm

RH-381-04
4.0 mm

RH-381-05
5.0 mm



Septum Elevators
 Septumelevatorien
 Elevadores para tabique nasal
 Élévateurs pour la cloison nasale
 Elevatori per il setto nasale

RH
RHINOLOGY INSTRUMENTS


McKenty Septum Elevator
 Overall Length 15 cm (6")

RH-455-16
Halle Septum Elevator
 3.5 mm
 Overall Length 16 cm (6 1/4")

RH-461-16
Joseph Septum Elevator
 4.0 mm
 Overall Length 16 cm (6 1/4")

Septum Elevators
 Septumelevatorien
 Elevadores para tabique nasal
 Elévateurs pour la cloison nasale
 Elevatori per il setto nasale



Daniel Periosteal Elevator
 Overall Length 16 cm (6 1/4")

RH-465-01
 Slightly Curved

RH-465-02
 Medium Curved

RH-465-03
 Strongly Curved



RH-471-19
 Sharp
 5.0 mm


RH-473-19
 Blunt
 6.0 mm

Freer Septum Elevator
 Overall Length 19 cm (7 1/2")



RH-475-18
Freer Septum Elevator
 Double Ended 5.0 mm, 5.0 mm
 Overall Length 18.5 cm (7 1/4")

Septum Elevators
 Septumelevatorien
 Elevadores para tabique nasal
 Élévateurs pour la cloison nasale
 Elevatori per il setto nasale

RH
RHINOLOGY INSTRUMENTS

RH-476-18

Freer-Yasargil Septum Elevator
 Double Ended, Sharp/Blunt
 3.0 mm, 3.0 mm
 Overall Length 18.5 cm (7 1/4")


RH-477-20

Freer Septum Elevator
 Double Ended, Sharp/Blunt
 6.0 mm, 5.0 mm
 Overall Length 20.5 cm (8")


RH-479-20

McKenty Septum Elevator
 4.0 mm, 4.0 mm
 Overall Length 20.5 cm (8")


RH-482-19

Ballenger-Hajek Septum Elevator
 6.0 mm, 6.0 mm
 Overall Length 19 cm (7 1/2")

Septum Elevators
 Septumelevatorien
 Elevadores para tabique nasal
 Élevateurs pour la cloison nasale
 Elevatori per il setto nasale



RH-501-23
Cottle Septum Elevator
 5.0 mm, 5.0 mm
 Overall Length 23.5 cm (9 1/4")



RH-503-18
Killian Septum Elevator
 4.0 mm, 4.0 mm
 Overall Length 18 cm (7")



RH-505-23
Killian Septum Elevator
 4.5 mm, 4.5 mm
 Overall Length 23 cm (9")



RH-507-21
Pierce Septum Elevator
 4.0 mm, 4.0 mm
 Overall Length 21.5 cm (8 1/2")

Periosteal Elevators
 Periostelevatorien
 Legras para periostio
 Elévateurs pour le périoste
 Elevatori per il periosto

RH
RHINOLOGY INSTRUMENTS


RH-509-21
Pennington Septum Elevator
 4.0 mm, 4.0 mm
 Overall Length 21.5 cm (8 1/2")



RH-511-19
Roger Septum Elevator
 4.0 mm, 4.0 mm
 Overall Length 19.5 cm (7 3/4")



RH-513-15
Traquair Septum Elevator
 1 x 3 mm (7 mm Long) &
 1 x 2 mm (5 mm Long)
 Overall Length 15.5 cm (6")



RH-541-22
Cottle Septum Elevator
 3.3 mm, 3.7 mm
 Overall Length 23 cm (9")

Periosteal Elevators
 Periostelevatorien
 Legras para periostio
 Élevateurs pour le périoste
 Elevatori per il periosto



RH-543-21
Masing Septum Elevator
 3.0 mm, 2.0 mm
 Overall Length 19.5 cm (7 3/4")



RH-545-21
Cottle Septum Elevator
 4.5 mm, 4.5 mm
 Overall Length 19.5 cm (7 3/4")



RH-547-21
Howarth Septum Elevator
 5.5 mm, 4.5 mm
 Overall Length 21.5 cm (8 1/2")



RH-551-16
Joseph Periosteal Elevator
 6.0 mm
 Overall Length 18 cm (7")

Septum Chisels, Septum Gouges
 Septummeißel, Septum-Hohlmeißel
 Cinceles para tabique nasal, Gubias para tabique nasal
 Elévateurs pour le périoste, Gouges pour la cloison nasale
 Elevatori per il periostio, Sgorbie per il setto nasale

RH
RHINOLOGY INSTRUMENTS


RH-600-16
Straight



RH-601-16
Curved



RH-604-16
West Septum Chisel
Overall Length 16 cm (6 1/4")



RH-605-16
West Septum Gouge
Overall Length 16 cm (6 1/4")



RH-608-15
Freer Septum Gouge
6.0 mm
Overall Length 15 cm (6")



Freer Septum Chisel
4.0 mm
Overall Length 16.5 cm (6 1/2")

Septum Gouges
 Septum-Hohlmeißel
 Gubias para tabique nasal
 Gouges pour la cloison nasale
 Sgorbie per il setto nasale



Ballenger Septum Gouge
 Overall Length 15 cm (6")



RH-610-02
 2.0 mm



RH-610-04
 4.0 mm



RH-610-06
 6.0 mm



RH-610-08
 8.0 mm



RH-613-16
Killian Septum Gouge
 Bayonet Shaped, 8.0 mm
 Overall Length 16 cm (6 1/4")



RH-615-16
Killian-Claus Septum Gouge
 Bayonet Shaped, 5.0 mm
 Overall Length 16 cm (6 1/4")



Osteotomes
Osteotome
Osteótomas
Ostéotomes
Osteotomi

RH
RHINOLOGY INSTRUMENTS


RH-620-02	RH-620-04	RH-620-07
2.0 mm	4.0 mm	7.0 mm
RH-620-09	RH-620-12	RH-620-14
9.0 mm	12.0 mm	14.0 mm

Cottle Nasal Osteotome
Straight
Overall Length 18.5 cm (7 1/4")



RH-621-02	RH-621-04	RH-621-06
2.0 mm	4.0 mm	6.0 mm
RH-621-09	RH-621-12	
9.0 mm	12.0 mm	

Cottle Nasal Osteotome
Curved
Overall Length 18.5 cm (7 1/4")

Chisels, Osteotomes
 Meißel, Osteotome
 Cinceles, Osteótomas
 Ciseaux-burins, Ostéotomes
 Scalpelli, Osteotomi



Cottle Nasal Chisel
 Straight, Cross Handle
 Overall Length 18.5 cm (7 1/4")

RH-626-02
 2.0 mm

RH-626-04
 4.0 mm

RH-626-07
 7.0 mm

RH-626-09
 9.0 mm

RH-626-12
 12.0 mm



RH-627-06
Cottle Nasal Chisel
 Curved, Cross Handle
 6.0 mm
 Overall Length 18.5 cm (7 1/4")



Cottle Nasal Osteotome
 Overall Length 18.5 cm (7 1/4")



RH-630-14
 14.0 mm


RH-630-16
 16.0 mm

RH-630-18
 18.0 mm

Osteotomes
Osteotome
Osteótomas
Ostéotomes
Osteotomi

RH
RHINOLOGY INSTRUMENTS


Cottle Nasal Osteotome
Cross Handle
Overall Length 18.5 cm (7 1/4")



RH-632-04
4.0 mm

RH-632-06
16.0 mm

RH-632-07
7.0 mm

RH-632-16
16.0 mm



Rubin Nasal Osteotome
Overall Length 18.5 cm (7 1/4")



RH-636-10
10.0 mm

RH-636-12
12.0 mm

RH-632-14
14.0 mm

RH-632-16
16.0 mm




McIndoe Nasal Chisel
Overall Length 15.5 cm (6")



Chisels, Osteotomes
Meißel, Osteotome
Cinceles, Osteótomos
Ciseaux-burins, Ostéotomes
Scalpelli, Osteotomi



RH-640-03
3.0 mm



RH-640-04
4.5 mm



Walter Nasal Osteotome
Overall Length 19 cm (7 1/2")



RH-642-03
3.0 mm

RH-642-04
4.0 mm

RH-642-07
7.0 mm

RH-642-09
9.0 mm

RH-642-12
12.0 mm



Halle Antrum Curette
Malleable
Overall Length 21 cm (8 1/4")



RH-643-01
4.8 mm



RH-643-02
5.8 mm



RH-643-03
6.8 mm

Chisels, Osteotomes
 Meißel, Osteotome
 Cinceles, Osteótomas
 Ciseaux-burins, Ostéotomes
 Scalpelli, Osteotomi

RH
RHINOLOGY INSTRUMENTS


Antrum Curette
 Oval
 Overall Length 19 cm (7 1/2")



RH-645-01
 7 x 11 mm



RH-645-02
 8 x 11 mm



RH-645-03
 9 x 13 mm



Antrum Curette
 Overall Length 19 cm (7 1/2")



RH-647-01
 Round
 Ø 4 mm



RH-647-02
 Oblong Small
 1.5 x 6 mm



RH-647-03
 Oblong Large
 3 x 9 mm



RH-649-01
 Angled
 Forward Cutting
 7 x 10 mm



RH-649-02
 Angled
 Backwards Cutting
 7 x 10 mm



Osteotomes
Osteotome
Osteótomas
Ostéotomes
Osteotomi



Fomon Nasal Chisel
Overall Length 16 cm (6 1/4")

RH-648-04	4.0 mm
RH-648-05	5.0 mm
RH-648-06	6.0 mm
RH-648-07	7.0 mm




Cinelli Nasal Osteotome
Overall Length 16 cm (6 1/4")

RH-650-10	10.0 mm
RH-650-12	12.0 mm
RH-650-14	14.0 mm
RH-650-16	16.0 mm


Osteotomes
Osteotome
Osteótomas
Ostéotomes
Osteotomi

RH
RHINOLOGY INSTRUMENTS


RH-654-07
Straight
7.0 mm



RH-655-07
Left
7.0 mm



RH-657-07
Right
7.0 mm



Cottle Nasal Chisel
Overall Length 18 cm (7")



RH-660-02
2.0 mm



RH-660-04
4.0 mm



RH-660-06
6.0 mm



RH-660-07
7.0 mm



RH-660-09
9.0 mm



RH-660-12
12.0 mm

Neivert-Anderson Nasal Osteotome
Overall Length 20.5 cm (8")

Osteotomes
Osteotome
Osteótomas
Ostéotomes
Osteotomi



RH-664-01
1.0 mm

RH-664-02
2.0 mm

RH-664-03
3.0 mm

RH-664-04
4.0 mm

RH-664-05
5.0 mm

RH-664-06
6.0 mm

RH-664-07
7.0 mm

RH-664-08
8.0 mm

RH-664-10
10.0 mm

RH-664-12
12.0 mm


Sheehan Nasal Osteotome
Overall Length 15 cm (6")

Rhinoplastic Chisels
 Nasenmeißel, Meißel
 Cinceles nasales, Cinceles
 Ciseaux-burins nasaux, Ciseaux-burins
 Scalpelli nasali, Scalpelli


RH
RHINOLOGY INSTRUMENTS


SURTEX CE


RH-666-05
 Straight
 5.0 mm




RH-667-05
 Right
 5.0 mm



RH-669-05
 Left
 5.0 mm



RH-672-01
 Parkes Nasal Osteotome
 Overall Length 16 cm (6 1/4")



Silver Nasal Chisel
 Overall Length 18 cm (7")

Nasal Saws
 Nasensägen
 Sierras nasales
 Scies nasales
 Seghe nasalì



RH-691-01
Joseph Nasal Saw
 Angled to Left
 Overall Length 17 cm (6 $\frac{3}{4}$ ")

RH-691-02
Joseph Nasal Saw
 Angled to Right
 Overall Length 17 cm (6 $\frac{3}{4}$ ")

Nasal Saws
Nasensägen
Sierras nasales
Scies nasales
Seghe nasalí

RH
RHINOLOGY INSTRUMENTS

RH-692-19
Joseph Nasal Saw

Straight, 40.0 mm
Overall Length 19 cm (7 1/2")


RH-693-01
Joseph Nasal Saw

Bayonet Shaped, 40.0 mm
Overall Length 19 cm (7 1/2")


RH-693-02
Joseph Nasal Saw

Bayonet Shaped, 40.0 mm
Overall Length 19 cm (7 1/2")


RH-698-00
Cottle Cartilage Crusher



Nasal Rasps
Nasenraspeln
Escofinas nasales
Râpes nasales
Raspe nasali



Miller-Colburn Nasal Rasp
Drawing Cut
Overall Length 18 cm (7")




RH-701-01
Fig. 1




RH-701-02
Fig. 2



RH-701-03
Fig. 3



RH-701-04
Fig. 4



RH-701-05
Fig. 5



Miller-Colburn Nasal Rasp
Cross Cut
Overall Length 18 cm (7")



RH-703-01
Fig. 1



RH-703-02
Fig. 2



RH-703-03
Fig. 3



RH-703-04
Fig. 4

Nasal Rasps
Nasenraspeln
Escopinas nasales
Râpes nasales
Raspe nasali


RH
RHINOLOGY INSTRUMENTS

RH-705-01

UltraCut™ TC Miller Nasal Rasp
Drawing Cut
Overall Length 19 cm (7 1/2")


RH-707-01

UltraCut™ TC Miller Nasal Rasp
Cross Cut
Overall Length 19 cm (7 1/2")


RH-710-18


Maltz Nasal Rasp
Drawing Cut, 40 x 9 mm
Overall Length 17.5 cm (7")

ULTRACut™
ULTRACut™

Nasal Rasps
Nasenraspeln
Escofinas nasales
Râpes nasales
Raspe nasali




RH-712-01
34 x 7 mm
Drawing Cut



RH-712-02
34 x 7 mm
Pushing Cut



RH-714-01
60 x 8 mm



RH-714-02
45 x 9 mm

Maltz Nasal Rasp
Overall Length 17 cm (6 3/4")

Joseph Nasal Rasp
Overall Length 16 cm (6 1/4")

Nasal Rasp
Nasenraspeln
Escopinas nasales
Râpes nasales
Raspe nasali

RH
RHINOLOGY INSTRUMENTS


Fomon Nasal Rasp
Overall Length 21 cm (8 1/4")




RH-720-01
30 x 8 mm
Fine




RH-720-02
30 x 8 mm
Coarse



Parkes Nasal Rasp
Drawing Cut
Overall Length 19.5 cm (7 3/4")



RH-722-01
35 x 8 mm



RH-722-02
35 x 10 mm



Gallaher Antrum Rasp
Overall Length 17 cm (6 3/4")

RH-724-01

Straight
45 x 6 mm

RH-724-02

Left
45 x 6 mm

RH-724-03

Right
45 x 6 mm


Nasal Rasps
Nasenraspeln
Escofinas nasales
Râpes nasales
Raspe nasali





RH-726-21
Cottle Nasal Rasp
30 x 8 mm
Overall Length 20.5 cm (8")



Aufrecht Glabella Nasal Rasp
Overall Length 21 cm (8 1/4")



RH-729-01
Pushing Cut
32 x 9 mm



RH-729-02
Drawing Cut
32 x 9 mm




RH-730-21
UltraCut™ TC Webster Glabella Nasal Rasp
Bi-Directional, (Push-Pull Cutting)
Overall Length 21 cm (8 1/4")

ULTRACut™

Nasal Rasps
Nasenraspeln
Escopinas nasales
Râpes nasales
Raspe nasal

RH
RHINOLOGY INSTRUMENTS


Lewis Nasal Rasp
Overall Length 19 cm (7 1/2")



RH-732-02
Drawing
15 x 8 mm



RH-734-01
Drawing Cut,
Fine Cutting Edge
of 15 x 8 mm


RH-734-02
Drawing Cut,
Coarse Cutting Edge
of 15 x 8 mm



RH-735-18
Wiener Nasal Rasp
Overall Length 18 cm (7")

ULTRACUT™

Bone Files, Nasal Rasps
 Knochenfeilen, Nasenraspeln
 Limas para huesos, Escofinas nasales
 Limes à os, Râpes nasales
 Lime per ossa, Raspe nasali



Diam-n-Dust™ Nasal Rasp
 Overall Length 19 cm (7 1/2")

RH-738-19
 Straight


RH-739-19
 Curved



UltraCut™ TC Maltz Nasal Rasp
 Overall Length 19.5 cm (7 3/4")

ULTRA*Cut*™

 **Diam-n-GRIP**



RH-742-01

RH-742-02

RH-742-03

RH-742-04

RH-742-05

RH-742-06


RH-742-07

RH-742-08

RH-742-09

RH-742-10


Nasal Rasps, Antrum Rasps
 Nasenraspeln, Antrumraspeln
 Escofinas nasales, Escofinas para el antró
 Râpes nasales, Râpes pour antre
 Raspe nasali, Raspe per seno mascellare

RH
RHINOLOGY INSTRUMENTS


RH-744-21
Fomon Nasal Rasp
 30 x 8 mm
 Overall Length 20.5 cm (8")

UltraCut™ Fomon Nasal Rasp
 Overall Length 20.5 cm (8")

ULTRACut™



RH-746-01
 30 x 8 mm
 Fig. 1 - 2

RH-746-02
 40 x 7 mm
 Fig. 3 - 4

RH-746-03
 40 x 7 mm
 Fig. 5 - 6

RH-746-04
 40 x 7 mm
 Fig. 7 - 8

RH-746-05
 40 x 7 mm
 Fig. 9 - 10

Alar Hooklets
 Nasenflügelhäkchen
 Ganchitos para alas de nariz
 Crochets pour les ailes du nez
 Uncini per le pinne nasali



RH-761-22
Lange Alar Hook
 Button End
 Overall Length 23 cm (9")



RH-763-17
Fomon Alar Hook
 Button End
 Overall Length 17 cm (6 3/4")



RH-765-16
Fomon Alar Hook
 Button End
 Overall Length 16 cm (6 1/4")



RH-767-20
Cottle-Neivert Alar Hook
 Overall Length 19.5 cm (7 3/4")

Alar Hooklets, Nasal Retractors
 Nasenflügelhäkchen, Nasenhaken
 Ganchitos para alas de nariz, Separadores nasales
 Crochets pour les ailes du nez, Écarteurs nasaux
 Uncini per le pinne nasali, Divaricatori nasali

RH
RHINOLOGY INSTRUMENTS


Cottle Alar Hook
 Overall Length 14.5 cm (5 3/4")



RH-771-01
 Sharp




RH-771-02
 Probe on
 Left Prong




RH-771-03
 Probe on
 Right Prong




RH-775-00
Cottle-Millard Finger Ring Retractor
 5 x 12 mm
 Overall Length 5.5 cm (2 1/4")




Converse Nasal Retractor
 Overall Length 9 cm (3 1/2")




RH-781-01
 18 x 9.5 mm



RH-781-02
 28 x 12.5 mm



RH-781-03
 41 x 13.5 mm



RH-781-04
 53 x 13.5 mm

Nasal Retractors
 Nasenhaken
 Separadores nasales
 Ecarteurs nasaux
 Divinatori nasal



RH-785-19
Aufricht Nasal Retractor
 40 x 12 mm
 Overall Length 19 cm (7 1/2")



RH-787-11
Converse Nasal Retractor
 44 x 11 mm
 Overall Length 10.5 cm (4")

Nasal Retractors
 Nasenhaken
 Separadores nasales
 Ecarteurs nasaux
 Divaricatori nasalii

RH
RHINOLOGY INSTRUMENTS


RH-789-16
Neivert Nasal Retractor
 Overall Length 17 cm (6 $\frac{3}{4}$ ")

RH-791-15
Cottle Alar Retractor
 Overall Length 16 cm (6 $\frac{1}{4}$ ")

RH-793-14
Aufrecht Nasal Retractor
 Blunt, 60 x 6.5 mm
 Overall Length 12 cm (4 $\frac{3}{4}$ ")

RH-801-20
Kressner Septum Forceps
 Overall Length 21 cm (8 $\frac{1}{4}$ ")

Septum Straightening Forceps

Septumrichtzangen

Pinzas para enderezar el tabique nasal

Pince pour le redressement de la cloison nasale

Pinze per raddrizzare il setto nasale


Cottle-Walsham Septum Straightening Forcep

Overall Length 23 cm (9")

RH-804-01
Straight

RH-804-02
Left

RH-804-03
Right

RH-807-23
Asch Septum Straightening Forcep

Overall Length 22.5 cm (8 3/4")


RH-811-15
Aiach Columella Suturing Forcep

5.0 mm

Overall Length 15.5 cm (6")

Septum Straightening Forceps

Septumrichtzangen

Pinzas para enderezar el tabique nasal

Pince pour le redressement de la cloison nasale


Pinze per raddrizzare il setto nasale

RH
RHINOLOGY INSTRUMENTS

RH-814-11
Ogura Skin Grasping Forcep

4.0 mm

Overall Length 11 cm (4 1/4")


RH-816-12
Derlacki Ossicle Holding Clamp

1.5 x 5.0 mm Vertical and

1.5 x 2.0 Horizontal Grooves

Overall Length 12.5 cm (5")


RH-817-14
Cottle Skin Grasping Forcep

1.7 mm

Overall Length 14 cm (5 1/2")



Antrum Curettes, Adenoid Curettes
 Antrumkürettten, Ringmesser
 Curetas para el antro, Cuchillos para adenoides
 Curettes pour l'antre, Couteaux à végétations adenoïdales
 Curette per antro, Coltelli per vegetazioni adenoidee



RH-820-11
Cottle Columella Forcep
 12.0 mm
 Overall Length 11 cm (4 1/4")



RH-822-11
Hildmann Cartilage Forcep
 Overall Length 11 cm (4 1/4")



RH-839-21
Faulkner Antrum Curette
 Overall Length 21 cm (8 1/4")

Adenoid Curettes
 Ringmesser
 Cuchillos para adenoides
 Couteaux à végétations adenoïdales
 Cotelli per vegetazioni adenoidee

RH
RHINOLOGY INSTRUMENTS


Coakley Antrum Curette
 Overall Length 17 cm (6 $\frac{3}{4}$ ")



Kuhn-Bolger Frontal Sinus Curette
 Overall Length 19 cm (7 $\frac{1}{2}$ ")

Adenoid Curettes

Ringmesser

Cuchillos para adenoides

Couteaux à végétations adenoïdales

Coltelli per vegetazioni adenoidee



Beckmann Adenoid Curette
Overall Length 22.5 cm (8 3/4")

RH-851-01 13.0 mm Fig. 1	RH-851-02 15.0 mm Fig. 2	RH-851-03 17.0 mm Fig. 3	RH-853-01 13.0 mm Fig. 1	RH-853-02 15.0 mm Fig. 2	RH-853-03 17.0 mm Fig. 3
RH-851-04 19.0 mm Fig. 4	RH-851-05 21.0 mm Fig. 5	RH-851-06 23.0 mm Fig. 6	RH-853-04 19.0 mm Fig. 4	RH-853-05 21.0 mm Fig. 5	RH-853-06 23.0 mm Fig. 6



Beckmann Adenoid Curette
Angled Shaft
Overall Length 22 cm (8 3/4")

Adenoid Curettes
 Ringmesser
 Cuchillos para adenoides
 Couteaux à végétations adenoïdales
 Cotelli per vegetazioni adenoidee

RH
RHINOLOGY INSTRUMENTS

RH-855-01

 11.0 mm
 Fig. 1

RH-855-02

 13.0 mm
 Fig. 2

RH-855-03

 15.0 mm
 Fig. 3

RH-855-04

 17.0 mm
 Fig. 4

RH-855-05

 19.0 mm
 Fig. 5

Barnhill Adenoid Curette
 Overall Length 22 cm (8 ¾")


RH-857-01

 10.0 mm
 Fig. 1

RH-857-02

 12.0 mm
 Fig. 2

RH-857-03

 14.0 mm
 Fig. 3

RH-857-04

 16.0 mm
 Fig. 4

St. Clair Thompson Adenoid Curette
 Overall Length 21 cm (8 ¼")

Adenoid Curettes




Ringmesser

Cuchillos para adenoides

Couteaux à végétations adenoïdales

Coltelli per vegetazioni adenoidee


LaForce Adenotome

 With Non-Perforated Blade
 Overall Length 25 cm (9 3/4")

RH-861-01
 Fig. 1

RH-861-02
 Fig. 2

RH-861-03
 Fig. 3

RH-861-11

Small

RH-861-12

Medium

RH-861-13


Large

LaForce Blade
 Non-Perforated


Adenoid Curettes
 Ringmesser
 Cuchillos para adenoides
 Couteaux à végétations adenoïdales
 Cotelli per vegetazioni adenoidee

RH
RHINOLOGY INSTRUMENTS



LaForce Adenotome
 With Non-Perforated Blade
 Overall Length 25 cm (9 ¾")



RH-863-01
 Fig. 1



RH-863-02
 Fig. 2



RH-863-03
 Fig. 3




LaForce Blade
 Perforated

RH-863-11
 Small

RH-863-12
 Medium

RH-863-13
 Large

Frontal Sinus Dilators
 Sinus-Dilatatoren
 Dilatadores para el seno
 Dilatateurs pour le sinus frontal
 Dilatatori per il seno frontale




Eicken-Killian Irrigation Cannula
 Overall Length 14 cm (5 ½")

- RH-881-20**
 \varnothing 2.0 mm
- RH-881-25**
 \varnothing 2.5 mm
- RH-881-30**
 \varnothing 3.0 mm
- RH-881-40**
 \varnothing 4.0 mm



Ritter-Halle Frontal Sinus Dilator
 Overall Length 14.5 cm (5 ¾")

- RH-891-25**
 \varnothing 2.5 mm
- RH-891-30**
 \varnothing 3.0 mm
- RH-891-40**
 \varnothing 4.0 mm



SURTEX
INSTRUMENTS

Crafted to Excellence

TL

TONSILECTOMY, LARYNGO-BRONCHOSCOPY

Tonsillektomie, Laryngo-Bronchoskopie

Tonsilectomía, Laringo-Bronchoscopia

Tonsillectomie, Laryngo-Bronchoscopie

Tonsillectomia, Laringo-Bronchoscopia

Tonsil Dissector / Retractors
 Tonsillenelevator / Retraktor
 Dissector / Separador para amígdalas
 Dissecteur / Rétracteur pour amygdalas
 Disettori / Retrattori per tonsille



TL-011-21
Fischer Tonsil Knife /
Dissector
 Overall Length 21.5 cm (8 1/2")



TL-013-23
Bruenings Tonsil Knife
 Overall Length 23 cm (9")



TL-015-19
Abraham Tonsil Knife
 Overall Length 19 cm (7 1/2")



TL-021-25
Carpenter Tonsil Knife
 4.4 mm, 8.8 mm
 Overall Length 25 cm (9 3/4")



TL-027-24
Carpenter Tonsil Dissector
 Overall Length 24 cm (9 1/2")

Tonsil Dissector / Retractors
 Tonsillenelevator / Retraktor
 Dissector / Separador para amígdalas
 Dissecteur / Rétracteur pour amygdales
 Disettori / Retrattori per tonsille

TL


TL-031-22
Hurd Tonsil Dissector / Retractor
 11.0, 14.0 mm
 Overall Length 22.5 cm (8 3/4")



TL-033-21
Mollison Tonsil Dissector/Retractor
 10. mm 10 mm
 Overall Length 21 cm (8 1/4")



TL-037-23
 10.0, 5.0 mm

TL-039-23
 16.0, 5.0 mm

Tonsil Elevator / Retractors
 Tonsil Aufzüge / Wundhaken
 Amígdalas Elevadores / Separadores
 Amygdale Ascenseurs / Ecarteurs
 Tonsille Ascensori / Divaricatori



TL-043-22
Gwynne-Evans Tonsil Elevator
 Overall Length 20 cm (8")



TL-045-23
Lopez-Reinke Tonsil Dissector
 With Suction
 Overall Length 23 cm (9")



Love Tonsil Retractor
 Overall Length 20 cm (8")

- TL-051-18**
18 mm
- TL-051-20**
20 mm
- TL-051-22**
22 mm

Tonsil Dissector / Retractors
 Tonsillenelevator / Retractor
 Dissector / Separador para amígdalas
 Dissecteur / Rétracteur pour amygdalas
 Disettori / Retrattori per tonsille

TL
TONSILECTOMY, LARYNGO-BRONCHOSCOPY INSTRUMENTS


TL-055-23
Nager Tonsil Retractor
 Overall Length 23 cm (9")



TL-063-17
Haslinger Tonsil Retractor
 Overall Length 17 cm (6 ¾")



TL-065-00
Roeder Tonsil Ligator
 Overall Length 26.5 cm (10 ½")



Tonsil Instruments
Tonsil Instrumente
Instrumentos de la amígdala
Instruments de amygdales
Strumenti Tonsille



Luniatschek Gauze Packer
Overall Length 17.5 cm (7")



Overall Length 20 cm (8")



White Tonsil Grasping Forcep

Tonsil Instruments
 Tonsil Instrumente
 Instrumentos de la amígdala
 Instruments de amygdales
 Strumenti Tonsille

TL

TONSILLECTOMY, LARYNGO-BRONCHOSCOPY INSTRUMENTS


TL-105-20
Biohmke Tonsil Grasping Forcep

1 x 2 Teeth

Overall Length 20 cm (8")

Schnidt Tonsil Forcep
TL-111-19

 Slightly Curved
 19 cm (7 1/2")

TL-113-19

 Strongly Curved
 18.5 cm (7 1/4")

TL-115-20
Wilson Tonsil Artery Forcep

Overall Length 20 cm (8")




Tonsil Instruments
 Tonsil Instrumente
 Instrumentos de la amígdala
 Instruments de amygdales
 Strumenti Tonsille



Colver Tonsil Grasping Forcep
 Overall Length 19 cm (7 1/2")

TL-116-19
 Straight



TL-119-22
Denis-Browne Tonsil Holding Forcep
 Curved
 Overall Length 22 cm (8 1/2")

TL-117-19
 Curved



TL-120-20
Museux Tonsil Grasping Forcep
 Overall Length 20 cm (8")

Tonsil Instruments
 Tonsil Instrumente
 Instrumentos de la amígdala
 Instruments de amygdales
 Strumenti Tonsille

TL


TL-150-28
Eves Tonsil Snare

 Without Ratchet
 Overall Length 28.5 cm (11 1/4")

TL-152-28
Eves Tonsil Snare

 With Ratchet
 Overall Length 28.5 cm (11 1/4")

TL-156-29
Bruenings Tonsil Snare

 With Double Ratchet
 Overall Length 29 cm (11 1/2")

Wire Loop for Tonsil Snare
 Pack of 100

TL-161-40
 \varnothing 40.0 mm


TL-161-50
 \varnothing 50.0 mm

TL-163-00
Snare Wire
 Pack of 100
 \varnothing 0.3 mm

Tonsil Instruments
 Tonsil Instrumente
 Instrumentos de la amígdala
 Instruments de amygdales
 Strumenti Tonsille




TL-173-00
Tyding Tonsil Snare
 With 1 Straight Tip



TL-181-00
Sluder-Ballenger Tonsillectome Complete
 With 3 Tip Ref: TL-183-01 to TL-183-03



TL-181-01
Sluder-Ballenger Tonsillectome Handle Only



TL-183-01

Fig. 1
 Size:18 mm

TL-183-02

Fig. 2
 Size:20 mm

TL-183-03

Fig. 3
 Size:23 mm

Sluder-Ballenger Tonsillectome Tip

Tonsil Instruments
 Tonsil Instrumente
 Instrumentos de la amígdala
 Instruments de amygdales
 Strumenti Tonsille

TL


TL-191-14
Corwin Tonsil Haemostat
 Overall Length 14 cm (5 1/2")



Hartmann Tonsil Punch
 Through Cutting
 Overall Length 12 cm (4 3/4")



TL-201-13
 \varnothing 13.0 mm



TL-201-15
 \varnothing 15.0 mm



TL-201-17
 \varnothing 17.0 mm



Schmeden Tonsil Punch
 Through Cutting
 Overall Length 12 cm (4 3/4")

TL-213-01
 6.5 mm


TL-213-02
 8.0 mm

TL-213-03
 9.5 mm

Tonsil Instruments
 Tonsil Instrumente
 Instrumentos de la amígdala
 Instruments de amygdales
 Strumenti Tonsille




TL-215-01
Fig. 1



TL-215-02
Fig. 2



TL-215-03
Fig. 3



JURASZ
 Adenoид Cutting Forcep
 Overall Length 21 cm (8 1/4")



TL-251-01
Fishbone Laryngeal Grasping Forcep
 Serrated Jaws Open Upwards
 Overall Length 22 cm (8 3/4")



TL-301-01
Huber Universal Handle

Tonsil Instruments
Tonsil Instrumente
Instrumentos de la amígdala
Instruments de amygdales
Strumenti Tonsille

TL

TONSILECTOMY, LARYNGO-BRONCHOSCOPY INSTRUMENTS



Bruenings Interchangeable Tube
Straight

TL-302-20
20 cm (8")

TL-302-30
30 cm (12")

TL-302-40
40 cm (16")

TL-303-20
Bruenings Interchangeable Tube
Curved
Overall Length 20 cm (8")

Microlaryngeal Instruments
 Mikrolaryngeale Instrumente
 Microlaryngeal Instrumentos
 Instruments de Microlaryngeal
 Strumenti Microlaryngeal



TL-310-01
 Fraenkel Broncho-Esophagoscopy Forcep Tip



TL-318-01
 Fraenkel Broncho-Esophagoscopy Forcep Tip



TL-323-01
 Scheinmann Broncho-Esophagoscopy Forcep Tip



TL-314-01
 Seiffert Broncho-Esophagoscopy Forcep Tip



TL-321-01
 Scheinmann Broncho-Esophagoscopy Forcep Tip



Microlaryngeal Instruments
 Mikrolaryngeale Instrumente
 Microlaryngeal Instrumentos
 Instruments de Microlaryngeal
 Strumenti Microlaryngeal

TL

TONSILLECTOMY, LARYNGO-BRONCHOSCOPY INSTRUMENTS



Bruenings Broncho-Esophagoscopy Forcep Tip



TL-326-02



TL-326-03



Bruenings Broncho-Esophagoscopy Forcep Tip



TL-328-02



TL-328-03



TL-334-01

Lange Broncho-Esophagoscopy Forcep Tip



Krause Broncho-Esophagoscopy Forcep Tip



TL-338-01



TL-338-02



TL-342-01

Krause Broncho-Esophagoscopy Forcep Tip

Microlaryngeal Instruments
 Mikrolaryngeale Instrumente
 Microlaryngeal Instrumentos
 Instruments de Microlaryngeal
 Strumenti Microlaryngeal



Cordes Broncho-Esophagoscopy Forcep Tip



TL-347-01



TL-347-02



TL-347-03



TL-347-04



TL-349-01

Schumacher Broncho-Esophagoscopy Forcep Tip



Cordes Broncho-Esophagoscopy Forcep Tip



TL-353-01



TL-353-02



TL-357-01

Cordes Broncho-Esophagoscopy Forcep Tip




TL-359-01

Rosenberg Broncho-Esophagoscopy Forcep Tip

Microlaryngeal Instruments
 Mikrolaryngeale Instrumente
 Microlaryngeal Instrumentos
 Instruments de Microlaryngeal
 Strumenti Microlaryngeal

TL

TONSILLECTOMY, LARYNGO-BRONCHOSCOPY INSTRUMENTS


*** Micro Laryngeal Scissors**

TL-361-01
 Straight, 25 cm

TL-361-02
 Straight Horizontal 25 cm

TL-361-03
 Angled Up, 25 cm

TL-361-04
 Angled Up Horizontal 25 cm

TL-361-05
 Curved Right, 25 cm

TL-361-06
 Curved Left, 25 cm

*** Micro Laryngeal Forceps**

TL-363-01
 Straight, 25 cm

TL-363-02
 Straight Horizontal, 25 cm

TL-363-03
 Angled Up Horizontal, 25 cm

TL-363-04
 Curved Right, 25 cm

TL-363-05
 Curved Left, 25 cm

TL-363-06
 Curved Vertical, 25 cm

*** Micro Laryngeal Forceps**

Cup Shaped, Ø 2.5 mm


TL-365-01
 Straight, 25 cm

TL-365-02
 Angled Up, 25 cm

TL-365-03
 Angled Up Horizontal, 25 cm

TL-365-04
 Curved Right, 25 cm

TL-365-05
 Curved Left 25 cm

TL-365-06
 Oval Shaped, 25 cm

Microlaryngeal Instruments
 Mikrolaryngeale Instrumente
 Microlaryngeal Instrumentos
 Instruments de Microlaryngeal
 Strumenti Microlaryngeal



* **Micro Laryngeal Forceps**

For Children
 Shaft Length 35 cm (13 3/4")



TL-367-01



TL-367-02



TL-367-03



TL-367-04
 \varnothing 3.0 mm



TL-367-05
 \varnothing 3.0 mm



TL-367-06
 \varnothing 3.0 mm



TL-367-07
 \varnothing 3.0 mm



TL-367-08
 \varnothing 3.0 mm



TL-367-09
 2×5 mm



TL-367-10
 3×5 mm



TL-367-11

* **Micro Laryngeal Forceps**

For Adults



TL-369-01
 Shaft Length 40 cm (15 3/4")



TL-369-02
 Shaft Length 40 cm (15 3/4")



TL-369-03
 Shaft Length 40 cm (15 3/4")



TL-369-04
 Shaft Length 40 cm (15 3/4")

TL-371-01

Shaft Length 50 cm (20")

TL-371-02

Shaft Length 50 cm (20")

TL-371-03

Shaft Length 50 cm (20")

TL-371-04

Shaft Length 50 cm (20")

TL-373-01

Shaft Length 60 cm (23")

TL-373-03

Shaft Length 60 cm (23")

TL-373-03

Shaft Length 60 cm (23")

TL-373-04

Shaft Length 60 cm (23")

Microlaryngeal Instruments
 Mikrolaryngeale Instrumente
 Microlaryngeal Instrumentos
 Instruments de Microlaryngeal
 Strumenti Microlaryngeal

TL

* **Micro Laryngeal Forceps**
 For Adults



TL-369-05

Shaft Length 40 cm (15 ¾")



TL-369-06

Shaft Length 40 cm (15 ¾")



TL-369-07

Shaft Length 40 cm (15 ¾")



TL-369-08

Shaft Length 40 cm (15 ¾")

TL-371-05

Shaft Length 50 cm (20")

TL-371-06

Shaft Length 50 cm (20")

TL-371-07

Shaft Length 50 cm (20")

TL-371-08

Shaft Length 50 cm (20")

TL-373-05

Shaft Length 60 cm (23")

TL-373-06

Shaft Length 60 cm (23")

TL-373-07

Shaft Length 60 cm (23")

TL-373-08

Shaft Length 60 cm (23")



TL-369-09

Shaft Length 40 cm (15 ¾")



TL-369-10

Shaft Length 40 cm (15 ¾")



TL-369-11

Shaft Length 40 cm (15 ¾")



TL-369-12

Shaft Length 40 cm (15 ¾")

TL-371-09

Shaft Length 50 cm (20")

TL-371-09

Shaft Length 50 cm (20")

TL-371-11

Shaft Length 50 cm (20")

TL-371-12

Shaft Length 50 cm (20")

TL-373-05

Shaft Length 60 cm (23")

TL-373-06

Shaft Length 60 cm (23")

TL-373-07

Shaft Length 60 cm (23")

TL-373-08

Shaft Length 60 cm (23")



TL-369-13

Shaft Length 40 cm (15 ¾")

TL-371-13

Shaft Length 50 cm (20")

TL-373-13

Shaft Length 60 cm (23")

Tonsillectomy Instruments
 Tonsil Instrumente
 Instrumentos de la amígdala
 Instruments de amygdales
 Strumenti Tonsille


TL-381-30

30 cm (12")

TL-381-40

 40 cm (15 $\frac{3}{4}$ ")

TL-381-50

 50 cm (19 $\frac{1}{2}$)

TL-381-60

 60 cm (23 $\frac{1}{2}$)

Chevalier-Jackson Square Punch

TL-383-30

30 cm (12")

TL-383-40

 40 cm (15 $\frac{3}{4}$ ")

TL-383-50


 50 cm (19 $\frac{1}{2}$)

TL-383-60

 60 cm (23 $\frac{1}{2}$)

Chevalier-Jackson Grasping Forcep


Tonsillectomy Instruments
 Tonsil Instrumente
 Instrumentos de la amígdala
 Instruments de amygdales
 Strumenti Tonsille

TL

TL-385-30

30 cm (12")

TL-385-40

40 cm (15 ¾")

TL-385-50

50 cm (19 ½")

TL-385-60

60 cm (23 ½")

Chevalier-Jackson Grasping Forcep

TL-387-30

30 cm (12")

TL-387-40

40 cm (15 ¾")

TL-387-50

50 cm (19 ½")

TL-387-60

60 cm (23 ½")

Chevalier-Jackson Cup Shaped Biopsy Forcep
 Straight


Tonsillectomy Instruments

Tonsil Instrumente

Instrumentos de la amígdala

Instruments de amygdales

Strumenti Tonsille



Tonsillectomy Instruments
Tonsil Instrumente
Instrumentos de la amígdala
Instruments de amygdales
Strumenti Tonsille

TL

**TL-393-30**

30 cm (12")

TL-393-40

40 cm (15 ¾")

TL-393-50

50 cm (19 ½")

TL-393-60


60 cm (23 ½")

Chevalier-Jackson Cup Shaped Biopsy Forcep

Curved Upwards

TONSILLECTOMY, LARYNGO-BRONCHOSCOPY INSTRUMENTS

Tonsillectomy Instruments
 Tonsil Instrumente
 Instrumentos de la amígdala
 Instruments de amygdales
 Strumenti Tonsille



Laryngeal Mirror
 Rhodium Plated With Fixed Handle
 Overall Length 18.5 cm (7 1/4")

TL-402-08

Fig. 00
 Ø 8 mm

TL-402-20

Fig. 5
 Ø 20 mm

TL-402-10

Fig. 0
 Ø 10 mm

TL-402-22

Fig. 6
 Ø 22 mm

TL-402-12

Fig. 1
 Ø 12 mm

TL-402-24

Fig. 7
 Ø 24 mm

TL-402-14

Fig. 2
 Ø 14 mm

TL-402-26

Fig. 8
 Ø 26 mm

TL-402-16

Fig. 3
 Ø 16 mm

TL-402-28

Fig. 9
 Ø 28 mm

TL-402-18

Fig. 4
 Ø 18 mm

TL-402-30

Fig. 10
 Ø 30 mm



Laryngeal Mirror
 Rhodium Plated Without Handle
 Overall Length 13.5 cm (5 1/4")

TL-406-08

Fig. 00
 Ø 8 mm

TL-406-20

Fig. 5
 Ø 20 mm

TL-406-10

Fig. 0
 Ø 10 mm

TL-406-22

Fig. 6
 Ø 22 mm

TL-406-12

Fig. 1
 Ø 12 mm

TL-406-24

Fig. 7
 Ø 24 mm

TL-406-14

Fig. 2
 Ø 14 mm

TL-406-26

Fig. 8
 Ø 26 mm

TL-406-16

Fig. 3
 Ø 16 mm

TL-406-28

Fig. 9
 Ø 28 mm

TL-406-18

Fig. 4
 Ø 18 mm

TL-406-30

Fig. 10
 Ø 30 mm

Tonsillectomy Instruments
 Tonsil Instrumente
 Instrumentos de la amígdala
 Instruments de amygdales
 Strumenti Tonsille

TL

TONSILLECTOMY, LARYNGO-BRONCHOSCOPY INSTRUMENTS



Michel Laryngeal Mirror
 With Handle , Moveable Mirror
 Overall Length 25 cm (9 3/4")

TL-413-12

Fig. 1
 Ø 12 mm

TL-413-14

Fig. 2
 Ø 14 mm

TL-413-16

Fig. 3
 Ø 16 mm

TL-413-18

Fig. 4
 Ø 18 mm



TL-420-10
Laryngeal Mirror Handle
 Overall Length 10 cm (4")



SURTEX
INSTRUMENTS

Crafted to Excellence

TM

TRACHEOSTOMY INSTRUMENTS

Tracheotomie

Traqueotomía

Trachéotomie

Tracheotomia



Tracheal Dilators
Trachealdilatatoren
Dilatadores traqueales
Dilatateurs pour la trachéotomie
Dilatatori per tracheotomia



TM-101-75

\varnothing 7.5 mm

TM-101-90

\varnothing 9.0 mm

Uckermann Tracheotomy Trocar
Overall Length 11.5 cm (4 1/2")



TM-151-14

Laborde Tracheal Dilator

Overall Length 14 cm (5 1/2")

Tracheal Dilators, Tracheal Hooks
 Trachealdilatatoren, Tracheahaken
 Dilatadores traqueales, Ganchos traqueales
 Dilatateurs pour la trachéotomie, Ecarteurs pour trachéotomie
 Dilatatori per tracheotomia, Uncini tracheali

TM

TRACHEOSTOMY INSTRUMENTS



TM-155-14
Trousseau Tracheal Dilator
 Overall Length 14.5 cm (5 $\frac{3}{4}$ ')



TM-161-07
Bose Tracheal Retractor
 Overall Length 6 cm (2 $\frac{1}{4}$ ')




TM-165-01
 Left

TM-165-02
 Right

Bose Tracheal Retractor
 Overall Length 16 cm (6 $\frac{1}{4}$ ')

Tracheal Hooks
 Tracheahaken
 Ganchos traqueales
 Ecarteurs pour trachéotomie
 Uncini tracheali



TM-171-16

Sharp



TM-181-21

Schoeborn Tracheal Retractor
 Overall Length 21 cm (8 1/4")



TM-183-13

Rose Tracheal Retractor
 Overall Length 13.5 cm (5 1/4")



TM-185-11

Luer Tracheal Retractor
 Overall Length 11 cm (4 1/2")

Goitre Enucleator
 Kropfsonden
 Sondas para bocio
 Spatules à goître
 Uncini tracheali

TM

TRACHEOSTOMY INSTRUMENTS



TM-189-17
Lukens Thymus Retractor
 Blunt
 Overall Length 17.5 cm (7")



TM-191-17
Jackson Tracheal Hook
 Overall Length 18 cm (7")




Kocher Goiter Enucleator
 Overall Length 16.5 cm (6 1/2")

TM-201-15
 Without Eye

TM-203-15
 With Eye



Goitre Grasping Forceps, Tracheal Tubes
 Kropffasszangen, Trachealkanülen
 Pinzas para coger el bocio, Cánulas traqueales
 Pinces à goître, Canules tracheales
 Pinze per gozzo, Cannule tracheali



Lahey Goiter Grasping Forcep

TM-250-15

2 x 2 Teeth
 15.5 cm (6")

TM-252-15

3 x 3 Teeth
 15.5 cm (6")

TM-252-20

3 x 3 Teeth
 20 cm (8")



Luer Tracheal Tube

TM-301-04

45 mm
 Ø 4.0
 Fig. 0

TM-301-05

40 mm
 Ø 5.0
 Fig. 1

TM-301-06

55 mm
 Ø 6.0
 Fig. 2

TM-301-07

60 mm
 Ø 7.0
 Fig. 3

TM-301-08

65 mm
 Ø 8.0
 Fig. 4

TM-301-09

70 mm
 Ø 9.0
 Fig. 5

TM-301-10

75 mm
 Ø 10.0
 Fig. 6

TM-301-11

80 mm
 Ø 11.0
 Fig. 7

TM-301-12

85 mm
 Ø 12.0
 Fig. 8

TM-301-13

90 mm
 Ø 13.0
 Fig. 9

Tracheal Tubes
 Trachealkanülen
 Cánulas traqueales
 Canules tracheales
 Cannule tracheali

TM

TRACHEOSTOMY INSTRUMENTS




Chevalier-Jackson Tracheal Tube

TM-321-04	TM-321-05	TM-321-06	TM-321-07
45 mm	40 mm	55 mm	60 mm
Ø 4.0	Ø 5.0	Ø 6.0	Ø 7.0
Fig. 0	Fig. 1	Fig. 2	Fig. 3

TM-321-08	TM-321-09	TM-321-10	TM-321-11
65 mm	70 mm	75 mm	80 mm
Ø 8.0	Ø 9.0	Ø 10.0	Ø 11.0
Fig. 4	Fig. 5	Fig. 6	Fig. 7

TM-321-12	TM-321-13
85 mm	90 mm
Ø 12.0	Ø 13.0
Fig. 8	Fig. 9



Speaking Tube

TM-351-08	TM-351-09	TM-351-10	TM-351-11
65 mm	70 mm	75 mm	80 mm
Fig. 4	Fig. 5	Fig. 6	Fig. 7

TM-351-12	TM-351-13
85 mm	90 mm
Fig. 8	Fig. 9

General Terms and Conditions of Sales

1. Definitions

“Company” means Surtex Instruments Ltd (registered in the United Kingdom under number 7035724).

“Contract” means the Order accepted by the Company and incorporating these conditions.

“Customer” means the company, firm or organisation who’s Order for the Goods is accepted by the Company.

“Goods” means the goods which the Company agrees to supply in accordance with these conditions and which are identified on the Order. “Order” means the order for the purchase of the Goods submitted by the Customer.

“Goods” means the goods agreed to be provided to the customs by Surtex Instrument Ltd.

“Quotation” means the quotation by Surtex Instruments stating the price at which and the time period during which it will supply the goods to the customer.

1. Any reference to a law is a reference to it as in force for the time being taking account of any amendment, extension, application or re-enactment and includes any subordinate legislation for the time being in force made under it.

2. Where a condition provides for a list of items following the word “including” or “includes”, then such list is not to be interpreted as being an exhaustive list.

2. General

2.1 Headings are purely for the ease of reference and do not form part of effect the interpretation of Conditions.

2.2 In the event of a conflict the contractual documents shall be applied in the following order:

a) the Quotation

b) the Conditions

2.3 No other terms and conditions or warrantied are to be incorporated into the Contract unless agreed expressly in writing with Surtex Instruments. In particular, the Customer’s terms and conditions are expressly excluded whether provided in full, referred to in other documents and whether before or after the Contract is formed.

2.4 Surtex Instruments may update or amend these Conditions from time to time to comply with the law or to meet the changing business requirements without notice to the customer.

3. Ordering

3.1 Each Order (whether submitted by telephone, fax, post, email or electronic ordering system) shall be deemed to be an offer by the Customer to purchase the Goods pursuant to these conditions. The Company may, at its complete discretion, refuse to accept any Order. The Contract will be formed upon confirmation of acceptance of the Customer’s purchase order by Surtex Instruments.

3.2 The Customer shall be responsible to the Company for ensuring the accuracy of the terms of any Order (including the description of the Goods ordered).

Samples and Free of Charge Goods

3.3 For the avoidance of doubt, where no order is placed and the Company provides Samples, to familiarize the Customer with the goods, as part of a sampling program or where the Company provides Free of Charge Goods to the Customer as part of a clinical trial or compassionate use program these conditions shall, to the extent applicable, apply to such Samples and Free of Charge Goods.

4. Prices

4.1 The price is exclusive of VAT and all other applicable taxes and duties, which the Customer shall be additionally liable to pay the Surtex Instruments.

4.2 All prices quoted are subject to change up to the date of delivery, subject to prior written notice by Surtex Instruments.

4.3 Surtex Instruments reserves the right, by giving notice to the Customer at any time before delivery, to increase the price of the Goods to reflect any increase in the cost of Surtex Instruments which is due to any factor beyond the reasonable control of Surtex Instruments.

4.4 Errors on Surtex Instruments website prices and availability of products and services are subject to change without notice errors will be corrected where discovered, and Surtex Instruments reserves the right to revoke any stated offer and to correct any errors, inaccuracies or omissions including after an order has been submitted and whether or not the order has been confirmed and your credit card charged. Surtex Instruments shall have the right to refuse or cancel any orders placed for product listed at the incorrect price. Surtex Instruments shall have the right to refuse or cancel any such orders whether or not the order has been confirmed and your credit card charged. If your credit card has already been charged for the purchase and your order is cancelled, Surtex Instruments will issue a credit to your credit card account in the amount of the charge. Individual bank policies will dictate when this amount is credited to your account. If you are not fully satisfied with your purchase, you may return it in accordance with Surtex Instruments’ Return Policy.

5. Delivery

5.1 Delivery of the Goods shall be made to the delivery address specified as stated on the Pro Forma Invoice or communicated separately whichever is greater after the date of the Surtex Instruments acceptance of the Order.

5.2 Surtex Instruments reserves the right to withhold delivery of any Goods if any sum due under any contract with the Customer is overdue or if, in the reasonable opinion of Surtex Instruments, the financial standing of the Customer has been impaired for any other reason.

5.3 Time for delivery shall not be of the essence and, Surtex Instruments shall not be liable for any loss or damage whatever due to any failure by the Surtex Instruments to deliver the Goods (or any of them) promptly or at all. The customer will have no right damages or to cancel the order for any failure to meet the delivery times stated.

5.4 The delivery of the Goods may happen by separate instalments. Each separate instalment shall be invoiced and paid for in accordance with the provisions of the Contract. Each instalment shall be a separate Contract and no cancellation or termination of any one Contract relating to an instalment shall entitle the Customer to repudiate or cancel any other Contractor instalment.

5.5 The Customer shall take delivery of the Goods promptly whenever they are tendered for delivery (which shall include signing the relevant delivery note) If, for any reason, the Customer fails to take such delivery of the Goods or otherwise causes or requests a delay in delivery:

a) The Goods shall be deemed to have been delivered;

b) Risk in the Goods shall pass to the Customer, and

c) Surtex Instruments shall be entitled to store or arrange for the storage of the Goods until actual delivery and charge the Customer for the reasonable costs (including insurance) of such storage and any additional charges incurred in delivering the Goods at a later date.

5.6 The Customer shall ensure that adequate and safe facilities and procedures exist at the Customer’s premises for the delivery of the Goods.

5.7 The Customer shall immediately notify Surtex Instruments Customer Support Centre if the Goods are not received within 10 days of the specified delivery date.

5.8 The customer outside of the UK shall be solely responsible for obtaining all import authorisations.

6. Title and Risk

6.1 The Goods shall be at the Customer’s risk at the time of delivery unless the parties have agreed in writing that the Goods are to be collected from the Surtex Instruments premises, in which case, the Goods shall be at the Customer’s risk when Surtex Instruments notifies the Customer that the Goods are available for collection.

1. Property in the Goods shall not pass from Surtex Instruments until:

a) The Customer has paid the price in full; and

b) No other sums whatsoever are due from the Customer to Surtex Instruments.

2. Until property in the Goods passes to the Customer in accordance with the condition 6.1, the Customer shall:
 - a) Hold the Goods on a fiduciary basis as bailed for Surtex Instruments
 - b) Store the Goods separately from all other goods in its possession and mark them in such a way that they are clearly identified as Surtex Instrument's property; and
 - c) Maintain the Goods in satisfactory condition and keep them insured on Surtex Instrument's behalf for their full price.
3. The Customer may sell at full market value or use the Goods only in the ordinary course of the Customer's business. Any such sale or dealing shall be a sale or use of Surtex Instrument's property on the Customer's own behalf and the Customer shall deal as principal when making such sales or dealings. Until property in the Goods passes from Surtex Instruments, the proceeds of sale or otherwise of the Goods, to the full extent of the price of the Goods, shall be held in trust for Surtex Instruments and shall not be mixed with other money or paid into any overdrawn bank account and shall be at all material times identified as Surtex Instruments money.
4. Surtex Instruments shall be entitled to recover the price notwithstanding that property in any of the Goods has not passed from Surtex Instruments.
5. Until such time as the property in the Goods passes from Surtex Instruments, if the Customer: a) Enters into liquidation or has a winding-up order made against it or has an administrator, receiver or manager appointed in respect of its assets; or b) Fails to make payment in accordance with these conditions, the Customer shall upon request deliver up Surtex Instruments such of the Goods as have not ceased to be in existence or resold. If the Customer fails to do so, Surtex Instruments may enter upon any premises owned. Occupied or controlled by the Customer where the Goods are situated and repossess the Goods.
6. The Customer shall not pledge or in any way charge by way of security for any indebtedness, any of the Goods which are the property Surtex Instruments. If the Customer does so, all sums whatever owing by the Customer to Surtex Instruments shall forthwith become due and payable.
7. On termination of the Contract, howsoever caused, Surtex Instruments rights contained in this condition 6 shall remain in effect.

7. Terms of Payment

- 7.1 Payment is due to Surtex Instruments within 30 days of invoice (the "Due Date"). In the case of payment other than by cash, payment shall only be deemed to have been made on receipt of cleared funds.
- 7.2 Each consignment or partial delivery shall be invoiced and paid separately.
- 7.3 No disputes arising under the contract no delays in the delivery period beyond the control of Surtex Instruments shall interfere with prompt payment in full by Customer. The Customer shall not in any circumstances or for any reason whatsoever be entitled to make any deduction or withhold any sum from the price by way of set-off or otherwise.
- 7.4 If the Customer fails to make the payment on the due date, then Surtex Instruments without prejudice to other rights and remedies shall be entitled to charge the Customer interest (both before and after judgment) on the amount unpaid at the rate of 4% per annum above the basic rate of HSBC plc during the delay and suspend all further deliveries without notice.
- 7.5 All payments payable to the Company under the Contract shall become due immediately on its termination, despite any other provision.

8. Warranties

- 8.1 Surtex Instruments will use reasonable endeavours to ensure that the Goods
 - a) Correspond to their description and any applicable specification;
 - b) Be of satisfactory quality and fit for purpose held out by Surtex Instruments or made to known to the customer by Surtex Instruments
 - c) Where applicable, be free from defects in design, materials and workmanship and remain so for a period of 12 months after delivery; and
 - d) Comply with all applicable statutory and regulatory requirements relating to the manufacture, labelling, packaging, storage, handling and delivery of Goods
- 8.3 In case of Goods not manufactured by Surtex Instruments
 - a) Surtex Instruments gives no assurance or warranty whatsoever that the sale or use of the Goods will not infringe the patent design right registered, design copyright or other industrial property rights of other person firm or Company and
 - b) The warranty given to the Customer shall be limited to the warranty (if any) which Surtex Instruments receives from the manufacturer or supplier of the Goods to Surtex Instruments.

9. Missing and Defective Goods

- 9.1 The Customer must inspect the Goods immediately on arrival at their premises.
- 9.2 In case of missing goods or defects apparent on inspection:
 - a) a written complaint must be made to Surtex Instruments by the customer within seven days of receipt of the Goods specifying the shortage and/or defect and
 - b) Surtex Instruments shall be given an opportunity to inspect the Goods and investigate any complaint before any use is made of the Goods.
- 9.3 In the case of defects not apparent on inspection:
 - a) The Customer shall send a written complaint to Surtex Instruments as soon as reasonably practicable after the defect is noticed and ensure no use is made of or alterations made to the Goods thereafter;
 - b) The Customer shall allow Surtex Instruments 30 days to inspect the Goods and shall take such steps as are necessary to enable Surtex Instruments to do so including returning of the Goods to its head office in the United Kingdom if requested.
- 9.4 In any potential case of a product being defective pursuant to Part 1 of the Consumer Protection Act 1987, the customer must send the complaint within 12 months of the date of delivery by Surtex Instruments or in the case of items not manufactured by Surtex Instruments, within the guarantee period specified by the manufacturer of such item.
- 9.5 Any complaints made outside of the periods specified in clause 9 will not be considered.
- 9.6 All allegations of total non-delivery of any consignment of the Goods must be made by notice in writing by the customer to Surtex Instruments within 10 days of the specified delivery date. Failure by the Customer to give such notice shall render the Customer liable for any temporary or permanent loss of the Goods and all additional costs and expenses incurred by Surtex Instruments in relation to such loss.

10. Liability

- 10.1 Save where Goods manufactured by Surtex Instruments are held to be defective as defined in Part 1 of the Consumer Protection Act 1987, Surtex Instruments sole obligation in the event of any missing or defective Goods, shall be to repair or replace the Goods.
- 10.2 Surtex Instruments shall not be liable in respect of consequential loss, loss of profits, damage to the property, loss of goodwill, loss of business opportunity, wasted expenditure, cost of mitigation arising out of or occasioned by any fault or defect in the Goods supplied or any other losses howsoever caused.
- 10.3 The maximum total liability under or arising out of in connection with this Agreement will not exceed the total of the Contract value.
- 10.4 Nothing in this Agreement shall exclude or limit either party's liability for death or personal injury resulting from negligence or their fraudulent acts, omissions or statements or for infringement of third party's intellectual property rights.
- 10.5 For the purpose of the Health and Safety at Work Act 1974, Surtex Instruments relies upon the testing examination and research carried out by the suppliers to it of propriety equipment and component parts of Goods.
- 10.6 This clause 10 shall survive the termination of the contract.

11. Return of Goods

- 11.1 In cases where Goods are not defective and the customer requests to return the Goods, Surtex Instruments shall be under no contractual obligation to accept the return except in the event of any error on its parts as to the amount or type of Goods delivered.

- 11.2 Any requests to return must be made to Surtex Customer Support Centre within ten working days of the date of actual delivery.
- 11.3 The return of Goods to correct an error in delivery or ordering will only be considered if the Customer Care Team is notified immediately or in any event within 48 hours of delivery to the Customer. Goods must be returned to Surtex Instruments within 5 days of delivery to the customer to correct an error in delivery or ordering.
- 11.4 For the avoidance of doubt, unless defective, the returns of Goods as detailed in sections 11.2, 11.4, 11.5 and 11.6 must be unused, undamaged and in the original packaging. The Goods must be the same batch as originally supplied by Surtex Instruments and are not falsified/counterfeit. All goods are individually assessed by Surtex Instruments before being authorised to return.
- 11.5 Return of Goods for any reason shall require prior consent from Surtex Instruments which shall be sent solely at the discretion of Surtex Instruments. To obtain consent, the Customer Care Team must be notified. Consideration to requests for a return may be given if:
- a) of a bona fide nature;
 - b) made within 10 working days of the actual delivery; and
 - c) unless defective, the Goods must be unused, undamaged and in the original undamaged packaging. The Goods must be the same batch as originally supplied by Surtex Instruments and not falsified/counterfeit.
- 11.6 On the customer contacting Surtex Instruments to request a return, the Customer will be issued with a returns merchandise authorisation number by Surtex Instruments. Goods returned without the prior written approval of Surtex Instruments may be returned to the customer or disposed of at the absolute discretion of the company.
- 11.7 Goods returned will be subjected to the following charges:
- a) for palletised items there will be a standard charge of £75.00 per pallet for any return uplifted by Surtex Instruments.
 - b) for returns of small parcel there will be a minimum £15.00 charge per parcel and a 30% re-stocking charge.
 - c) the event of return resulting from an error on the part of Surtex Instruments in relation to the amount or type of Goods delivered, no charges under this clause 11.7 will apply.
- 11.8 In the event, that any Goods are subject to be called by Surtex Instruments for the purpose of recall, the customer shall provide to Surtex Instruments as much information as it is reasonably possible regarding the whereabouts of the Goods.

12. Cancellation Policy

- 12.1 In the event of cancellation before delivery and if returned after delivery where products are not defective, there will be a re-stocking charge of not less than 30% of the invoice value charged by Surtex Instruments.
- 12.2 Cancellation of an order will only be agreed by Surtex Instruments at its sole discretion, on condition that all costs and expenses incurred by Surtex Instruments up to the time of cancellation and loss of profits and other loss or damage resulting to Surtex Instruments of such cancellation will be reimbursed by the Customer to Surtex Instruments forthwith.

13. Specifications

Surtex Instruments reserve the right to alter or change the specification and dimension of the Goods supplied with reasonable limits having regard to the nature of the Goods. These do not materially affect the quality or fitness for purpose of Goods. Dimensions specified by Surtex Instruments are to be treated as an approximate only unless the Customer specifically states in writing that exact dimensions are required. Our website and catalogues display the product images as accurately as possible, However, Surtex Instruments cannot guarantee the finish/colour you see matches the actual product as the display of the finish and/or colour depends in part, upon your monitor and its settings.

14. Marking and Instructions

- 14.1 No name, mark, numbering, colouring, appearance or log on the goods will be obscured, removed or concealed by the Customer.
- 14.2 The Customer will not repackage or alter the presentation of Goods and it will not assist, cause, or enable any other party to do any of the said acts or deal in Goods or packaging which have been subject to such acts.
- 14.3 The Customer acknowledges that it is under a duty to pass on to its customers (where appropriate) all instructions, information and warnings supplied by Surtex Instruments with the Goods.

15. Confidential Information

- 15.1 A party shall treat the Contract as private and confidential and neither party shall publish or disclose any particulars without the prior written consent of the other with the exceptions of 15.2 below.
- 15.2 The obligations expressed in 15.1 above shall not apply to any information which:
- a) Is or subsequently comes into the public domain other than by breach of this clause;
 - b) The disclosure of which is required to be disclosed by law, any governmental or regulatory authority or by a court of competent jurisdiction;
 - c) Comes into the possession of one of the parties prior to its disclosure by the other, or which is acquired lawfully and in good faith from an independent third party.
- 15.3 On the conclusion or termination of the Contract, both parties shall destroy all copies of confidential information obtained from the other and the obligation relating to the confidentiality shall continue for a period of 10 years.
- 15.4 The Customer shall not, without the prior written consent of Surtex Instruments, advertise or make known to any third that Surtex Instruments provided Goods to the Customer.

16. Insolvency

- 16.1 Without limiting its other rights or remedies, Surtex Instruments shall be entitled to cancel the contract, in whole or in part, by notice in writing if:
- a) The Customer makes a voluntary arrangement with its creditors or (being an individual firm) becomes bankrupt or (being a company) becomes subject to an administration order or goes into liquidation (otherwise than for purpose of amalgamation or reconstruction);
 - b) An encumbrancer takes possession, or a receiver is appointed, of any of the property or assets of the Customer.
 - c) The Customer ceases, or threatens to cease, to carry on business; d) Surtex Instruments reasonably apprehends that any of the events mentioned above is about to occur in relation to the Customer and notifies the Customer accordingly.

17. Sub-Contractors

- 17.1 Surtex Instruments may assign the Contractor any part of it to any company, firm or organization.
- 17.2 The Customer shall not be entitled to assign the Contractor any part of it without the prior written consent of Surtex Instruments (such consent not to be unreasonably withheld or delayed).

18. Force Majeure

- 18.1 Surtex Instruments shall not be responsible for any failure to fulfil its obligations under these conditions and/or in relation to the Goods if such failure is caused by an event which is beyond its reasonable control including acts of God, war, armed conflict, acts of terrorism, riot, strikes, industrial action (including action involving Surtex Instruments employees), fire, flood, any natural disaster, any defaults by Surtex Instruments suppliers or subcontractors, difficulties obtaining workmen or materials, and breakdown of machinery.
- 18.2 If such event continues for more than one month, Surtex Instruments may terminate the Contract without liability.

19. Waiver

19.1 waiver of any right or remedy under the contract or law is only effective if given in writing and shall not be deemed a waiver of any subsequent breach or default. No failure or delay by a party to exercise any right or remedy provided under the Contract or by law shall constitute a waiver of that or any other right or remedy. No single or partial exercise of such right or remedy shall prevent or restrict the further exercise of that or any other right or remedy.

20. Severance

20.1 any provision or part-provision of the Contract is or becomes invalid, illegal or unenforceable, it shall be deemed modified to the minimum extent necessary to make it valid, legal and enforceable whilst retaining the risk and reward between the parties. If such modification is not possible, the relevant provision or part-provision under this clause shall not affect the validity and enforceability of the rest of the Contract.

21. Variation

21.1 as set out in these Conditions, no variation of the Contract, including the introduction for any additional terms and conditions, shall become effective unless it is agreed in writing and signed by Surtex Instruments.

22. General

22.1 rights and remedies of Surtex Instruments set out in these conditions shall be in addition and without prejudice to any other rights and remedies which may be available to Surtex Instruments at common law or under the statute.

23. Governing Law and Jurisdiction

23.1 The contract and any dispute or claim arising out of or in connection with it or its subject matter or formation (including non-contractual disputes or claims) shall be governed by and construed in accordance with the law of England and Wales.

23.2 Each party irrevocably agrees that the courts of England and Wales shall have exclusive jurisdiction to settle any dispute or claim arising out of or in connection with his agreement or its subject matter or formation (including non-contractual disputes or claims).

23.3 The schedules to the Uniform Law of International Sales Act 1967 shall not in any circumstances apply to the contract.

24. Entire Agreement

24.1 This contract constitutes the whole agreement and understanding between the parties in relation to its subject matter and supersedes any previous arrangement, understanding, representations, (including any information or data set out in any sales and marketing literature), or agreement between them relating to the subject matter of this Agreement.

24.2 Each party acknowledges that, in entering into this Agreement, it has not relied on any statement, assurance or warranty (whether made negligently or innocently) other than those expressly set out in this Agreement and has no rights and remedies in respect thereof.

25. Third-Party Rights

25.1 term of this agreement shall be enforceable under the Contracts (Rights of Third Parties) Act 1999 by any legal entity who is not a party to this Agreement.

26. Notices

26.1 communications between the parties about the Contract shall be in writing and delivered by pre-paid first-class post or by hand to the party's registered office, principal place of business or such other address as may at the relevant time have been notified pursuant to this provision to the party giving notice. In the absence of evidence of earlier receipt, any notices shall be deemed to be duly served:

- a) If sent by first class post, two working days after posting; or
- b) If delivered personally, when left at the relevant address.

INSTRUCTIONS FOR REPROCESSING OF REUSABLE SURGICAL INSTRUMENTS

Please read all information contained in this insert. The use of a device for a task other than that for which it is intended, as well as improper, ineffective and insufficient maintenance can greatly reduce the life of a device and will invalidate the device's warranty. Incorrect handling and care as well as misuse can lead to premature wear or can cause hazards to patients and users.


PRODUCTS:

Scalpel Handles, knives	Raspatories	Bone Holding & Cutting Forceps
Scissors	Rhinology	Elevators
Forceps	Hooks	Vaginal Speculums
Dressing Forceps and Sponge Forceps	Curette	Bone Chisel & Gouges
Probes, cotton applicators, spatulas	Needle Holders	Obstetrics
Bone surgery	Tissue & Dissecting Forceps	

INTENDED USE:	Surtex produces the standard surgical instruments, designed for routine use in general surgical procedures. Nevertheless, it is the responsibility of the treating surgeon to select the appropriate instruments for specific procedures or operative applications. Additionally, the surgeon is accountable for ensuring that operating room personnel are adequately trained and informed, and that they possess the requisite experience to use the instruments effectively.
WARNINGS:	<ul style="list-style-type: none"> • All types of surgical instruments are supplied non-sterile. • Every instrument must undergo cleaning, disinfection, and sterilization processes prior to and following each usage. • Utilization of strong alkaline or acidic cleaning agents is discouraged.
Limitations on Reprocessing:	<ul style="list-style-type: none"> • Reusable surgical instruments are designed to withstand a maximum of 20 re-sterilization cycles. • Signs of wear or damage incurred during surgical procedures typically signify the end of a product's life.

INSTRUCTIONS	
Point of use:	<ul style="list-style-type: none"> • Use a disposable cloth or paper towel to wipe away any excess debris. • Ensure that all instruments are subjected to a high-pressure water rinse
Containment and transportation:	<ul style="list-style-type: none"> • Transporting contaminated instruments in a sealed and clean container is advisable. • It is best practice to reprocess the instrument within 2 hours after use.
Preparation for decontamination:	Personal protective equipment, including gloves, a water-resistant apron, a face shield or safety glasses, should be worn in accordance with the manufacturer's instructions for the cleaning and disinfectant agents used.
Cleaning: Manually	<p>Prior to manual cleaning and subsequent manual disinfection: Do Pre-cleaning (see pre-cleaning above).</p> <p>Equipments: A multienzymatic detergent (for example, ECOLAB's Sekusept® MultiEnzyme (www.ecolab.com), soft-bristled brushes, a detergent container, and tap water ($20\pm2^\circ\text{C}$) at least meeting drinking water standards as per RKI guidelines.</p> <ol style="list-style-type: none"> 1) The cleaner is prepared in accordance to the manufacturer's instructions (e.g. Sekusept® MultiEnzyme 2% solution). 2) Soak the device into cleaning solution. 3) Exposure time in accordance to the manufacturer's instructions (at least 15 minutes for Sekusept® MultiEnzyme 2% solution). 4) While the instrument is submerged, scrub the exterior surfaces and any crevices with a brush to remove residue. 5) With the brush, clean the instrument's internal surfaces, grooves, and blind holes at least 15 times to eliminate surface contaminants. 6) During exposure time while device is soaked brush the canal inside the surgical instruments for at least 15 times, and brush all other hard-to-reach areas and accessible areas inside the device to remove all visible dirt with a soft brush. <p>Pay attention to the critical and hard-to-reach areas where it is not possible to assess the cleaning efficacy for the devices by means of visual inspection.</p> <ol style="list-style-type: none"> 7) Take the instrument out of the cleaning solution. 8) Rinse the instrument under running tap water for one minute, ensuring all surfaces are thoroughly washed. 9) Inspect the instrument for cleanliness. If contamination is visible, repeat the cleaning steps.
Disinfection: Manually	<p>Equipment: Disinfectant on peroxide basis (e.g. Sekusept® aktiv of ECOLAB (www.ecolab.com), disinfectant tank, tap water ($20\pm2^\circ\text{C}$) of at least drinking water quality (in accordance with RKI guideline), ultrasonic bath (35 kHz), deionized water (in accordance with KRINKO/BfArM-recommendation free of facultative pathogen microorganisms), fluff-free wipe or towel or air gun with compressed medical air (in accordance with European Pharmacopoeia), a sterile 20-ml syringe</p> <ol style="list-style-type: none"> 1) Prepare the disinfectant solution in an ultrasonic cleaner as per the manufacturer's instructions (e.g., a 2% solution of Sekusept® aktiv). 2) Fully submerge the instrument in the disinfectant. 3) Adhere to the manufacturer's recommended exposure time (at least 10 minutes for a 2% solution of Sekusept® aktiv). 4) Use a 20-ml syringe filled with disinfectant to flush the instrument's channels, blind holes, and other hard-to-reach areas at least five times. <ol style="list-style-type: none"> a. Pay attention to areas where visual inspection is not feasible for verifying cleanliness. 5) After 5 minutes exposure time, activate the ultrasonic cleaner for an additional 5 minutes. 6) Switch off the ultrasonic cleaner. 7) Remove the device from the disinfectant. 8) Rinse the device under running tap water for at least 1 minute. 9) Perform a secondary rinse using deionized water. 10) Wipe the device with single use, fluff free wipe or towel or dry by an air gun with compressed air. Cleaning and disinfection: Automated If visible soil remains in

	<p>If visible soil remains in gaps, holes, joining surfaces, threads and other hard to reach areas: Start with Pre-cleaning (see pre-cleaning above) and continue with manual cleaning.</p> <p>Equipment: Washer disinfector in accordance with DIN EN ISO 15883-1 + -2 with thermal program (temperature 90-93°C), mildly alkaline cleaner detergent (e.g. Sekumatic® MultiClean 200 ml), fluff free wipe or towel or air gun connected to compressed air of medical grade (in accordance with European Pharmacopoeia).</p> <ol style="list-style-type: none"> 1) Place all components of the instruments onto a compatible tray within a load carrier module designated for the washer-disinfector, or connect them with their hollow channels to a load carrier module designed for minimally invasive surgery instruments. 2) Once the door is securely closed, initiate the thermal cycle, with program parameters outlined in the following table. 3) Upon completion of the cycle, retrieve the instruments. 4) Inspect the instruments to ensure they are fully dry. If any moisture remains, use a single-use, lint-free cloth or towel to dry them, or alternatively, apply an air gun with compressed air for drying. 																																			
Cleaning and disinfection: Automated	<table border="1"> <thead> <tr> <th>Program Step</th><th>Water</th><th>Dosage</th><th>Time</th><th>Temperature</th></tr> </thead> <tbody> <tr> <td>Pre-Rinse</td><td>Cold</td><td></td><td>5 minutes</td><td></td></tr> <tr> <td>Dosage cleaner</td><td></td><td>According to manufacturer's instructions</td><td></td><td>According to manufacturer's instructions</td></tr> <tr> <td>Cleaning</td><td>Deionized water</td><td></td><td>10 minutes</td><td>55°C</td></tr> <tr> <td>Rinse</td><td>Deionized water</td><td></td><td>2 minutes</td><td></td></tr> <tr> <td>Disinfection</td><td>Deionized water</td><td></td><td>Ao value >30001 (e.g. 5 minutes, 90°C)</td><td></td></tr> <tr> <td>Drying</td><td></td><td></td><td>15 minutes</td><td>Up to 120°C</td></tr> </tbody> </table> <p>Public authorities may adopt other implementing provisions (parameters for disinfection performance) within their field of responsibility</p>	Program Step	Water	Dosage	Time	Temperature	Pre-Rinse	Cold		5 minutes		Dosage cleaner		According to manufacturer's instructions		According to manufacturer's instructions	Cleaning	Deionized water		10 minutes	55°C	Rinse	Deionized water		2 minutes		Disinfection	Deionized water		Ao value >30001 (e.g. 5 minutes, 90°C)		Drying			15 minutes	Up to 120°C
Program Step	Water	Dosage	Time	Temperature																																
Pre-Rinse	Cold		5 minutes																																	
Dosage cleaner		According to manufacturer's instructions		According to manufacturer's instructions																																
Cleaning	Deionized water		10 minutes	55°C																																
Rinse	Deionized water		2 minutes																																	
Disinfection	Deionized water		Ao value >30001 (e.g. 5 minutes, 90°C)																																	
Drying			15 minutes	Up to 120°C																																
Maintenance, Inspection and Testing:	<p>Equipment: Silicone-free steam sterilizable lubricant allowed for use on medical products according to the manufacturer (e.g. Sterilite®, spray of B. BRAUN Aesculap or other fabricates).</p> <ul style="list-style-type: none"> Under adequate lighting, scrutinize all parts of the instrument for any remaining soil. If soils detected, repeat the manual cleaning steps (including pre-cleaning) and either the automated or manual disinfection processes until the components are visibly clean. Verify the instrument for completeness, any damage, and proper functionality. Examine elongated components of the device for any signs of bending or distortion. Do not proceed with any damaged instrument. Apply a silicone-free steam sterilization-compatible lubricant to all movable parts of the instrument. 																																			
Packaging:	<p>Packaging should be in a sterile barrier system compliant with DIN EN ISO 11607 or DIN 58953 standards. (For validation purposes, double-layered sterilization paper wraps were employed)</p>																																			
Sterilization:	<p>Equipment: Sterilization should be performed using a sterilizer that adheres to DIN EN 13060 and/or DIN EN 285 standards.</p> <ul style="list-style-type: none"> Employ a pre-vacuum sterilization process at temperatures ranging from 134°C to 137°C, with a minimum sterilization duration of 3 to 3.5 minutes. Sterilize the device in its disassembled state. <ol style="list-style-type: none"> 1) Place the packaged device in the sterilizer chamber. 2) Activate the sterilization cycle. 3) Once the cycle is complete, carefully extract the device and allow it to cool to room temperature. 4) Ensure the integrity of the package by checking that the seal is secure and the 																																			
Storage:	<p>Adhere to the storage conditions and shelf life guidelines established by the processor.</p> <ul style="list-style-type: none"> Instruments should be stored in dry, clean and humidity-free environment without direct solar irradiation. Do not store them in places where chemicals may emit corrosive fumes or where fluctuating temperature or humidity may cause vapors to condense on the instruments. 																																			
Additional Information:	<p>Processing of medical devices should only be carried out using validated procedures.</p>																																			
Symbols used on Labelling:	<table> <tbody> <tr> <td></td> <td>Manufacturer</td> <td></td> <td>Non-sterile</td> </tr> <tr> <td> EC REP</td> <td>Authorized Representative in the European Community</td> <td></td> <td>Consult instructions for use</td> </tr> <tr> <td> REF</td> <td>Catalog Number</td> <td></td> <td>Keep away from rain</td> </tr> <tr> <td> LOT</td> <td>Lot Number</td> <td></td> <td>Keep away from sunlight</td> </tr> <tr> <td> !</td> <td>See instructions for use</td> <td>SERIAL</td> <td>Serial Number</td> </tr> </tbody> </table>		Manufacturer		Non-sterile	 EC REP	Authorized Representative in the European Community		Consult instructions for use	 REF	Catalog Number		Keep away from rain	 LOT	Lot Number		Keep away from sunlight	 !	See instructions for use	SERIAL	Serial Number															
	Manufacturer		Non-sterile																																	
 EC REP	Authorized Representative in the European Community		Consult instructions for use																																	
 REF	Catalog Number		Keep away from rain																																	
 LOT	Lot Number		Keep away from sunlight																																	
 !	See instructions for use	SERIAL	Serial Number																																	



ENT Surgery Sets



SURTEX
INSTRUMENTS

Crafted to Excellence

Basic ENT Surgery Set

Codes	Description	Quantity
SI-012-15	Allis Intestinal and Tissue Grasping Forcep 4 x 5 Teeth 15.5 cm - 6" Stainless Steel	1
TR-151-15	Baron-Schuknecht Suction Tube With Finger Cutt Off Working Length 75 mm - Ø 1.5 mm Stainless Steel	1
TR-151-20	Baron-Schuknecht Suction Tube With Finger Cutt Off Working Length 75 mm - Ø 2.0 mm Stainless Steel	1
OT-660-01	Billeau Ear Loop 16 cm - 6 1/4" Working End 2.0 mm Stainless Steel	1
OT-660-02	Billeau Ear Loop 16 cm - 6 1/4" Working End 3.0 mm Stainless Steel	1
OT-660-03	Billeau Ear Loop 16 cm - 6 1/4" Working End 4.0 mm Stainless Steel	1
OT-570-10	Buck Ear Curette Fig. 0 - Straight - Blunt 17 cm - 6 3/4" Ø 2.6 mm Stainless Steel	1
OT-570-11	Buck Ear Curette Fig. 1 - Straight - Blunt 17 cm - 6 3/4" Ø 3.0 mm Stainless Steel	1
OT-570-12	Buck Ear Curette Fig. 2 - Straight - Blunt 17 cm - 6 3/4" Ø 3.7 mm Stainless Steel	1
OT-572-10	Buck Ear Curette Fig. 0 - Straight - Sharp 17 cm - 6 3/4" Ø 2.6 mm Stainless Steel	1
OT-572-11	Buck Ear Curette Fig. 1 - Straight - Sharp 17 cm - 6 3/4" Ø 3.0 mm Stainless Steel	1
OT-572-12	Buck Ear Curette Fig. 2 - Straight - Sharp 17 cm - 6 3/4" Ø 3.7 mm Stainless Steel	1
OT-677-01	Day Ear Hook Long 17 cm - 6 3/4" Stainless Steel	1
OT-677-02	Day Ear Hook Short 17 cm - 6 3/4" Stainless Steel	1
OT-255-14	Hartmann Alligator Forcep 14 cm - 5 1/2" Stainless Steel	1
OT-261-25	Hartmann Alligator Forcep Cup Shaped 8.5 cm - 3 1/4" Ø 2.5 mm Stainless Steel	1
SS-440-11	Iris Scissor Straight 11.5 cm - 4 1/2" Stainless Steel	1
SS-441-11	Iris Scissor Curved 11.5 cm - 4 1/2" Stainless Steel	1
OT-211-14	Gruenwald (Lucae) Ear Forcep Bayonet Shaped 14.5 cm - 5 3/4" Stainless Steel	1
OT-283-08	McGee Micro Alligator Forcep Serrated-Straight 8 cm - 3" Jaw Size 4.0 x 0.8 mm Stainless Steel	1
OT-293-08	Mini-McGee Micro Alligator Forcep Serrated-Straight 8 cm - 3" Jaw Size 3.5 x 0.6 mm Stainless Steel	1
OT-311-08	Micro Alligator Forcep Straight - Cup Shaped 8 cm - 3" Cup Size - Jaw Size 1.0 x 0.9 mm - 4.0 mm Stainless Steel	1
OT-255-08	Hartmann Alligator Forcep 8.5 cm - 3 1/4" Stainless Steel	1
OT-555-01	Shapleigh Ear Curette Fig. 1 15 cm - 6" Stainless Steel	1
OT-555-02	Shapleigh Ear Curette Fig. 2 15 cm - 6" Stainless Steel	1
SS-456-11	Strabismus Scissor Straight 11.5 cm - 4 1/2" Stainless Steel	1
SS-457-11	Strabismus Scissor Curved 11.5 cm - 4 1/2" Stainless Steel	1
RH-023-01	Mod. Vienna Nasal Speculum Fig. 1 Blade Length 27 mm Stainless Steel	1
RH-023-02	Mod. Vienna Nasal Speculum Fig. 2 Blade Length 31 mm Stainless Steel	1
RH-023-03	Mod. Vienna Nasal Speculum Fig. 3 Blade Length 35 mm Stainless Steel	1
OT-201-11	Troeltsch (Wilde) Ear Forcep 11 cm / 4 3/8" Stainless Steel	1

Basic Ear Set

Codes	Description	Quantity
SL-012-03	Scalpel Handle No. 3 Solid, cm-Graduation 12 cm - 4 3/4" Stainless Steel	1
SS-861-14	UltraCut TC Mayo Dissecting Scissor Curved 14.5 cm - 5 3/4" Stainless Steel	1
SS-921-14	UltraCut TC Metzenbaum-Fine Dissecting Scissor - Slender Pattern Curved 14.5 cm - 5 3/4" Stainless Steel	1
SS-923-18	UltraCut TC Metzenbaum-Fine Dissecting Scissor - Slender Pattern Curved 18 cm - 7" Stainless Steel	1
AF-160-12	Halsted-Mosquito Haemostatic Forcep Straight 12.5 cm - 5" Stainless Steel	6
AF-200-14	Kelly Haemostatic Forcep Straight 14.5 cm - 5 3/4" Stainless Steel	2
SU-586-13	UltraGripX TC Halsey Needle Holder 13 cm - 5" Stainless Steel	1
SS-991-12	UltraCut TC Universal Wire Cutting Scissor Angled - One Toothing Cutting Edge 12.5 cm - 5" Stainless Steel	1
SI-012-15	Allis Intestinal and Tissue Grasping Forcep 4 x 5 Teeth 15.5 cm - 6" Stainless Steel	2
RT-716-13	Weitlaner-Wullstein Self Retaining Retractor 3 x 3 Sharp Prongs 13 cm - 5" Stainless Steel	1
RH-475-18	Freer Septum Elevator Double Ended - Sharp/Blunt 18.5 cm - 7 1/4" Blade Size 1 - Blade Size 2 5.0 mm - 5.0 mm Stainless Steel	1
RT-034-15	Joseph Mucus Hook Sharp 15.5 cm - 6" Stainless Steel	2
RT-677-10	Jansen Self Retaining Retractor 3 x 3 Blunt Prongs - With Adjustment Screw 10 cm - 4" Stainless Steel	1
AF-773-08	Backhaus Towel Clamp 9 cm - 3 1/2" Stainless Steel	4
TR-101-07	Frazier Suction Tube 7 Charr. 19.5 cm - 7 3/4" Ø 2.3 mm Stainless Steel	1
TR-101-09	Frazier Suction Tube 9 Charr. 19.5 cm - 7 3/4" Ø 3.0 mm Stainless Steel	1
FR-050-12	Adson Dressing Forcep 12 cm - 4 3/4" Stainless Steel	1
FR-250-12	Adson Dissecting Forcep 1 x 2 Teeth 12 cm - 4 3/4" Stainless Steel	1
RT-251-01	Langenbeck Retractor 22.5 cm - 8 3/4" Blade Size 30 x 11 mm Stainless Steel	1
RT-251-03	Langenbeck Retractor 22.5 cm - 8 3/4" Blade Size 30 x 16 mm Stainless Steel	1
RT-162-16	Senn-Miller Retractor Sharp 16 cm - 6 1/4" Blade Size 20 x 7 mm Stainless Steel	2
RT-672-10	Allport Wound Spreader 4 x 4 Sharp Prongs 9.5 cm - 3 3/4" Stainless Steel	1
SS-457-11	Strabismus Scissor Curved 11.5 cm - 4 1/2" Stainless Steel	1
FR-772-01	Graefe Fixation Forcep Without Lock 11 cm - 4 1/4" Stainless Steel	1
OT-213-16	Jansen Ear Forcep Bayonet Shaped 16 cm - 6 1/4" Stainless Steel	1
OT-572-02	Buck Ear Curette Fig. 00 - Straight - Sharp 17 cm - 6 3/4" Ø 2.3 mm Stainless Steel	1
OT-572-10	Buck Ear Curette Fig. 0 - Straight - Sharp 17 cm - 6 3/4" Ø 2.6 mm Stainless Steel	1
OT-010-01	Hartmann Ear Specula Set of 3 Ref:- OT-012-00 to OT-012-02 Stainless Steel	1

Basic Nasal Surgery Set

Codes	Description	Quantity
SL-010-03	Bard-Parker Scalpel Handle No. 3 Solid 12 cm - 4 3/4" Stainless Steel	1
SL-030-07	Scalpel Handle No. 7 Solid 16 cm - 6 1/4" Stainless Steel	1
SS-041-14	Mayo Dissecting Scissor Curved 14.5 cm - 5 3/4" Stainless Steel	1
SS-201-14	Metzenbaum-Fino Delicate Dissecting Scissor Curved - Blunt/Blunt Slender Pettern 14.5 cm - 5 3/4" Stainless Steel	1
SS-215-18	Metzenbaum-Fino Delicate Dissecting Scissor Curved - Sharp/sharp Slender Pettern 18 cm - 7" Stainless Steel	2
SS-441-11	Iris Scissor Curved 11.5 cm - 4 1/2" Stainless Steel	2
FR-250-12	Adson Dissecting Forcep 1 x 2 Teeth 12 cm - 4 3/4" Stainless Steel	2
AF-160-12	Halsted-Mosquito Haemostatic Forcep Straight 12.5 cm - 5" Stainless Steel	6
AF-161-12	Halsted-Mosquito Haemostatic Forcep Curved 12 cm - 4 3/4" Stainless Steel	6
AF-773-08	Backhaus Towel Clamp 9 cm - 3 1/2" Stainless Steel	6
RT-162-16	Senn-Miller Retractor Sharp 16 cm - 6 1/4" Blade Size 20 x 7 mm Stainless Steel	1
RT-163-16	Senn-Miller Retractor Blunt 16 cm - 6 1/4" Blade Size 20 x 7 mm Stainless Steel	1
RT-579-21	Parker-Langebeck (US Army) Retractor Set 21.5 cm Stainless Steel	1
RT-677-10	Jansen Self Retaining Retractor 3 x 3 Blunt Prongs - With Adjustment Screw 10 cm - 4" Stainless Steel	1
TR-101-08	Frazier Suction Tube 8 Charr. 19.5 cm - 7 3/4" Ø 2.6 mm Stainless Steel	2
TR-251-27	Yankauer Suction Tube Complete Central Bore 2.0 mm 31 cm - 12 1/4" Stainless Steel	1
SU-586-13	UltraGripX TC Halsey Needle Holder 13 cm - 5" Stainless Steel	1
OR-504-19	FiberGrip Vickers Mallet With Interchangeable Nylon Discs 19 cm - 7 1/2" Head Ø 25.0 mm - Weight 300 Grams Stainless Steel	1
SI-012-15	Allis Intestinal and Tissue Grasping Forcep 4 x 5 Teeth 15.5 cm - 6" Stainless Steel	2
FR-772-01	Graefe Fixation Forcep Without Lock 11 cm - 4 1/4" Stainless Steel	1
OT-211-14	Gruenwald (Lucae) Ear Forcep Bayonet Shaped 14.5 cm - 5 3/4" Stainless Steel	3
RH-033-02	Hartmann-Halle Nasal Speculum Fig. 2 15.5 cm - 6" Blade Length 31 mm Stainless Steel	1
RH-101-15	Troeltsch Nasal Tampon Forcep 14.5 cm - 5 3/4" Stainless Steel	1
RT-034-15	Joseph Mucus Hook Sharp 15.5 cm - 6" Stainless Steel	2
RT-074-02	Joseph Mucosa Hook 16.5 cm - 6 1/2" Width 2.0 mm Stainless Steel	2
RT-074-07	Joseph Mucosa Hook 16.5 cm - 6 1/2" Width 7.0 mm Stainless Steel	2

Major Nasal Set

Codes	Description	Quantity
AF-773-08	Backhaus Towel Clamp 9 cm - 3 1/2" Stainless Steel	6
RH-380-03	Ballanger Swivel Knife Straight 21 cm - 8 1/4" Cutting Edge Size 3.0 mm Stainless Steel	1
RH-820-11	Cottle Columella Forcep 11 cm - 4 1/4" Width 12.0 mm Stainless Steel	1
RH-771-01	Cottle Alar Hook Sharp 14.5 cm - 5 3/4" Stainless Steel	1
RH-367-15	Cottle Rhinoplastic Knife 14 cm - 5 1/2" Stainless Steel	1
RH-621-04	Cottle Osteotome Curved 18.5 cm - 7 1/4" Blade Width 4.0 mm Stainless Steel	1
RT-080-14	Cottle Alar Hook 4 Semi Sharp Prongs 14 cm - 5 1/2" Stainless Steel	1
RT-120-08	Desmarres Lid Retractor 16 cm - 6 1/4" Blade Width 8.0 mm Stainless Steel	1
SS-357-15	Fomon Tonsil Scissor 15 cm - 6" Stainless Steel	1
RH-355-15	Fomon Rhinoplastic Knife 16.5 cm - 6 1/2" Stainless Steel	1
SS-189-14	Fomon Nasal Scissor Curved 14 cm - 5 1/2" Stainless Steel	1
RH-765-16	Fomon Alar Hook Button Ended 16 cm - 6 1/4" Stainless Steel	1
TR-101-08	Frazier Suction Tube 8 Charr. 19.5 cm - 7 3/4" Ø 2.6 mm Stainless Steel	2
RH-475-18	Freer Septum Elevator Double Ended - Sharp/Blunt 18.5 cm - 7 1/4" Blade Size 1 - Blade Size 2 5.0 mm - 5.0 mm Stainless Steel	1
RT-012-02	Gillies Skin Hook Fig. 2 18 cm - 7" Stainless Steel	1
FR-772-01	Graefe Fixation Forcep Without Lock 11 cm - 4 1/4" Stainless Steel	1
AF-160-12	Halsted-Mosquito Haemostatic Forcep Straight 12.5 cm - 5" Stainless Steel	6
AF-161-12	Halsted-Mosquito Haemostatic Forcep Curved 12 cm - 4 3/4" Stainless Steel	6
RH-127-16	Hartmann Nasal Polypus Forcep 15 cm - 6" Stainless Steel	1
SS-440-11	Iris Scissor Straight 11.5 cm - 4 1/2" Stainless Steel	1
RT-677-10	Jansen Self Retaining Retractor 3 x 3 Blunt Prongs - With Adjustment Screw 10 cm - 4" Stainless Steel	1
RH-551-16	Joseph Periosteal Elevator 18 cm - 7" Blade Width 6.5 mm Stainless Steel	1
RT-034-15	Joseph Mucus Hook Sharp 15.5 cm - 6" Stainless Steel	2
RT-074-02	Joseph Mucosa Hook 16.5 cm - 6 1/2" Width 2.0 mm Stainless Steel	2
RT-074-05	Joseph Mucosa Hook 16.5 cm - 6 1/2" Width 5.0 mm Stainless Steel	1
RT-074-07	Joseph Mucosa Hook 16.5 cm - 6 1/2" Width 7.0 mm Stainless Steel	2
PR-500-16	Farrell Cotton Applicator Round Spiral End 16 cm - 6 1/4" Working End Ø 1.1 Stainless Steel	2
OT-211-14	Gruenwald (Lucae) Ear Forcep Bayonet Shaped 14.5 cm - 5 3/4" Stainless Steel	3
OR-504-19	FiberGrip Vickers Mallet With Interchangeable Nylon Discs 19 cm - 7 1/2" Head Ø 25.0 mm - Weight 300 Grams Stainless Steel	1
SS-041-14	Mayo Dissecting Scissor Curved 14.5 cm - 5 3/4" Stainless Steel	1
SS-201-14	Metzenbaum-Fino Delicate Dissecting Scissor Curved - Blunt/Blunt Slender Pettern 14.5 cm - 5 3/4" Stainless Steel	1
SS-141-18	Metzenbaum-Nelson Dissecting Scissor Curved - Blunt/Blunt 18 cm - 7" Stainless Steel	2
RH-023-01	Mod. Vienna Nasal Speculum Fig. 1 Blade Length 27 mm Stainless Steel	1
RH-023-02	Mod. Vienna Nasal Speculum Fig. 2 Blade Length 31 mm Stainless Steel	1
RH-789-16	Neivert Nasal Retractor 17 cm - 6 3/4" Stainless Steel	1
SL-010-03	Bard-Parker Scalpel Handle No. 3 Solid 12 cm - 4 3/4" Stainless Steel	1
SL-030-07	Scalpel Handle No. 7 Solid 16 cm - 6 1/4" Stainless Steel	1
RT-162-16	Senn-Miller Retractor Sharp 16 cm - 6 1/4" Blade Size 20 x 7 mm Stainless Steel	1

Major Nasal Set

Codes	Description	Quantity
RT-163-16	Senn-Miller Retractor Blunt 16 cm - 6 1/4" Blade Size 20 x 7 mm Stainless Steel	1
RH-664-04	Sheehan Osteotome 15 cm - 6" Blade Width 4.0 mm Stainless Steel	1
RH-664-08	Sheehan Osteotome 15 cm - 6" Blade Width 8.0 mm Stainless Steel	1
RH-664-10	Sheehan Osteotome 15 cm - 6" Blade Width 10.0 mm Stainless Steel	1
RH-664-12	Sheehan Osteotome 15 cm - 6" Blade Width 12.0 mm Stainless Steel	1
SS-456-11	Strabismus Scissor Straight 11.5 cm - 4 1/2" Stainless Steel	1
SS-457-11	Strabismus Scissor Curved 11.5 cm - 4 1/2" Stainless Steel	1
FR-078-17	Taylor Dressing Forcep Bayonet Shaped - With dissector End 18.5 cm - 7 1/4" Stainless Steel	2
FR-962-12	UltraGripX TC Adson Dissecting Forcep 1 x 2 Teeth 12 cm - 4 3/4" Stainless Steel	2
SU-594-18	UltraGripX TC Crile-Wood Needle Holder 18 cm - 7" Stainless Steel	1
SU-586-13	UltraGripX TC Halsey Needle Holder 13 cm - 5" Stainless Steel	1
RH-101-15	Troeltsch Nasal Tampon Forcep 14.5 cm - 5 3/4" Stainless Steel	1
RT-579-21	Parker-Langebeck (US Army) Retractor Set 21.5 cm Stainless Steel	1
TR-251-27	Yankauer Suction Tube Complete Central Bore 2.0 mm 31 cm - 12 1/4" Stainless Steel	1
AF-269-18	Adson Haemostatic Forcep Curved 18.5 cm - 7 1/4" Stainless Steel	1

Minor Nasal Surgery Set

Codes	Description	Quantity
OR-504-19	FiberGrip Vickers Mallet With Interchangeable Nylon Discs 19 cm - 7 1/2" Head √ø 25.0 mm - Weight 300 Grams Stainless Steel	1
SL-010-03	Bard-Parker Scalpel Handle No. 3 Solid 12 cm - 4 3/4" Stainless Steel	1
SL-030-07	Scalpel Handle No. 7 Solid 16 cm - 6 1/4" Stainless Steel	1
RT-162-16	Senn-Miller Retractor Sharp 16 cm - 6 1/4" Blade Size 20 x 7 mm Stainless Steel	1
RT-163-16	Senn-Miller Retractor Blunt 16 cm - 6 1/4" Blade Size 20 x 7 mm Stainless Steel	1
SU-586-13	UltraGripX TC Halsey Needle Holder 13 cm - 5" Stainless Steel	1
RH-101-15	Troeltsch Nasal Tampon Forcep 14.5 cm - 5 3/4" Stainless Steel	1
RT-579-21	Parker-Langebeck (US Army) Retractor Set 21.5 cm Stainless Steel	1
TR-251-27	Yankauer Suction Tube Complete Central Bore 2.0 mm 31 cm - 12 1/4" Stainless Steel	1

Nasal Fracture Reduction Set

Codes	Description	Quantity
OT-213-16	Jansen Ear Forcep Bayonet Shaped 16 cm - 6 1/4" Stainless Steel	1
RH-839-21	Faulkner Antrum Curette 21 cm - 8 1/4" Stainless Steel	1
RH-023-01	Mod. Vienna Nasal Speculum Fig. 1 Blade Length 27 mm Stainless Steel	1
RH-807-23	Asch Septum Straightening Forcep 22.5 cm - 8 3/4" Stainless Steel	1
RH-804-02	Cottle-Walsham Septum Straightening Forcep Left 23 cm - 9" Stainless Steel	1
RH-804-03	Cottle-Walsham Septum Straightening Forcep Right 23 cm - 9" Stainless Steel	1
OM-607-19	Boise Nasal Fracture Elevator 19.5 cm - 7 3/4" Stainless Steel	1

Nasal Polypectomy Set

Codes	Description	Quantity
RH-011-20	Thudichum Nasal Specula Set of 4 Ref. RH-013-01 to RH-013-04 Stainless Steel	1
RH-101-15	Troeltsch Nasal Tampon Forcep 14.5 cm - 5 3/4" Stainless Steel	1
RH-091-26	Krause Polypus Snare 26 cm - 10 1/4" Shaft Ø 3.0 mm Stainless Steel	1
RH-097-12	Snare Wire Bundle with 12 Pieces 10.5 cm / 4" Ø 0.3 mm Stainless Steel	1
TR-111-15	Fergusson Suction Tube 15 Charr. 18 cm - 7" Ø 5.0 mm Stainless Steel	1

Basic Throat Surgery Set

Codes	Description	Quantity
AF-346-18	Birkett Haemostatic Forcep Straight 18.5 cm - 7 1/4" Stainless Steel	2
SF-350-25	Foerster Sponge Holding Forcep Straight 25 cm - 9 3/4" Stainless Steel	2
AF-315-01	Negus Haemostatic Forcep Fig. 1 19.5 cm - 7 3/4" Stainless Steel	2
OM-159-19	Bruenings Tongue Depressor 19 cm - 7 1/2" Blade Width 1 - Blade Width 2 19 mm - 15 mm Stainless Steel	1
RH-475-18	Freer Septum Elevator Double Ended - Sharp/Blunt 18.5 cm - 7 1/4" Blade Size 1 - Blade Size 2 5.0 mm - 5.0 mm Stainless Steel	1
TL-037-23	Henke Tonsil Elevator 22.5 cm - 8 3/4" Blade Width 10.5 mm - 5.0 mm Stainless Steel	1
OR-623-19	Kraemer Perosteal Respiratory / Elevator 19 cm - 7 1/2" Width 10 mm Stainless Steel	1
PR-015-14	Probe 14.5 cm - 5 3/4" Stainless Steel	1
PR-100-14	Butterfly Probe / Grooved Director 14 cm - 5 1/2" Stainless Steel	1
RT-716-13	Weitlaner-Wullstein Self Retaining Retractor 3 x 3 Sharp Prongs 13 cm - 5" Stainless Steel	1
RT-712-16	Weitlaner Self Retaining Retractor 3 x 4 Sharp Prongs 16 cm - 6 1/4" Stainless Steel	1
SS-341-21	Killian Bone Shear 21 cm - 8 1/4" Stainless Steel	1
TR-111-09	Fergusson Suction Tube 9 Charr. 18 cm - 7" \varnothing 3.0 mm Stainless Steel	1
TR-111-12	Fergusson Suction Tube 12 Charr. 18 cm - 7" \varnothing 4.0 mm Stainless Steel	1
TR-111-15	Fergusson Suction Tube 15 Charr. 18 cm - 7" \varnothing 5.0 mm Stainless Steel	1
TL-105-20	Blohmke Tonsil Grasping Forcep 1 x 2 Teeth 20 cm - 8" Stainless Steel	1
RH-241-02	Hartmann Nasal Cutting Forcep (Conchotome) Fig. 2 12 cm - 4 3/4" \varnothing 7.0 mm Stainless Steel	1
HW-300-50	Medicine Cup Capacity 50 cc Stainless Steel	1
TR-251-27	Yankauer Suction Tube Complete Central Bore 2.0 mm 31 cm - 12 1/4" Stainless Steel	12
RT-100-01	Wound Retractor 1 Sharp Prong 16.5 cm - 6 1/2" Stainless Steel	2
RT-100-02	Wound Retractor 2 Sharp Prongs - Small Curve 16.5 cm - 6 1/2" Width 4.2 mm Stainless Steel	2
RT-210-02	Volkmann Retractor 2 Sharp Prongs 22 cm - 8 3/4" Stainless Steel	2
RT-210-04	Volkmann Retractor 4 Sharp Prongs 22 cm - 8 3/4" Stainless Steel	2
RT-210-06	Volkmann Retractor 6 Sharp Prongs 22 cm - 8 3/4" Stainless Steel	2
RT-211-02	Volkmann Retractor 2 Blunt Prongs 22 cm - 8 3/4" Stainless Steel	2
RT-211-04	Volkmann Retractor 4 Blunt Prongs 22 cm - 8 3/4" Stainless Steel	2
RT-211-06	Volkmann Retractor 6 Blunt Prongs 22 cm - 8 3/4" Stainless Steel	1
RT-335-20	Shoenborn Retractor 20 cm - 8" Blade Size 13 x 6 mm Stainless Steel	1
RT-251-02	Langenbeck Retractor 22.5 cm - 8 3/4" Blade Size 30 x 14 mm Stainless Steel	1
RT-281-20	Kocher Retractor 23 cm - 9" Blade Size 61 x 20 mm Stainless Steel	1
FR-012-14	Slender Pattern Dressing Forcep 14.5 cm - 5 3/4" Stainless Steel	1
FR-016-14	Fine Pattern Dressing Forcep 14.5 cm - 5 3/4" Stainless Steel	1
FR-200-14	Standard Pattern Dissecting Forcep 1 x 2 Teeth 14.5 cm - 5 3/4" Stainless Steel	1
FR-214-14	Fine Pattern Dissecting Forcep 1 x 2 Teeth 14.5 cm - 5 3/4" Stainless Steel	1
FR-010-20	Standard Pattern Dressing Forcep 20 cm - 8" Stainless Steel	1
FR-200-20	Standard Pattern Dissecting Forcep 1 x 2 Teeth 20 cm - 8" Stainless Steel	1
FR-590-16	DeBakey Atrauma Forcep 16 cm - 6 1/4" Tip Size 1.5 mm Stainless Steel	1
AF-311-02	Schnidt Tonsil Haemostatic Forcep Gently Curved 19 cm - 7 1/2" Stainless Steel	2
SL-010-03	Bard-Parker Scalpel Handle No. 3 Solid 12 cm - 4 3/4" Stainless Steel	1
SL-010-04	Bard-Parker Scalpel Handle No. 4 Solid 13.5 cm - 5 3/8" Stainless Steel	1
TR-705-10	Redon Guide Needle 10 Charr. - Lancet Tip 19.5 cm - 7 3/4" Tip Size 3.3 mm Stainless Steel	1
SU-875-20	Deschamps Ligature Needle Blunt for Right Hand Medium 20 cm - 8" Stainless Steel	1
SU-877-20	Deschamps Ligature Needle Blunt for Left Hand Medium 20 cm - 8" Stainless Steel	1
SU-586-13	UltraGripX TC Halsey Needle Holder 13 cm - 5" Stainless Steel	1
SU-594-15	UltraGripX TC Crile-Wood Needle Holder 15 cm - 6" Stainless Steel	1
SU-602-20	UltraGripX TC Hegar Needle Holder 20 cm - 8" Stainless Steel	1
OT-211-14	Gruenwald (Lucae) Ear Forcep Bayonet Shaped 14.5 cm - 5 3/4" Stainless Steel	2
FR-050-12	Adson Dressing Forcep 12 cm - 4 3/4" Stainless Steel	1
FR-250-12	Adson Dissecting Forcep 1 x 2 Teeth 12 cm - 4 3/4" Stainless Steel	1
FR-258-12	Adson-Brown Dissecting Forcep 7 x 7 Teeth 12 cm - 5" Stainless Steel	2
AF-161-12	Halsted-Mosquito Haemostatic Forcep Curved 12 cm - 4 3/4" Stainless Steel	20
AF-163-12	Halsted-Mosquito Haemostatic Forcep Curved - 1 x 2 Teeth 12 cm - 4 3/4" Stainless Steel	20
AF-773-11	Backhaus Towel Clamp 11.5 cm - 4 1/2" Stainless Steel	20
AF-220-14	Rochester-Pean Haemostatic Forcep Straight 14.5 cm - 5 3/4" Stainless Steel	12
AF-220-16	Rochester-Pean Haemostatic Forcep Straight 16.5 cm - 6 1/2" Stainless Steel	4
AF-230-14	Ochsner-Kocher Haemostatic Forcep Straight - 1 x 2 Teeth 14.5 cm - 5 3/4" Stainless Steel	12
SI-015-15	Allis Intestinal and Tissue Grasping Forcep 5 x 6 Teeth 15.5 cm - 6" Stainless Steel	2
GY-616-16	Braun Tenaculum Forcep 16 cm - 6 1/4" Stainless Steel	2
SS-251-15	Reynolds Dissecting Scissor Curved 15.5 cm - 6" Stainless Steel	1
SS-131-14	Metzenbaum Dissecting Scissor Curved - Blunt/Blunt 14.5 cm - 5 3/4" Stainless Steel	2
SS-141-18	Metzenbaum-Nelson Dissecting Scissor Curved - Blunt/Blunt 18 cm - 7" Stainless Steel	2
SS-931-14	UltraCut- \ddot{o} TC Metzenbaum-Lahey Dissecting Scissor - Slender Pattern Curved 14.5 cm - 5 3/4" Stainless Steel	1
SS-453-10	Stevens Tenotomy Scissor Curved - Blunt/Blunt 10 cm - 4" Stainless Steel	1
SS-105-15	Joseph Dissecting Scissor / Opeating Scissor Curved 15 cm - 6" Stainless Steel	2
SS-423-09	Very Delicate Operating Scissor Curved 9 cm - 3 1/2" Stainless Steel	1
SS-004-14	Operating Scissor Straight - Sharp/Blunt 14.5 cm - 5 3/4" Stainless Steel	1
SS-005-14	Operating Scissor Curved - Sharp/Blunt 14.5 cm - 5 3/4" Stainless Steel	1

Tonsillectomy & Adenoidectomy

Codes	Description	Quantity
SL-030-07	Scalpel Handle No. 7 Solid 16 cm - 6 1/4" Stainless Steel	1
SS-141-18	Metzenbaum-Nelson Dissecting Scissor Curved - Blunt/Blunt 18 cm - 7" Stainless Steel	1
AF-269-18	Adson Haemostatic Forcep Curved 18.5 cm - 7 1/4" Stainless Steel	1
AF-773-13	Backhaus Towel Clamp 13 cm - 5" Stainless Steel	2
SF-350-18	Foerster Sponge Holding Forcep Straight 18 cm - 7" Stainless Steel	2
TR-251-27	Yankauer Suction Tube Complete Central Bore 2.0 mm 31 cm - 12 1/4" Stainless Steel	1
SU-650-20	UltraGripX TC DeBakey Needle Holder 20 cm - 8" Stainless Steel	1
SI-012-15	Allis Intestinal and Tissue Grasping Forcep 4 x 5 Teeth 15.5 cm - 6" Stainless Steel	2
RH-851-01	Beckmann Adenoid Curette Fig. 1 22.5 cm - 8 3/4" Width 13.0 mm Stainless Steel	1
RH-851-02	Beckmann Adenoid Curette Fig. 2 22.5 cm - 8 3/4" Width 15.0 mm Stainless Steel	1
RH-851-03	Beckmann Adenoid Curette Fig. 3 22.5 cm - 8 3/4" Width 17.0 mm Stainless Steel	1
RH-851-04	Beckmann Adenoid Curette Fig. 4 22.5 cm - 8 3/4" Width 19.0 mm Stainless Steel	1
RH-851-05	Beckmann Adenoid Curette Fig. 5 22.5 cm - 8 3/4" Width 21.0 mm Stainless Steel	1
RH-851-06	Beckmann Adenoid Curette Fig. 6 22.5 cm - 8 3/4" Width 23.0 mm Stainless Steel	1
RH-853-01	Beckmann Adenoid Curette Angled Shaft - Fig. 1 22 cm - 8 3/4" Width 13.0 mm Stainless Steel	1
RH-853-02	Beckmann Adenoid Curette Angled Shaft - Fig. 2 22 cm - 8 3/4" Width 15.0 mm Stainless Steel	1
RH-853-03	Beckmann Adenoid Curette Angled Shaft - Fig. 3 22 cm - 8 3/4" Width 17.0 mm Stainless Steel	1
RH-853-04	Beckmann Adenoid Curette Angled Shaft - Fig. 4 22 cm - 8 3/4" Width 19.0 mm Stainless Steel	1
RH-853-05	Beckmann Adenoid Curette Angled Shaft - Fig. 5 22 cm - 8 3/4" Width 21.0 mm Stainless Steel	1
RH-853-06	Beckmann Adenoid Curette Angled Shaft - Fig. 6 22 cm - 8 3/4" Width 23.0 mm Stainless Steel	1
OM-080-00	McIvor Mouth Gag Complete With 3 Tongue Depressors Ref:- OM-083-01 to OM-083-03 Stainless Steel	1
OM-147-01	Wieder Tongue Depressor 14 cm - 5 1/2" Blade Width 36 mm Stainless Steel	1
OM-147-02	Wieder Tongue Depressor 15 cm - 6" Blade With 43 mm Stainless Steel	1
TL-011-21	Fischer Tonsil Knife / Dissector 21.5 cm - 8 1/2" Stainless Steel	1
TL-031-22	Hurd Tonsil Dissector / Retractor 22.5 cm - 8 3/4" Blade Width 11.0 mm - 14.0 mm Stainless Steel	1
TL-101-18	White Tonsil Grasping Forcep 18 cm - 7" Stainless Steel	1
AF-311-02	Schnidt Tonsil Haemostatic Forcep Gently Curved 19 cm - 7 1/2" Stainless Steel	1
TL-150-28	Eves Tonsil Snare Without Ratchet 28.5 cm - 11 1/4" Stainless Steel	1

Tympanoplasty And Stapedectomy Set

Codes	Description	Quantity
OT-731-03	Lempert Periosteal Elevator 17 cm - 6 3/4" Width 3 mm Stainless Steel	1
OT-415-15	Plester Sickle Knife 15 cm - 6" Stainless Steel	1
OT-417-15	Sexton Myringotomy Knife 15 cm - 6" Stainless Steel	1
OT-409-15	Plester Flap Knife 15 cm - 6" Stainless Steel	1
OT-675-04	Wagener Ear Hook Fig. 4 14 cm - 5 1/2" Stainless Steel	1
OT-505-03	Fisch Micro Dissector Angled to Right 15.5 cm - 6" Stainless Steel	1
OT-453-16	Barbara Micro Ear Pick 16 cm - 6 1/4" Stainless Steel	1
OT-701-16	Politzer Furuncle Knife 14.5 cm - 5 3/4" Stainless Steel	1
OT-463-10	Barbara Micro Ear Needle Angled 90° 16 cm - 6 1/4" Tip Size 1.0 mm Stainless Steel	1
OT-520-01	House Micro Ear Curette 14.5 cm - 5 3/4" Cup Size 1 / Cup Size 2 1.0 mm - 1.2 mm Stainless Steel	1
OT-520-02	House Micro Ear Curette 14.5 cm - 5 3/4" Cup Size 1 / Cup Size 2 1.5 mm - 1.8 mm Stainless Steel	1
PR-040-10	Stacke (Silver) Probe 10 cm - 4" Tip v̄o 1 mm Stainless Steel	1
OT-311-08	Micro Alligator Forcep Straight - Cup Shaped 8 cm - 3" Cup Size - Jaw Size 1.0 x 0.9 mm - 4.0 mm Stainless Steel	1
OT-303-08	Micro Alligator Forcep Serrated-Straight 8 cm - 3" Jaw Size 4.0 x 0.6 mm Stainless Steel	1
OT-365-01	House-Dieter Malleus Nipper Up Cutting 8 cm - 3" Jaw Opening 1.3 mm Stainless Steel	1
OT-331-08	Bellucci Micro Scissor Straight - Very Fine 8 cm - 3" Blade Size 4.0 x 0.8 mm Stainless Steel	1
TR-101-06	Frazier Suction Tube 6 Charr. 19.5 cm - 7 3/4" Ø 2.0 mm Stainless Steel	1
TR-101-08	Frazier Suction Tube 8 Charr. 19.5 cm - 7 3/4" Ø 2.6 mm Stainless Steel	1
TR-101-09	Frazier Suction Tube 9 Charr. 19.5 cm - 7 3/4" Ø 3.0 mm Stainless Steel	1
TR-141-01	Wullstein Suction Handle 14 cm - 5 1/2" Stainless Steel	1
TR-185-05	Verhoeven Suction Cannula With Luer Cone 10 cm - 4" Diameter 0.5 mm Ø Stainless Steel	1
TR-185-07	Verhoeven Suction Cannula With Luer Cone 10 cm - 4" Diameter 0.7 mm Ø Stainless Steel	1
TR-185-10	Verhoeven Suction Cannula With Luer Cone 10 cm - 4" Diameter 1.0 mm Ø Stainless Steel	1
TR-185-12	Verhoeven Suction Cannula With Luer Cone 10 cm - 4" Diameter 1.2 mm Ø Stainless Steel	1
RT-732-13	Plester Self Retaining Retractor 2 x 2 Sharp Prongs 13 cm - 5" Stainless Steel	1
RT-735-01	Plester Self Retaining Retractor Left 10.5 cm - 4" Stainless Steel	1
RT-735-02	Plester Self Retaining Retractor Right 10.5 cm - 4" Stainless Steel	1
OT-421-04	Rosen Circular Cut Knife Fig. 4 15 cm - 6" v̄o 2.0 mm Stainless Steel	1
OT-463-03	Barbara Micro Ear Needle Angled 90° 16 cm - 6 1/4" Tip Size 0.3 mm Stainless Steel	1
OT-463-10	Barbara Micro Ear Needle Angled 90° 16 cm - 6 1/4" Tip Size 1.0 mm Stainless Steel	1
OT-571-02	Buck Ear Curette Fig. 00 - Angled - Blunt 17 cm - 6 3/4" Ø 2.3 mm Stainless Steel	1
OT-571-11	Buck Ear Curette Fig. 1 - Angled - Blunt 17 cm - 6 3/4" Ø 3.0 mm Stainless Steel	1
OT-283-08	McGee Micro Alligator Forcep Serrated-Straight 8 cm - 3" Jaw Size 4.0 x 0.8 mm Stainless Steel	1
AF-160-12	Halsted-Mosquito Haemostatic Forcep Straight 12.5 cm - 5" Stainless Steel	1
TR-101-07	Frazier Suction Tube 7 Charr. 19.5 cm - 7 3/4" Ø 2.3 mm Stainless Steel	1
OT-421-03	Rosen Circular Cut Knife Fig. 3 15 cm - 6" Ø 2.6 mm Stainless Steel	1
OT-421-07	Rosen Circular Cut Knife Fig. 7 15 cm - 6" Ø 1.0 mm Stainless Steel	1

Tympanoplasty And Stapedectomy Set

Codes	Description	Quantity
OT-505-01	Fisch Micro Dissector Angled to Right 15.5 cm - 6" Stainless Steel	1
OT-505-02	Fisch Micro Dissector Angled to Left 15.5 cm - 6" Stainless Steel	1

Myringotomy Set

Codes	Description	Quantity
SS-440-11	Iris Scissor Straight 11.5 cm - 4 1/2" Stainless Steel	1
TR-151-10	Baron-Schuknecht Suction Tube With Finger Cutt Off Working Length 75 mm Ø 1.0 mm Stainless Steel	1
TR-151-15	Baron-Schuknecht Suction Tube With Finger Cutt Off Working Length 75 mm Ø 1.5 mm Stainless Steel	1
TR-151-20	Baron-Schuknecht Suction Tube With Finger Cutt Off Working Length 75 mm Ø 2.0 mm Stainless Steel	1
OT-223-12	Hartmann Ear Polypus Forcep Fine Model 12 cm - 4 3/4" Stainless Steel	1
OT-201-13	Troeltsch (Wilde) Ear Forcep 12.5 cm - 5" Stainless Steel	1
OT-704-16	Politzer Tympanum Needle 17.5 cm - 7" Stainless Steel	1
OT-570-10	Buck Ear Curette Fig. 0 - Straight - Blunt 17 cm - 6 3/4" Ø 2.6 mm Stainless Steel	1
OT-570-11	Buck Ear Curette Fig. 1 - Straight - Blunt 17 cm - 6 3/4" Ø 3.0 mm Stainless Steel	1
OT-572-10	Buck Ear Curette Fig. 0 - Straight - Sharp 17 cm - 6 3/4" Ø 2.6 mm Stainless Steel	1
OT-572-11	Buck Ear Curette Fig. 1 - Straight - Sharp 17 cm - 6 3/4" Ø 3.0 mm Stainless Steel	1
OT-010-01	Hartmann Ear Specula Set of 3 Ref:- OT-012-00 to OT-012-02 Stainless Steel	1
PR-500-16	Farrell Cotton Applicator Round Spiral End 16 cm - 6 1/4" Working End Ø 1.1 Stainless Steel	1
PR-520-16	Farrell Cotton Applicator Triangular Shaped End 16 cm - 6 1/4" Working End Ø 1.1 mm Stainless Steel	1

Middle Ear Set

Codes	Description	Quantity
TR-151-10	Baron-Schuknecht Suction Tube With Finger Cutt Off Working Length 75 mm - Ø 1.0 mm Stainless Steel	1
TR-151-15	Baron-Schuknecht Suction Tube With Finger Cutt Off Working Length 75 mm - Ø 1.5 mm Stainless Steel	1
TR-151-20	Baron-Schuknecht Suction Tube With Finger Cutt Off Working Length 75 mm - Ø 2.0 mm Stainless Steel	1
SL-010-03	Bard-Parker Scalpel Handle No. 3 Solid 12 cm - 4 3/4" Stainless Steel	2
SS-041-17	Mayo Dissecting Scissor Curved 17 cm - 6 3/4" Stainless Steel	1
SS-456-11	Strabismus Scissor Straight 11.5 cm - 4 1/2" Stainless Steel	1
SS-457-11	Strabismus Scissor Curved 11.5 cm - 4 1/2" Stainless Steel	1
SS-551-12	Universal Wire Cutting Scissor Angled - One Toothing Cutting Edge 12 cm - 4 3/4" Stainless Steel	1
SS-921-14	UltraCut TC Metzenbaum-Fine Dissecting Scissor - Slender Pattern Curved 14.5 cm - 5 3/4" Stainless Steel	1
SS-923-18	UltraCut TC Metzenbaum-Fine Dissecting Scissor - Slender Pattern Curved 18 cm - 7" Stainless Steel	1
FR-250-12	Adson Dissecting Forcep 1 x 2 Teeth 12 cm - 4 3/4" Stainless Steel	1
AF-160-12	Halsted-Mosquito Haemostatic Forcep Straight 12.5 cm - 5" Stainless Steel	6
AF-200-14	Kelly Haemostatic Forcep Straight 14.5 cm - 5 3/4" Stainless Steel	2
AF-773-08	Backhaus Towel Clamp 9 cm - 3 1/2" Stainless Steel	8
RT-162-16	Senn-Miller Retractor Sharp 16 cm - 6 1/4" Blade Size 20 x 7 mm Stainless Steel	1
RT-251-01	Langenbeck Retractor 22.5 cm - 8 3/4" Blade Size 30 x 11 mm Stainless Steel	1
RT-251-02	Langenbeck Retractor 22.5 cm - 8 3/4" Blade Size 30 x 14 mm Stainless Steel	1
RT-677-10	Jansen Self Retaining Retractor 3 x 3 Blunt Prongs - With Adjustment Screw 10 cm - 4" Stainless Steel	1
TR-101-08	Frazier Suction Tube 8 Charr. 19.5 cm - 7 3/4" Ø 2.6 mm Stainless Steel	1
TR-101-10	Frazier Suction Tube 10 Charr. 19.5 cm - 7 3/4" Ø 3.3 mm Stainless Steel	1
SU-588-12	UltraGripX-TC Webster Needle Holder Smooth Jaw 13 cm - 5" Stainless Steel	1
SU-586-13	UltraGripX TC Halsey Needle Holder 13 cm - 5" Stainless Steel	1
SI-012-15	Allis Intestinal and Tissue Grasping Forcep 4 x 5 Teeth 15.5 cm - 6" Stainless Steel	1
FR-772-01	Graefe Fixation Forcep Without Lock 11 cm - 4 1/4" Stainless Steel	1
OT-291-08	Mini-McGee Micro Alligator Forcep Smooth-Straight 8 cm - 3" Jaw Size 3.5 x 0.6 mm Stainless Steel	1
OT-293-08	Mini-McGee Micro Alligator Forcep Serrated-Straight 8 cm - 3" Jaw Size 3.5 x 0.6 mm Stainless Steel	1
OT-305-08	Micro Alligator Forcep Serrated-Bent Upwards 8 cm - 3" Jaw Size 4.0 x 0.6 mm Stainless Steel	1
OT-337-08	Bellucci Micro Scissor Left - Very Fine 8 cm - 3" Blade Size 4.0 x 0.8 mm Stainless Steel	1
OT-335-08	Bellucci Micro Scissor Right - Very Fine 8 cm - 3" Blade Size 4.0 x 0.8 mm Stainless Steel	1
OT-331-08	Bellucci Micro Scissor Straight - Very Fine 8 cm - 3" Blade Size 4.0 x 0.8 mm Stainless Steel	1
OT-365-01	House-Dieter Malleus Nipper Up Cutting 8 cm - 3" Jaw Opening 1.3 mm Stainless Steel	1
OT-405-15	Wullstein Sickle Knife 15.5 cm - 6" Stainless Steel	1
OT-421-04	Rosen Circular Cut Knife Fig. 4 15 cm - 6" Ø 2.0 mm Stainless Steel	1
OT-450-16	Barbara Micro Ear Needle Straight 16 cm - 6 1/4" Stainless Steel	1
OT-455-16	Rosen Micro Ear Pick 16 cm - 6 1/4" Stainless Steel	1
OT-520-01	House Micro Ear Curette 14.5 cm - 5 3/4" Cup Size 1 / Cup Size 2 1.0 mm - 1.2 mm Stainless Steel	1
OT-520-02	House Micro Ear Curette 14.5 cm - 5 3/4" Cup Size 1 / Cup Size 2 1.5 mm - 1.8 mm Stainless Steel	1
OT-520-03	House Micro Ear Curette 14.5 cm - 5 3/4" Cup Size 1 / Cup Size 2 1.0 mm - 1.2 mm Stainless Steel	1
OT-520-04	House Micro Ear Curette 14.5 cm - 5 3/4" Cup Size 1 / Cup Size 2 1.5 mm - 1.8 mm Stainless Steel	1
OT-052-01	Farrir Ear Speculum Fig. 1 - Black 3.6 cm / 1 1/2" Ø 5.5 mm Stainless Steel	1
OT-052-03	Farrir Ear Speculum Fig. 3 - Black 3.6 cm / 1 1/2" Ø 7.5 mm Stainless Steel	1
OT-052-05	Farrir Ear Speculum Fig. 5 - Black 3.6 cm / 1 1/2" Ø 9.5 mm Stainless Steel	1
OT-052-07	Farrir Ear Speculum Fig. 7 - Black 3.6 cm / 1 1/2" Ø 11.5 mm Stainless Steel	1

Middle Ear Set

Codes	Description	Quantity
OT-211-14	Gruenwald (Lucae) Ear Forcep Bayonet Shaped 14.5 cm - 5 3/4" Stainless Steel	1
OT-573-02	Buck Ear Curette Fig. 00 - Angled - Sharp 17 cm - 6 3/4" Ø 2.3 mm Stainless Steel	1
OT-573-11	Buck Ear Curette Fig. 1 - Angled - Sharp 17 cm - 6 3/4" Ø 3.0 mm Stainless Steel	1
OT-660-01	Billeau Ear Loop 16 cm - 6 1/4" Working End 2.0 mm Stainless Steel	1
OT-660-02	Billeau Ear Loop 16 cm - 6 1/4" Working End 3.0 mm Stainless Steel	1
OT-660-03	Billeau Ear Loop 16 cm - 6 1/4" Working End 4.0 mm Stainless Steel	1
RH-475-18	Freer Septum Elevator Double Ended - Sharp/Blunt 18.5 cm - 7 1/4" Blade Size 1 - Blade Size 2 5.0 mm - 5.0 mm Stainless Steel	1
RT-034-15	Joseph Mucus Hook Sharp 15.5 cm - 6" Stainless Steel	2
RT-074-07	Joseph Mucosa Hook 16.5 cm - 6 1/2" Width 7.0 mm Stainless Steel	1

Mastoid Surgery Set

Codes	Description	Quantity
SF-350-18	Foerster Sponge Holding Forcep Straight 18 cm - 7" Stainless Steel	2
AF-773-08	Backhaus Towel Clamp 9 cm - 3 1/2" Stainless Steel	4
SL-030-07	Scalpel Handle No. 7 Solid 16 cm - 6 1/4" Stainless Steel	2
FR-270-15	Gillies Dissecting Forcep 1 x 2 Teeth 15.5 cm - 6" Stainless Steel	2
FR-060-15	McIndoe Dressing Forcep 15 cm - 6" Stainless Steel	2
SS-041-14	Mayo Dissecting Scissor Curved 14.5 cm - 5 3/4" Stainless Steel	1
SS-440-11	Iris Scissor Straight 11.5 cm - 4 1/2" Stainless Steel	1
SS-430-10	Iris Scissor Straight 10 cm - 4" Stainless Steel	1
AF-183-14	Providence-Hospital Haemostatic Forcep Curved 14 cm - 5 1/2" Stainless Steel	6
SU-026-13	Halsey Needle Holder Serrated Jaws 13.5 cm - 5 1/4" Stainless Steel	2
OR-309-04	Partsch Bone Gouge 13.5 cm - 5 1/4" Blade Width 4 mm Stainless Steel	1
OR-309-05	Partsch Bone Gouge 13.5 cm - 5 1/4" Blade Width 5 mm Stainless Steel	1
OR-309-08	Partsch Bone Gouge 13.5 cm - 5 1/4" Blade Width 8 mm Stainless Steel	1
OR-302-10	Bone Gouge 13.5 cm - 5 1/4" Blade Width 10 mm Stainless Steel	1
OR-302-12	Bone Gouge 13.5 cm - 5 1/4" Blade Width 12 mm Stainless Steel	1
OT-201-13	Troeltsch (Wilde) Ear Forcep 12.5 cm - 5" Stainless Steel	1
OT-225-14	Hartmann Ear Polypus Forcep Regular Model 14 cm - 5 1/2" Stainless Steel	1
OT-251-08	Hartmann-Wullstein Alligator Forcep 8.5 cm - 3 1/4" Stainless Steel	1
RH-013-06	Thudichum Nasal Speculum Fig. 6 22 mm Stainless Steel	1
OR-476-16	Williger Bone Mallet 16.5 cm - 6 1/2" Head Vø 20.0 mm - Weight 140 Grams Stainless Steel	1
RT-012-01	Gillies Skin Hook Fig. 1 18 cm - 7" Stainless Steel	1
OR-563-16	Jansen Bone Curette 16 cm - 6 1/4" Stainless Steel	1
RT-735-01	Plester Self Retaining Retractor Left 10.5 cm - 4" Stainless Steel	1
RT-735-02	Plester Self Retaining Retractor Right 10.5 cm - 4" Stainless Steel	1
RH-547-21	Howarth Periosteal Elevator 21.5 cm - 8 1/2" Stainless Steel	1
RT-210-04	Volkmann Retractor 4 Sharp Prongs 22 cm - 8 3/4" Stainless Steel	2
SI-015-15	Allis Intestinal and Tissue Grasping Forcep 5 x 6 Teeth 15.5 cm - 6" Stainless Steel	2
FR-200-14	Standard Pattern Dissecting Forcep 1 x 2 Teeth 14.5 cm - 5 3/4" Stainless Steel	1
SS-004-14	Operating Scissor Straight - Sharp/Blunt 14.5 cm - 5 3/4" Stainless Steel	1
OR-558-01	Williger Bone Curette Double Ended - Oval/Oval - Fig. 0/Fig. 1 13.5 cm - 5 1/4" Stainless Steel	1
RT-732-13	Plester Self Retaining Retractor 2 x 2 Sharp Prongs 13 cm - 5" Stainless Steel	1
RH-475-18	Freer Septum Elevator Double Ended - Sharp/Blunt 18.5 cm - 7 1/4" Blade Size 1 - Blade Size 2 5.0 mm - 5.0 mm Stainless Steel	1
TR-101-10	Frazier Suction Tube 10 Charr. 19.5 cm - 7 3/4" Vø 3.3 mm Stainless Steel	1
TR-701-08	Redon Guide Needle 8 Charr. - Knife Tip 19.5 cm - 7 3/4" Tip Size 2.7 mm Stainless Steel	1
SU-060-14	Mayo-Hegar Needle Holder 14.5 cm - 5 3/4" Stainless Steel	1
SS-004-11	Operating Scissor Straight - Sharp/Blunt 12 cm - 4 3/4" Stainless Steel	1

Ear Microsurgery Set

Codes	Description	Quantity
RH-033-01	Hartmann-Halle Nasal Speculum Fig. 1 15.5 cm - 6" Blade Length 27 mm Stainless Steel	1
RH-451-05	McKenty Septum Elevator 15 cm - 6" Blade Size 5.0 mm Stainless Steel	1
OT-731-03	Lempert Periosteal Elevator 17 cm - 6 3/4" Width 3 mm Stainless Steel	1
FR-010-11	Standard Pattern Dressing Forcep 11.5 cm - 4 1/2" Stainless Steel	1
FR-012-14	Slender Pattern Dressing Forcep 14.5 cm - 5 3/4" Stainless Steel	1
SS-005-11	Operating Scissor Curved - Sharp/Blunt 12 cm - 4 3/4" Stainless Steel	1
RT-732-13	Plester Self Retaining Retractor 2 x 2 Sharp Prongs 13 cm - 5" Stainless Steel	1
RT-735-01	Plester Self Retaining Retractor Left 10.5 cm - 4" Stainless Steel	1
RT-735-02	Plester Self Retaining Retractor Right 10.5 cm - 4" Stainless Steel	1
RT-140-04	Wound Retractor 4 Sharp Prongs - Large Curve 16.5 cm - 6 1/2" Width 14.0 mm Stainless Steel	1
SL-010-03	Bard-Parker Scalpel Handle No. 3 Solid 12 cm - 4 3/4" Stainless Steel	2
SS-005-14	Operating Scissor Curved - Sharp/Blunt 14.5 cm - 5 3/4" Stainless Steel	1
AF-220-14	Rochester-Pean Haemostatic Forcep Straight 14.5 cm - 5 3/4" Stainless Steel	1
SU-586-13	UltraGripX TC Halsey Needle Holder 13 cm - 5" Stainless Steel	1
AF-160-12	Halsted-Mosquito Haemostatic Forcep Straight 12.5 cm - 5" Stainless Steel	1

Ear Microsurgery Set

Codes	Description	Quantity
AF-161-12	Halsted-Mosquito Haemostatic Forcep Curved 12 cm - 4 3/4" Stainless Steel	2
AF-773-11	Backhaus Towel Clamp 11.5 cm - 4 1/2" Stainless Steel	6
SS-551-12	Universal Wire Cutting Scissor Angled - One Toothing Cutting Edge 12 cm - 4 3/4" Stainless Steel	1
OT-012-02	Hartmann Ear Speculum Fig. 2 3.6 cm / 1 1/2" vø 5.0 mm Stainless Steel	1
OT-012-03	Hartmann Ear Speculum Fig. 3 3.6 cm / 1 1/2" vø 6.0 mm Stainless Steel	1
OT-012-04	Hartmann Ear Speculum Fig. 4 3.6 cm / 1 1/2" vø 7.0 mm Stainless Steel	1
OT-223-12	Hartmann Ear Polypus Forcep Fine Model 12 cm - 4 3/4" Stainless Steel	1
OT-255-08	Hartmann Alligator Forcep 8.5 cm - 3 1/4" Stainless Steel	1
OT-251-08	Hartmann-Wullstein Alligator Forcep 8.5 cm - 3 1/4" Stainless Steel	1
OT-253-08	Hartmann-Wullstein Alligator Forcep Cup Shaped 8.5 cm - 3 1/4" vø 1.0 mm Stainless Steel	1
OT-261-20	Hartmann Alligator Forcep Cup Shaped 8.5 cm - 3 1/4" vø 2.0 mm Stainless Steel	1
OT-331-08	Bellucci Micro Scissor Straight - Very Fine 8 cm - 3" Blade Size 4.0 x 0.8 mm Stainless Steel	1
OT-365-01	House-Dieter Malleus Nipper Up Cutting 8 cm - 3" Jaw Opening 1.3 mm Stainless Steel	1
OT-361-01	McGee Wire Bending Forcep Bent Downwards 8 cm - 3" Jaw Size 3.5 x 0.8 mm Stainless Steel	1
TR-185-07	Verhoeven Suction Cannula With Luer Cone 10 cm - 4" Diameter 0.7 mm vø Stainless Steel	4
TR-185-10	Verhoeven Suction Cannula With Luer Cone 10 cm - 4" Diameter 1.0 mm vø Stainless Steel	4
TR-185-12	Verhoeven Suction Cannula With Luer Cone 10 cm - 4" Diameter 1.2 mm vø Stainless Steel	4
TR-185-15	Verhoeven Suction Cannula With Luer Cone 10 cm - 4" vø 1.5 mm Stainless Steel	4
TR-180-00	House Suction Adapter With Finger Cut Off Stainless Steel	1
TR-101-06	Frazier Suction Tube 6 Charr. 19.5 cm - 7 3/4" vø 2.0 mm Stainless Steel	1
TR-101-07	Frazier Suction Tube 7 Charr. 19.5 cm - 7 3/4" vø 2.3 mm Stainless Steel	1
TR-101-09	Frazier Suction Tube 9 Charr. 19.5 cm - 7 3/4" vø 3.0 mm Stainless Steel	1
OT-673-03	Lucae Ear Hook Large 14 cm - 5 1/2" Stainless Steel	1
OT-572-02	Buck Ear Curette Fig. 00 - Straight - Sharp 17 cm - 6 3/4" vø 2.3 mm Stainless Steel	1
OT-572-10	Buck Ear Curette Fig. 0 - Straight - Sharp 17 cm - 6 3/4" vø 2.6 mm Stainless Steel	1
OT-473-10	McGee Micro Ear Needle For Footplate 15.5 cm - 6" Tip Size 1.0 mm Stainless Steel	1
OT-520-01	House Micro Ear Curette 14.5 cm - 5 3/4" Cup Size 1 / Cup Size 2 1.0 mm - 1.2 mm Stainless Steel	1
OT-509-05	Rosen Micro Dissector 15.5 cm - 6" Stainless Steel	1
OT-461-06	Barbara Micro Ear Needle Angled 45° 16 cm - 6 1/4" Tip Size 0.6 mm Stainless Steel	1
OT-405-15	Wullstein Sickle Knife 15.5 cm - 6" Stainless Steel	1
OT-459-03	Barbara Micro Ear Needle Angled 25° 16 cm - 6 1/4" Tip Size 0.3 mm Stainless Steel	1
OT-459-06	Barbara Micro Ear Needle Angled 25° 16 cm - 6 1/4" Tip Size 0.6 mm Stainless Steel	1
OT-535-18	Zoellner Application Instrument Stainless Steel	1
OT-461-10	Barbara Micro Ear Needle Angled 45° 16 cm - 6 1/4" Tip Size 1.0 mm Stainless Steel	1
OT-409-15	Plester Flap Knife 15 cm - 6" Stainless Steel	1
OT-421-04	Rosen Circular Cut Knife Fig. 4 15 cm - 6" vø 2.0 mm Stainless Steel	1
OT-415-15	Plester Sickle Knife 15 cm - 6" Stainless Steel	1
OT-453-16	Barbara Micro Ear Pick 16 cm - 6 1/4" Stainless Steel	1
OT-463-03	Barbara Micro Ear Needle Angled 90° 16 cm - 6 1/4" Tip Size 0.3 mm Stainless Steel	1
OT-463-06	Barbara Micro Ear Needle Angled 90° 16 cm - 6 1/4" Tip Size 0.6 mm Stainless Steel	1
OT-463-10	Barbara Micro Ear Needle Angled 90° 16 cm - 6 1/4" Tip Size 1.0 mm Stainless Steel	1

Ear Foreign Body Removal Set

Codes	Description	Quantity
OT-211-14	Gruenwald (Lucae) Ear Forcep Bayonet Shaped 14.5 cm - 5 3/4" Stainless Steel	1
OT-237-10	Quire Foreign Body Lever Forcep 10.5 cm / 4" Stainless Steel	1
OT-255-08	Hartmann Alligator Forcep 8.5 cm - 3 1/4" Stainless Steel	1
RT-012-02	Gilles Skin Hook Fig. 2 18 cm - 7" Stainless Steel	1
OT-660-01	Billeau Ear Loop 16 cm - 6 1/4" Working End 2.0 mm Stainless Steel	1
OT-660-02	Billeau Ear Loop 16 cm - 6 1/4" Working End 3.0 mm Stainless Steel	1
OT-010-03	Hartmann Ear Specula Set of 4 Ref:- OT-012-01 to OT-012-04 Stainless Steel	1

Ear Polyp Set

Codes	Description	Quantity
OT-223-12	Hartmann Ear Polypus Forcep Fine Model 12 cm - 4 3/4" Stainless Steel	1
OT-161-18	Krause Polypus Snare Open Tip 18.5 cm - 7 1/4" Stainless Steel	1
OT-170-30	Snare Wire 100 cm - 39 1/2" Diameter 0.3 mm Stainless Steel	1
TR-111-05	Fergusson Suction Tube 5 Charr. 18 cm - 7" vø 1.7 mm Stainless Steel	1
OT-201-11	Troeltsch (Wilde) Ear Forcep 11 cm / 4 3/8" Stainless Steel	1

OTO-RHINO-Laryngological Surgery Set

Codes	Description	Quantity
SI-012-15	Allis Intestinal and Tissue Grasping Forcep 4 x 5 Teeth 15.5 cm - 6" Stainless Steel	1
TR-151-15	Baron-Schuknecht Suction Tube With Finger Cutt Off Working Length 75 mm - Ø 1.5 mm Stainless Steel	1
TR-151-20	Baron-Schuknecht Suction Tube With Finger Cutt Off Working Length 75 mm - Ø 2.0 mm Stainless Steel	1
OT-660-01	Billeau Ear Loop 16 cm - 6 1/4" Working End 2.0 mm Stainless Steel	1
OT-660-02	Billeau Ear Loop 16 cm - 6 1/4" Working End 3.0 mm Stainless Steel	1
OT-660-03	Billeau Ear Loop 16 cm - 6 1/4" Working End 4.0 mm Stainless Steel	1
OT-570-10	Buck Ear Curette Fig. 0 - Straight - Blunt 17 cm - 6 3/4" Ø 2.6 mm Stainless Steel	1
OT-570-11	Buck Ear Curette Fig. 1 - Straight - Blunt 17 cm - 6 3/4" Ø 3.0 mm Stainless Steel	1
OT-570-12	Buck Ear Curette Fig. 2 - Straight - Blunt 17 cm - 6 3/4" Ø 3.7 mm Stainless Steel	1
OT-572-10	Buck Ear Curette Fig. 0 - Straight - Sharp 17 cm - 6 3/4" Ø 2.6 mm Stainless Steel	1
OT-572-11	Buck Ear Curette Fig. 1 - Straight - Sharp 17 cm - 6 3/4" Ø 3.0 mm Stainless Steel	1
OT-572-12	Buck Ear Curette Fig. 2 - Straight - Sharp 17 cm - 6 3/4" Ø 3.7 mm Stainless Steel	1
OT-677-01	Day Ear Hook Long 17 cm - 6 3/4" Stainless Steel	1
OT-677-02	Day Ear Hook Short 17 cm - 6 3/4" Stainless Steel	1
OT-255-14	Hartmann Alligator Forcep 14 cm - 5 1/2" Stainless Steel	1
OT-261-25	Hartmann Alligator Forcep Cup Shaped 8.5 cm - 3 1/4" Ø 2.5 mm Stainless Steel	1
SS-440-11	Iris Scissor Straight 11.5 cm - 4 1/2" Stainless Steel	1
SS-441-11	Iris Scissor Curved 11.5 cm - 4 1/2" Stainless Steel	1
OT-211-14	Gruenwald (Lucae) Ear Forcep Bayonet Shaped 14.5 cm - 5 3/4" Stainless Steel	1
OT-283-08	McGee Micro Alligator Forcep Serrated-Straight 8 cm - 3" Jaw Size 4.0 x 0.8 mm Stainless Steel	1
OT-293-08	Mini-McGee Micro Alligator Forcep Serrated-Straight 8 cm - 3" Jaw Size 3.5 x 0.6 mm Stainless Steel	1
OT-311-08	Micro Alligator Forcep Straight - Cup Shaped 8 cm - 3" Cup Size - Jaw Size 1.0 x 0.9 mm - 4.0 mm Stainless Steel	1
OT-255-08	Hartmann Alligator Forcep 8.5 cm - 3 1/4" Stainless Steel	1
OT-555-01	Shapleigh Ear Curette Fig. 1 15 cm - 6" Stainless Steel	1
OT-555-02	Shapleigh Ear Curette Fig. 2 15 cm - 6" Stainless Steel	1
SS-456-11	Strabismus Scissor Straight 11.5 cm - 4 1/2" Stainless Steel	1
SS-457-11	Strabismus Scissor Curved 11.5 cm - 4 1/2" Stainless Steel	1
RH-023-01	Mod. Vienna Nasal Speculum Fig. 1 Blade Length 27 mm Stainless Steel	1
RH-023-02	Mod. Vienna Nasal Speculum Fig. 2 Blade Length 31 mm Stainless Steel	1
RH-023-03	Mod. Vienna Nasal Speculum Fig. 3 Blade Length 35 mm Stainless Steel	1
OT-201-11	Troeltsch (Wilde) Ear Forcep 11 cm / 4 3/8" Stainless Steel	1
OT-225-14	Hartmann Ear Polypus Forcep Regular Model 14 cm - 5 1/2" Stainless Steel	1
OT-229-16	Tilley Ear Polypus Forcep 16 cm - 6 1/4" Stainless Steel	1
OT-800-00	Ear Syringe Complete With Protection Disk and 2 Tips Capacity 50 ml Stainless Steel	1
OT-751-02	Citelli Mastoid Rongeur 17 cm - 6 3/4" Bite Size 2 mm Stainless Steel	1
RH-804-01	Cottle-Walsham Septum Straightning Forcep Straight 23 cm - 9" Stainless Steel	1
AF-773-08	Backhaus Towel Clamp 9 cm - 3 1/2" Stainless Steel	6
SL-010-03	Bard-Parker Scalpel Handle No. 3 Solid 12 cm - 4 3/4" Stainless Steel	2
RT-251-01	Langenbeck Retractor 22.5 cm - 8 3/4" Blade Size 30 x 11 mm Stainless Steel	2
RT-220-01	Kocher Retractor 1 Sharp Prong 22 cm - 8 3/4" Stainless Steel	2
RT-210-04	Volkmann Retractor 4 Sharp Prongs 22 cm - 8 3/4" Stainless Steel	2
TM-181-21	Schoenborn Tracheal Retractor 21 cm - 8 1/4" Stainless Steel	1
AF-161-12	Halsted-Mosquito Haemostatic Forcep Curved 12 cm - 4 3/4" Stainless Steel	12
AF-220-24	Rochester-Pean Haemostatic Forcep Straight 24.5 cm - 9 3/4" Stainless Steel	4
AF-161-12	Halsted-Mosquito Haemostatic Forcep Curved 12 cm - 4 3/4" Stainless Steel	6
FR-012-14	Slender Pattern Dressing Forcep 14.5 cm - 5 3/4" Stainless Steel	1
FR-210-14	Slender Pattern Dissecting Forcep 1 x 2 Teeth 14.5 cm - 5 3/4" Stainless Steel	1
SS-005-14	Operating Scissor Curved - Sharp/Blunt 14.5 cm - 5 3/4" Stainless Steel	1
SS-105-15	Joseph Dissecting Scissor / Opreating Scissor Curved 15 cm - 6" Stainless Steel	1
SS-131-14	Metzenbaum Dissecting Scissor Curved - Blunt/Blunt 14.5 cm - 5 3/4" Stainless Steel	1
SS-041-17	Mayo Dissecting Scissor Curved 17 cm - 6 3/4" Stainless Steel	1
RH-475-18	Freer Septum Elevator Double Ended - Sharp/Blunt 18.5 cm - 7 1/4" Blade Size 1 - Blade Size 2 5.0 mm - 5.0 mm Stainless Steel	1
OR-949-18	Beyer Bone Rongeur Compound Action 18 cm - 7" Stainless Steel	1
SS-341-21	Killian Bone Shear 21 cm - 8 1/4" Stainless Steel	1
TM-252-20	Lahey Goitre Grasping Forcep 3 x 3 Teeth 20 cm - 8" Stainless Steel	1
SU-855-20	Kocher Ligature Needle 20.5 cm - 8" Stainless Steel	2
SU-877-20	Deschamps Ligature Needle Blunt for Left Hand Medium 20 cm - 8" Stainless Steel	2
TR-251-27	Yankauer Suction Tube Complete Central Bore 2.0 mm 31 cm - 12 1/4" Stainless Steel	1
SU-050-15	Crile-Wood Needl Holder 15 cm - 6" Stainless Steel	1
TM-321-06	Chevalier-Jackson Tracheal Tube Fig. 2 55 mm Ø 6.0 mm Stainless Steel	1
TM-321-08	Chevalier-Jackson Tracheal Tube Fig. 4 65 mm Ø 8.0 mm Stainless Steel	1
TM-321-10	Chevalier-Jackson Tracheal Tube Fig. 6 75 mm Ø 10.0 mm Stainless Steel	1
RT-101-01	Wound Retractor 1 Blunt Prong 16.5 cm - 6 1/2" Stainless Steel	2
TR-701-12	Redon Guide Needle 12 Charr. - Knife Tip 19.5 cm - 7 3/4" Tip Size 4.0 mm Stainless Steel	1
SF-350-20	Foerster Sponge Holding Forcep Straight 20.5 cm - 8" Stainless Steel	4
SI-015-19	Allis Intestinal and Tissue Grasping Forcep 5 x 6 Teeth 19 cm - 7 1/2" Stainless Steel	2
RT-172-18	Cushing Retractor / Saddle Hook 20.5 cm - 8" Blade Size 18 mm Stainless Steel	2
RT-561-00	Roux Retractor Set of Fig. 1, Fig. 2 and Fig. 3 Stainless Steel	1
RT-301-01	Zenker Retractor 25 cm - 9 3/4" Blade Size 35 x 8 mm Stainless Steel	1
RT-301-02	Zenker Retractor 26.5 cm - 10 1/2" Blade Size 60 x 12 mm Stainless Steel	1
AF-421-02	Overholt-Geissendorfer Dissecting and Ligature Forcep Fig. 2 20 cm - 8" Stainless Steel	1

OTO-RHINO-Laryngological Surgery Set

Codes	Description	Quantity
AF-421-03	Overholt-Geissendorfer Dissecting and Ligature Forcep Fig. 3 21.5 cm - 8 1/2" Stainless Steel	1
AF-421-04	Overholt-Geissendorfer Dissecting and Ligature Forcep Fig. 4 22 cm - 8 3/4" Stainless Steel	1
AF-421-07	Overholt-Geissendorfer Dissecting and Ligature Forcep Fig. 7 22 cm - 8 3/4" Stainless Steel	1
AF-269-18	Adson Haemostatic Forcep Curved 18.5 cm - 7 1/4" Stainless Steel	1
AF-471-18	Gemini Dissecting and Ligature Forcep Curved 18 cm - 7" Stainless Steel	1
AF-471-23	Gemini Dissecting and Ligature Forcep Curved 23 cm - 9" Stainless Steel	1
SS-004-16	Operating Scissor Straight - Sharp/Blunt 16.5 cm - 6 1/2" Stainless Steel	1
SS-002-18	Operating Scissor Straight - Blunt/Blunt 18.5 cm - 7 1/4" Stainless Steel	1
SS-457-11	Strabismus Scissor Curved 11.5 cm - 4 1/2" Stainless Steel	1
SU-060-18	Mayo-Hegar Needle Holder 18.5 cm - 7 1/4" Stainless Steel	2
FR-050-12	Adson Dressing Forcep 12 cm - 4 3/4" Stainless Steel	1
FR-250-12	Adson Dissecting Forcep 1 x 2 Teeth 12 cm - 4 3/4" Stainless Steel	1
SI-015-15	Allis Intestinal and Tissue Grasping Forcep 5 x 6 Teeth 15.5 cm - 6" Stainless Steel	2
AF-773-08	Backhaus Towel Clamp 9 cm - 3 1/2" Stainless Steel	4
AF-211-14	Crile Haemostatic Forcep Curved 14 cm - 5 1/2" Stainless Steel	4
SF-351-18	Foerster Sponge Holding Forcep Curved 17.5 cm - 7" Stainless Steel	1
AF-160-12	Halsted-Mosquito Haemostatic Forcep Straight 12.5 cm - 5" Stainless Steel	3
AF-161-12	Halsted-Mosquito Haemostatic Forcep Curved 12 cm - 4 3/4" Stainless Steel	3
RT-040-02	Cushing Nerve Hook Blunt - Fig.2 19.5 cm - 7 3/4" Stainless Steel	2
TM-151-14	Laborde Tracheal Dilator 14 cm - 5 1/2" Stainless Steel	1
SS-040-14	Mayo Dissecting Scissor Straight 14.5 cm - 5 3/4" Stainless Steel	1
SS-131-14	Metzenbaum Dissecting Scissor Curved - Blunt/Blunt 14.5 cm - 5 3/4" Stainless Steel	1
TM-171-16	Iterson Tracheal Hook Sharp 17 cm - 6 3/4" Stainless Steel	1
TR-251-27	Yankauer Suction Tube Complete Central Bore 2.0 mm 31 cm - 12 1/4" Stainless Steel	1
AF-221-14	Rochester-Pean Haemostatic Forcep Curved 14 cm - 5 1/2" Stainless Steel	2
SL-010-03	Bard-Parker Scalpel Handle No. 3 Solid 12 cm - 4 3/4" Stainless Steel	1
SS-542-13	Spencer Stich Scissor 13 cm - 5" Stainless Steel	1
SU-594-15	UltraGripX TC Crile-Wood Needle Holder 15 cm - 6" Stainless Steel	1
RT-710-11	Weitlaner Self Retaining Retractor 2 x 3 Sharp Prongs 11 cm - 4 1/4" Stainless Steel	1

Sinus Surgery Set

Codes	Description	Quantity
OR-603-18	Molt Periosteal Raspatory / Elevator 18 cm - 7" Stainless Steel	1
OR-941-01	Blumenthal Bone Rongeur 15 cm - 6" Stainless Steel	1
OR-440-05	Bone Tamper 15.5 cm - 6" vø 5.0 mm Stainless Steel	1
DI-850-20	Castroviejo Caliper Straight 8 cm - 3" Measuring Range 20 mm Stainless Steel	1
RH-703-01	Miller-Colburn Bone File Fig. 1 - Cross Cut 18 cm - 7" Stainless Steel	1
RH-703-03	Miller-Colburn Bone File Fig. 3 - Cross Cut 18 cm - 7" Stainless Steel	1
RH-620-04	Cottle Osteotome Straight 18.5 cm - 7 1/4" Blade Width 4.0 mm Stainless Steel	1
RH-620-07	Cottle Osteotome Straight 18.5 cm - 7 1/4" Blade Width 7.0 mm Stainless Steel	1
RH-620-09	Cottle Osteotome Straight 18.5 cm - 7 1/4" Blade Width 9.0 mm Stainless Steel	1
SS-101-16	Kelly Dissecting Scissor / Opreating Scissor Curved 16 cm - 6 1/4" Stainless Steel	1
RH-281-19	Jansen-Middleton Septum Forcep Through Cutting Jaws 19 cm - 7 1/2" Stainless Steel	1
SS-417-16	Locklin Gum Scissor Angled - One Toothing Cutting Edge 16 cm - 6 1/4" Stainless Steel	1
OR-460-17	Gerzog Bone Mallet 19 cm - 7 1/2" Head vø 25.0 mm - Weight 225 Grams Stainless Steel	1
OR-975-15	Cleveland Bone Cutting Forcep 14.5 cm - 5 3/4" Stainless Steel	1
OM-121-11	Martin-Simplex Lip and Cheek Retractor 10.5 cm / 4" Stainless Steel	1
OM-451-01	Tatum Sinus Lift Curette Fig. 1 18 cm - 7" Blade Width 2.6 mm - 2.6 mm Stainless Steel	1
OM-453-02	Tatum Sinus Lift Curette Fig. 2 17.5 cm - 7" Blade Width 2.6 mm - 2.6 mm Stainless Steel	1
OM-455-03	Tatum Sinus Lift Curette Fig. 3 18 cm - 7" Blade Width 4.1 mm - 4.1 mm Stainless Steel	1
OM-457-04	Tatum Sinus Lift Curette Fig. 4 18.5 cm - 7 1/4" Blade Width 4.1 mm - 4.1 mm Stainless Steel	1
OM-459-05	Tatum Sinus Lift Curette Fig. 5 18.5 cm - 7 1/4" Blade Width 6.7 mm - 6.7 mm Stainless Steel	1
OM-461-06	Tatum Sinus Lift Curette Fig. 6 19 cm - 7 1/2" Blade Width 6.7 mm - 6.7 mm Stainless Steel	1
OM-463-07	Tatum Sinus Lift Curette Fig. 7 16 cm - 6 1/4" Blade Width 2.8 mm - 6.2 mm Stainless Steel	1
SU-594-15	UltraGripX TC Crile-Wood Needle Holder 15 cm - 6" Stainless Steel	1
TR-111-09	Fergusson Suction Tube 9 Charr. 18 cm - 7" vø 3.0 mm Stainless Steel	1

FESS (Functional Endoscopic Sinus Surgery) SET

Codes	Description	Quantity
TR-101-09	Frazier Suction Tube 9 Charr. 19.5 cm - 7 3/4" √ø 3.0 mm Stainless Steel	1
PR-500-14	Farrell Cotton Applicator Round Spiral End 14 cm - 5 1/2" Working End √ø 0.9 mm Stainless Steel	1
SF-350-18	Foerster Sponge Holding Forcep Straight 18 cm - 7" Stainless Steel	1
RH-475-18	Freer Septum Elevator Double Ended - Sharp/Blunt 18.5 cm - 7 1/4" Blade Size 1 - Blade Size 2 5.0 mm - 5.0 mm Stainless Steel	1
RH-735-18	Wiener Nasal Rasp 18 cm - 7" Stainless Steel	1
RH-841-01	Coakley Antrum Curette 17 cm - 6 3/4" Stainless Steel	1
RH-841-04	Coakley Antrum Curette 17 cm - 6 3/4" Stainless Steel	1
RH-841-05	Coakley Antrum Curette 17 cm - 6 3/4" Stainless Steel	1
TL-015-19	Abraham Tonsil Knife 19 cm - 7 1/2" Stainless Steel	1
RH-645-01	Antrum Curette Oval Fig. 1 19 cm - 7 1/2" Stainless Steel	1
RH-645-02	Antrum Curette Oval Fig. 2 19 cm - 7 1/2" Stainless Steel	1
RH-645-03	Antrum Curette Oval Fig. 3 19 cm - 7 1/2" Stainless Steel	1
RH-647-02	Antrum Curette Oblong Small 19 cm - 7 1/2" Stainless Steel	1
RH-649-01	Antrum Curette Angled Forward Cutting 19 cm - 7 1/2" Tip Size 7 x 10 mm Stainless Steel	1
RH-649-02	Antrum Curette Angled Backward Cutting 19 cm - 7 1/2" Tip Size 7 x 10 mm Stainless Steel	1
RH-541-22	Cottle Periosteal Elevator 23 cm - 9" Blade Size 1 - Blade Size 2 3.3 mm - 3.7 mm Stainless Steel	1
RH-033-01	Hartmann-Halle Nasal Speculum Fig. 1 15.5 cm - 6" Blade Length 27 mm Stainless Steel	1
RH-033-02	Hartmann-Halle Nasal Speculum Fig. 2 15.5 cm - 6" Blade Length 31 mm Stainless Steel	1
TR-010-00	Universal Trocar With 4 Trocar Tips Diameter 1.7 mm, 3.0 mm, 4.5 mm and 6.0 mm 14.5 cm - 5 3/4" Stainless Steel	1
RH-881-30	Eicken-Killian Irrigation Cannula 14 cm - 5 1/2" √ø 3.0 mm Stainless Steel	1
RH-881-40	Eicken-Killian Irrigation Cannula 14 cm - 5 1/2" √ø 4.0 mm Stainless Steel	1
TR-493-20	Tonsil Needle Tapered Tip - With Luer Lock Connection Needle √ø 0.65 mm Stainless Steel	1
OT-266-08	Struempler Alligator Forcep Angled 45° Upwards 8.5 cm - 3 1/4" Stainless Steel	1
RH-181-12	Takahashi Ethmoid Forcep Straight 12 cm - 4 3/4" Stainless Steel	1
RH-105-21	Gruenwald (Jansen) Nasal Tampon Forcep 21.5 cm - 8 1/2" Stainless Steel	1
RH-180-21	Sinus Forcep Cup Shaped Jaws Straight 12 cm - 4 3/4" Bite Size Ø 3 mm Stainless Steel	1
RH-231-02	Struycken Nasal Cutting Forcep Fig. 2 12 cm - 4 3/4" Bite Size 1.2 mm Stainless Steel	1
RH-309-01	Antrum Punch Backward Cutting 14 cm - 5 1/2" Bite size 4 mm x 5 mm Stainless Steel	1
RH-309-02	Antrum Punch Forward Cutting 14 cm - 5 1/2" Bite size 4 mm x 5 mm Stainless Steel	1
RH-201-01	Weil-Blakesley Nasal Cutting Forcep Straight - Fig. 1 12 cm - 4 3/4" Bite Size 3.0 mm Stainless Steel	1
RH-201-02	Weil-Blakesley Nasal Cutting Forcep Straight - Fig. 2 12 cm - 4 3/4" Bite Size 3.5 mm Stainless Steel	1
RH-203-01	Weil-Blakesley Nasal Cutting Forcep Angled 45° - Fig. 1 12 cm - 4 3/4" Bite Size 3.0 mm Stainless Steel	1
RH-203-02	Weil-Blakesley Nasal Cutting Forcep Angled 45° - Fig. 2 12 cm - 4 3/4" Bite Size 3.5 mm Stainless Steel	1
RH-205-02	Weil-Blakesley Nasal Cutting Forcep Angled 90° - Fig. 2 12 cm - 4 3/4" Bite Size 3.5 mm Stainless Steel	1
RH-251-25	Blakesley Nasal Cutting Forcep Through Cut Straight Bite Size 2.5 mm Shaft Length 12 cm - 4 3/4" Stainless Steel	1
RH-251-35	Blakesley Nasal Cutting Forcep Through Cut Straight Bite Size 3.5 mm Shaft Length 12 cm - 4 3/4" Stainless Steel	1
RH-257-35	Blakesley Nasal Cutting Forcep Through Cut Angled 45° Upwards Bite Size 3.5 mm Shaft Length 12 cm - 4 3/4" Stainless Steel	1
RH-234-00	Stammberger Nasal Cutting Forcep Retrograde 360° Rotatable 12 cm - 4 3/4" Bite Size 2.0 mm Stainless Steel	1
RH-281-19	Jansen-Middleton Septum Forcep Through Cutting Jaws 19 cm - 7 1/2" Stainless Steel	1
SI-012-15	Allis Intestinal and Tissue Grasping Forcep 4 x 5 Teeth 15.5 cm - 6" Stainless Steel	1
FR-250-12	Adson Dissecting Forcep 1 x 2 Teeth 12 cm - 4 3/4" Stainless Steel	1
AF-161-12	Halsted-Mosquito Haemostatic Forcep Curved 12 cm - 4 3/4" Stainless Steel	2
AF-160-12	Halsted-Mosquito Haemostatic Forcep Straight 12.5 cm - 5" Stainless Steel	2
AF-211-14	Crile Haemostatic Forcep Curved 14 cm - 5 1/2" Stainless Steel	2
AF-210-14	Crile Haemostatic Forcep Straight 14.5 cm - 5 3/4" Stainless Steel	2
AF-773-08	Backhaus Towel Clamp 9 cm - 3 1/2" Stainless Steel	4
SS-361-18	Heymann Nasal Scissor 18.5 cm - 7 1/4" Stainless Steel	1
RH-261-11	Nasal Sinus Scissor Straight 12 cm - 4 3/4" Stainless Steel	1
RH-261-13	Nasal Sinus Scissor Curved Right 12 cm - 4 3/4" Stainless Steel	1
RH-261-14	Nasal Sinus Scissor Curved Left 12 cm - 4 3/4" Stainless Steel	1
SS-440-11 SC	SuperEdge Iris Scissor Straight 11.5 cm - 4 1/2" Stainless Steel	1
SS-441-11 SC	SuperEdge-ø Iris Scissor Curved 11.5 cm - 4 1/2" Stainless Steel	1
SU-594-15	UltraGripX TC Crile-Wood Needle Holder 15 cm - 6" Stainless Steel	1
SU-584-13	UltraGripX-ø TC Halsey Needle Holder Smooth Jaw 13 cm - 5" Stainless Steel	1
RH-271-15	Cottle Masing Needle Holder Smooth Jaws Shaft Size 15 cm - 6" Stainless Steel	1
OT-683-02	Kuhn Frontal Sinus Seeker Double Ended Fig. 2, Both Sides Curved 90° One Tip Straight Other Tip Reverse Angle 22 cm - 8 1/2" Stainless Steel	1
OT-681-19	Maxillary Ostium Sinus Seeker With Ball Probes Double Ended (√ø 1.2 and 2.0 mm) 19 cm - 7 1/2" Stainless Steel	1
SL-010-03	Bard-Parker Scalpel Handle No. 3 Solid 12 cm - 4 3/4" Stainless Steel	1

Caldwell Luc Set

Codes	Description	Quantity
OR-470-19	Cottle Bone Mallet Plain - Convex 19 cm - 7 1/2" Head Ø 30.0 mm - Weight 300 Grams Stainless Steel	1
SL-030-07	Scalpel Handle No. 7 Solid 16 cm - 6 1/4" Stainless Steel	1
AF-773-13	Backhaus Towel Clamp 13 cm - 5" Stainless Steel	2
SF-350-18	Foerster Sponge Holding Forcep Straight 18 cm - 7" Stainless Steel	1
SU-588-12	UltraGripX-ö TC Webster Needle Holder Smooth Jaw 13 cm - 5" Stainless Steel	1
NS-406-01	Love-Gruenwald Laminectomy Rongeur Straight 18 cm - 7" Bite Size 3 x 10 mm Stainless Steel	1
SI-012-15	Allis Intestinal and Tissue Grasping Forcep 4 x 5 Teeth 15.5 cm - 6" Stainless Steel	1
RH-201-01	Weil-Blakesley Nasal Cutting Forcep Straight - Fig. 1 12 cm - 4 3/4" Bite Size 3.0 mm Stainless Steel	1
RH-201-03	Weil-Blakesley Nasal Cutting Forcep Straight - Fig. 3 12 cm - 4 3/4" Bite Size 4.5 mm Stainless Steel	1
RH-201-05	Weil-Blakesley Nasal Cutting Forcep Straight - Fig. 5 12 cm - 4 3/4" Bite Size 5.5 mm Stainless Steel	1
RH-303-01	Hajek-Kofler Sphenoid Punch Up Cutting 14 cm - 5 1/2" Bite size 3.5 x 3.5 mm Stainless Steel	1
RH-305-05	Kerrison Sphenoid Punch Up Cutting 9 cm - 3 1/2" Bite size 5.0 x 5.0 mm Stainless Steel	1
RH-461-16	Joseph Septum Elevator 16 cm - 6 1/4" Blade Size 4.0 mm Stainless Steel	1
RH-600-16	Freer Septum Chisel Straight 16.5 cm - 6 1/2" Blade Width 4.0 mm Stainless Steel	1
RH-610-02	Ballenger Septum Gouge 15 cm - 6" Blade Width 2.0 mm Stainless Steel	1
RH-610-04	Ballenger Septum Gouge 15 cm - 6" Blade Width 4.0 mm Stainless Steel	1
RH-610-06	Ballenger Septum Gouge 15 cm - 6" Blade Width 6.0 mm Stainless Steel	1
RH-841-01	Coakley Antrum Curette 17 cm - 6 3/4" Stainless Steel	1
RH-841-02	Coakley Antrum Curette 17 cm - 6 3/4" Stainless Steel	1
RH-841-03	Coakley Antrum Curette 17 cm - 6 3/4" Stainless Steel	1
RH-841-04	Coakley Antrum Curette 17 cm - 6 3/4" Stainless Steel	1
RH-841-05	Coakley Antrum Curette 17 cm - 6 3/4" Stainless Steel	1
RH-841-06	Coakley Antrum Curette 17 cm - 6 3/4" Stainless Steel	1

Thyroid Surgery Set

Codes	Description	Quantity
RT-162-16	Senn-Miller Retractor Sharp 16 cm - 6 1/4" Blade Size 20 x 7 mm Stainless Steel	2
RT-353-01	Middeldorf Retractor 21.5 cm - 8 1/2" Blade Size 15 x 17 mm Stainless Steel	2
RT-357-22	Greene Retractor 22.5 cm - 8 3/4" Blade Size 22 x 22 mm Stainless Steel	2
RT-174-12	Cushing Retractor / Saddle Hook 24 cm - 9 1/2" Blade Size 12 mm Stainless Steel	2
RT-805-13	Gelpi Self Retaining Retractor 13 cm - 5" Stainless Steel	2
RT-712-13	Weitlaner Self Retaining Retractor 3 x 4 Sharp Prongs 13.5 cm - 5 1/4" Stainless Steel	1
RT-271-19	Lahey Retractor 19.5 cm - 7 3/4" Blade Size 31 x 6 mm Stainless Steel	2
FR-070-12	Semken Dressing Forcep 12.5 cm - 5" Stainless Steel	2
FR-050-12	Adson Dressing Forcep 12 cm - 4 3/4" Stainless Steel	1
FR-590-16	DeBakey Atrauma Forcep 16 cm - 6 1/4" Tip Size 1.5 mm Stainless Steel	2
FR-206-14	Dissecting Forcep 4 x 5 Teeth 14.5 cm - 5 3/4" Stainless Steel	2
FR-252-12	Mini-Adson Dissecting Forcep 1 x 2 Teeth 12 cm - 4 3/4" Stainless Steel	2
SS-861-14	UltraCut TC Mayo Dissecting Scissor Curved 14.5 cm - 5 3/4" Stainless Steel	1
SS-861-17	UltraCut TC Mayo Dissecting Scissor Curved 17 cm - 6 3/4" Stainless Steel	1
SS-905-18	UltraCut TC Metzenbaum Dissecting Scissor Curved 18 cm - 7" Stainless Steel	2
SS-852-14	UltraCut-ö TC Operating Scissor Straight - Sharp/Blunt 14.5 cm - 5 3/4" Stainless Steel	1
AF-161-12	Halsted-Mosquito Haemostatic Forcep Curved 12 cm - 4 3/4" Stainless Steel	12
AF-211-14	Crile Haemostatic Forcep Curved 14 cm - 5 1/2" Stainless Steel	6
AF-221-16	Rochester-Pean Haemostatic Forcep Curved 16 cm - 6 1/4" Stainless Steel	2
SI-015-15	Allis Intestinal and Tissue Grasping Forcep 5 x 6 Teeth 15.5 cm - 6" Stainless Steel	6
AF-491-02	Mixer-Baby Dissecting and Ligature Forcep Strongly Curved 13 cm - 5" Stainless Steel	2
AF-471-18	Gemini Dissecting and Ligature Forcep Curved 18 cm - 7" Stainless Steel	2
AF-230-16	Ochsner-Kocher Haemostatic Forcep Straight - 1 x 2 Teeth 16.5 cm - 6 1/2" Stainless Steel	2
TM-250-15	Lahey Goitre Grasping Forcep 2 x 2 Teeth 15.5 cm - 6" Stainless Steel	2
AF-773-13	Backhaus Towel Clamp 13 cm - 5" Stainless Steel	1
SU-038-14	Hegar-Baumgartner Needle Holder 14.5 cm - 5 3/4" Stainless Steel	2
TR-101-08	Frazier Suction Tube 8 Charr. 19.5 cm - 7 3/4" ø 2.6 mm Stainless Steel	12
SL-010-03	Bard-Parker Scalpel Handle No. 3 Solid 12 cm - 4 3/4" Stainless Steel	1
SL-010-04	Bard-Parker Scalpel Handle No. 4 Solid 13.5 cm - 5 3/8 " Stainless Steel	1

ALPHABETICAL INDEX

	Pages
Broad Butterfly Forceps	344-345
Brown Hair Hook	389
Brimson Tickle Wiping Holder	368
Breiley Dressing Forceps	368
Burley Dressing Forceps	368
Burke O-Lipman Novette	421
Burkhardt Retractor	198-199
Burrone No. 5 Shear	141
Bryant Dressing Forceps	144-145
Buchs Retractor	139
Bulford Hill Spreader	143
Bullock Retractor	143
Bullock Suture	143-144
Buncke's Forceps	140-141
C	432
Caine Haemostatic Forceps	389
Carey-Roth Retractor	389
Carrel Ligation Clasper	389
Catgut Retractor	389
Catgut Surgical Suture	389
Gastric Airway Nasotracheal Vacuole Clamp	389
Gastric Retractor	389
Gastronome Airway Spreader	389
Gastronome Auto Needle Holder	389
Gastronome Forceps	389
Gastronome Dressing Forceps	389
Gauze Dressing Forceps	389
Chirurgic Appretation Forceps	389
Chirurgic Ligature Needle	389



SURTEX
INSTRUMENTS

Crafted to Excellence

Alphabetical Index

A	Pages
Abraham Tonsil Knife	166
Adams Napkin Holder	1
Adson Ear Forcep	52
Aiach Columella Suturing Forcep	157
Andrew Tongue Depressor	11
Asch Septum Straightening Forcep	157
Aufrecht Glabella Rasp	149
Aufrecht Nasal Retractor,	155-156
Austin Attic Dissector,	73

B	Pages
Ballenger Septum Gouge	133
Ballenger Swivel Knife	125
Ballenger-Hajek Septum Elevator	128
Barbara Micro Ear Needle	67-68
Barnhill Adenoid Curette	162
Barth Ear Scoop	84
Bauer Intra Oral Retractor	15
Bauer Osteotome	32
Beckmann Adenoid Curette	161
Bellucci Micro Scissor	61
Beyer Mastoid Rongeur	92
Billeau Ear Loop	85
Blakesley Nasal Cutting Forcep	115
Beales Elevator	73
Blohmke Tonsil Grasping Forcep	171
Boise Nasal Fracture Elevator	31
Bruenings Broncho-Esophagoscopy Forcep Tip	179
Bruenings Ethmoid Forcep,	108
Bruenings Interchangeable Tube	177
Bruenings Tongue Depressor	12
Bruenings Tonsil Snare	173
Buck Ear Curette	82-83
Buck Foreign Body Forcep	55
Buckingham Chisel	75

C	Pages
Carpenter Tonsil Dissector	166
Cawood-Minnesota Cheek Retractor	10
Chevalier	184-187
Chevalier-Jackson Biopsy Forcep	185-187
Chevalier-Jackson Grasping Forcep	184-185
Chevalier-Jackson Punch	184, 186
Cinelli Osteotome	139
Circular Cutting Punch	113
Citelli Mastoid Rongeur	91
Coakley Antrum Curette	159-160
Collin Tongue Holding Forcep	13
Colver Tonsil Grasping Forcep	172
Converse Nasal Retractor	154-155
Converse Rhinoplasty Knife	123
Cordes Broncho-Esophagoscopy Forcep Tip	180
Corwin Tonsil Haemostat	175
Cottle Alar Hook	154
Cottle Alar Retractor	156
Cottle Cartilage Crusher	144
Cottle Nasal Chisel	135
Cottle Nasal Osteotome	134-136, 140
Cottle Nasal Rasp	149
Cottle Nasal Speculum	97-99
Cottle Rhinoplasty Knife	124
Cottle Septum Elevator	129-131
Cottle Skiving Grasping Forcep	158
Cottle-Masing Needle Hodler	119
Cottle-Millard Finger Ring Retractor	154

C	Pages
Cottle-Neivert Alar Hook	153
Cottle-Walsham Septum Straightening Forcep	157
Craig Septum Forcep	106
Craig-Dominick Ethmoid Forcep	107

D	Pages
Daniel Periosteal Elevator	127
Dautrey-Munro Osteotome	33-34
Davis-Boyle Mouth Gag	6
Day Ear Hook	87
Deerlacki Ossicle Holding Clamp	158
Denhart Mouth Gag	2
Denis-Browne Tonsil Holding Forcep	172
Dingmann Bone Hook	32
Dingmann Mouth Gag	7
Dingmann Otobrader	76
Dingmann Tongue Depressor	7
Dingmann Zygoma Elevator	31
Dolwick-Type TMJ Condylar Retractor	35
Dolwick-Type TMJ Disc Clamp	35
Dolwick-Type TMJ Elevator	35
Draffin Bipod	9
Duckbill Alligator Forcep	58
Dunn-Dautrey Elevator	34

E	Pages
Ear Syringe	92
Eicken-Killian Irrigation Cannula	165
Epker-Bauer Osteotome	33
Eves Tonsil Snare	173

F	Pages
Farror Ear Speculum	42
Faulkner Antrum Curette	159
Fergusson-Ackland Mouth Gag	3
Fisch Footplate Perafortor	66
Fisch Measuring Rod	75
Fisch Micro Dissector	73
Fisch Tenotome	64
Fischer Tonsil Knife / Dissector	166
Fishbone Laryngeal Grasping Forcep	176
Fomon Alar Hook	153
Fomon Chisel	138
Fomon Nasal Rasp	148, 152
Fomon Rhinoplasty Knife	123
Fraenkel Broncho-Esophagoscopy Forcep Tip	178
Freer Mucosa Knife	124
Freer Septum Chisel	132
Freer Septum Elevator	127-128
Freer Septum Knife	125
Freer-Yasargil Septum Elevator	128
Fuller Introducing Forcep	62

G	Pages
Gallaher Antrum Rasp	148
Gillies Zygoma Elevator	31
Gross Ear Scoop / Hook	84
Gruenwald (Jansen) Nasal Tampon Forcep	102
Gruenwald (Lucae) Ear Forcep	52-53

Alphabetical Index

Pages	
G	
Guilford-Wright Duckbill Elevator	73
Gwynne-Evans Tonsil Elevator,	168
H	
	Pages
Hajek Sphenoid Punch	121
Hajek-Claus Antrum Punch	122
Hajek-Kofler Sphenoid Punch	121
Halle Antrum Curette	137
Halle Septum Elevator	126
Hargis Anterior Border Stripper	22
Hargis Bone Holding Forcep	23
Hargis Medial Ramus Retractor	21-22
Hargis Periform Rim Retractor	23
Hargis Perosteal Elevator	22
Hargis Pterygoid Osteotome	23
Hartmann Alligator Forcep	56-57, 59
Hartmann Ear Catheter	49
Hartmann Ear Hook,	86
Hartmann Ear Polypus Forcep	54
Hartmann Ear Speculum	39-41
Hartmann Nasal Cutting Forcep (Conchotome)	114
Hartmann Nasal Polypus Forcep	104-105
Hartmann Tongue Depressor	11
Hartmann Tonsil Punch	175
Hartmann Tuning Fork	44-45, 47
Hartmann-Halle Nasal Speculum	95
Hartmann-Hoffmar Alligator Punch Forcep	58
Hartmann-Weingartner Alligator Forcep	58
Hartmann-Wullstein Alligator Forcep	55
Haslinger Tonsil Retractor	169
Heister Mouth Gag	1
Henahan Retractor	14
Henke Tonsil Elevator	167
Heuwieser Antrum Grasping Forcep	108
Heymann Septum Forcep	106
Heymann-Knight Nasal Polypus Forcep	104
Hildmann Cartilage Forcep	159
Hofer Tonsil Retractor	169
House Micro Ear Curette	76-77
House Micro Measuring Rod	74
House-Dieter Malleus Nipper	62
Howarth Periosteal Elevator	131
Huber Universal Handle	176
Hughes Elevator	73
Hurd Tonsil Dissector / Retractor	167
I	
	Pages
Iterson Tracheal Hook	192
J	
	Pages
Jansen (Gruenwald) Nasal Tampon Forcep	102
Jansen Ear Forcep, 53Jansen Ear Scoop	84
Jansen-Middleton Septum Forcep	120
Jennings Mouth Gag	4
Joseph Bone File	147
Joseph Nasal Saw	143-144
Joseph Rhinoplasty Knife	123-124
Joseph Septum Elevator	126
Jurasz Adenoid Cutting Forcep	176
K	
	Pages
Kerrison Sphenoid Punch	122
Killian Nasal Polypus Forcep	105
Killian Nasal Speculum	96-97
Killian Septum Elevator	129
Killian-Claus Septum Gouge	133
Kilner-Doughty Mouth Gag	5
Kilner-Doughty Tongue Depressor	5
Knight Septum Forcep	106
Koenig Tongue Depressor	12
Koenig Tuning Fork	48
Krause Broncho-Esophagoscopy Forcep Tip	179
Krause Ear Polypus Snare	49
Krause-Nasal Polypus Snare	101
Krause-Voss Ear Polypus Snare	49
Krause-Voss Nasal Polypus Snare	101
Kressner Septum Forcep	156
Kuhn Frontal Sinus Seeker	88
Kuhn-Bolger Frontal Sinus Curette	160
L	
	Pages
LaForce Adenotome	163-164
Lange Alar Hook	153
Lange Broncho-Esophagoscopy Forcep Tip	179
Lange-Wilde Ear Polypus Snare	50
Lange-Wilde Nasal Polypus Snare	101
Langenbeck Ear Curette	81
Langenbeck Finger Protector	1
Laryngeal Mirror	188
Laryngeal Mirror Handle	189
Lempert Ear Scoop	84-85
Lempert Periosteal Elevator	91
Levasseur-Merrill Intra Oral Retractor	15
Lewis Nasal Rasp	150
Littauer Ear Polypus Forcep	54
Littauer Nasal Polypus Forcep	105
Lopez-Reinke Tonsil Dissector	168
Love Tonsil Retractor	168
Lower Lateral Morselizer	120
Lubet-Barbon Nasal Polypus Forcep	104
Luc Ethmoid Forcep	107
Lucae (Gruenwald) Ear Forcep	52-53
Lucae Ear Hook, 86Lucae Tuning Fork	46-47
Lucae Tympanum Needle	90
Luniatschek Gauze Packer	170
M	
	Pages
Maguaran Bipod Plate	9
Maltz Nasal Rasp	146-147, 151
Martin Nasal Speculum	95-96
Martin-Simplex Lip and Cheek Retractor	10
Masing Septum Elevator	131
Mauder Mouth Gag	1
Maxillary Ostium Sinus Seeker	87
Mayo Tongue Depressor	12
McGee Micro Alligator Forcep	59
McGee Micro Ear Needle	69-70
McGee Wire Bending Forcep	62
McIvor Mouth Gag	7
McKenty Septum Elevator	126, 128
McKesson Mouth Prop	9
Michel Laryngeal Mirror	189
Micro Alligator Forcep	60
Micro Ear Needle	68-69, 71
Micro Ear Pick	71-72
Micro Laryngeal Forcep	182

Alphabetical Index

M Pages

Micro Laryngeal Scissor	181
Micro Laryngeal Scissors	187
Micro Measuring Rod	75
Miller Bone File	146
Miller-Colburn Bone File	145
Mini-McGee Micro Alligator Forcep	59
Minnesota Cheek Retractor	10
Mod. Vienna Nasal Speculum	94
Moldenhauer Foreign Body Lever	90
Mollison Tonsil Dissector / Retractor	167
Molt Mouth Gag	2
Museux Tonsil Grasping Forcep	172

N Pages

Nager Tonsil Retractor	169
Negus Chest Support For Mouth Gag	4
Negus Knot Pusher	170
Negus Tongue Depressor	8
Neivert Nasal Retractor	156
Neivert-Anderson Osteotome	140
Neuro Tuning Fork	48

O Pages

Obwegeser Channel Retractor	21
Obwegeser Chin Retractor	20
Obwegeser J-Stripper	25
Obwegeser Lip and Cheek Retractor	10
Obwegeser Maxillary Awl	28
Obwegeser Mouth Gag	1
Obwegeser Osteotome	39
Obwegeser Periosteal Elevator	24
Obwegeser Ramus Retractor	18
Obwegeser Ramus Retractor (Nasal Spine Retractor)	19
Obwegeser Ramus Stripper	25
Obwegeser Septonasal Osteotome	36
Obwegeser Soft Tissue Retractor	16-17
Obwegeser Split Osteotome	37
Obwegeser Wedge Osteotome	37
Obwegeser Wire Loop Forcep	29
Obwegeser Zygomatic Arch Awl	28
Ogura Skin Grasping Forcep	158

P Pages

Parkes Nasal Rasp	148
Parkes Osteotome	142
Pennington Septum Elevator	130
Petri Intra Oral Retractor	15
Pierce Septum Elevator	129
Plester Sickle Knife	64
Plester Wire Guide and Aplicator	66
Politzer Furuncle Knife	89
Politzer Tympanum Needle	89

Q Pages

Quire Foreign Body Forcep	55
---------------------------	----

R Pages

Richnau-Holmgren Ear Speculum	43
Ritter-Halle Frontal Sinus Dilator	165
Roeder Tonsil Retractor	169
Roger Periosteal Elevator	130
Rosen Circular Cut Knife	65
Rosen Micro Curette	74
Rosen Micro Dissector	74
Rosen Micro Ear Periosteal Elevator	73
Rosenberg Broncho-Esophagoscopy Forcep Tip	180
Roser-Koenig Mouth Gag	3
Rowe Disimpaction Forcep	30
Rowe Zygoma Elevator	31
Rubin Nasal Osteotome	136
Rubin Septum Forcep	120
Russel-Davis Tongue Depressor	8
Rydel-Seiffer Tuning Fork	47

S Pages

Scheinmann Broncho-Esophagoscopy Forcep Tip,	178
Schmeden Nasal Cutting Forcep,	114
Schmeden Tonsil Punch,	175
Schnidt Tonsil Forcep,	171
Schukenchit Rotary Knife,	63
Schumacher Broncho-Esophagoscopy Forcep Tip,	180
Seiffert Broncho-Esophagoscopy Forcep Tip,	178
Sexton Myringotomy Knife,	64
Sexton Tympanum Knife,	65
Shambaugh Micro Ear Needle,	72
Shapleigh Ear Curette,	81
Shea Exploring Hook,	72
Shea Micro Ear Curette,	77
Shea Micro Scissor,	62
Sheehan Nasal Osteotome,	141
Siegle Otoscope,	43
Sigmoid Notch Extra Oral Retractor,	14
Silver Nasal Chisel,	142
Simplex (Martin) Lip and Cheek Retractor,	10
Sluder-Ballenger Tonsillectome,	174
Smith Ramus Separator,	29
Smith Sagittal Split Separator,	30
Snare Wire,	50, 101, 173
St. Clair Thompson Adenoid Curette,	162
Stammberger Nasal Cutting Forcep,	112
Standard Tongue Depressor,	12
Steinhauer Reposition Instrument,	26
Steinhauer Zygomatic Retractor,	20
Storage Rack for Speculum,	43
Struempel Alligator Forcep,	58
Struycken Nasal Cutting Forcep,	111

T Pages

Tabb Mucosa Knife	63
Takahashi Ethmoid Forcep	109
Tatum Sinus Lift Curette	27
Terry Ramus Pusher	26
Thudichum Nasal Speculum	93-94
Tieck-Halle Nasal Specula	96
Tilley Ear Polypus Forcep	54
Tilley-Henckel Nasal Cutting Forcep	111
Tobold Tongue Depressor	11
Traquair Periosteal Elevator	130
Troeltsch (Wilde) Ear Forcep	51-52
Troeltsch Nasal Tampon Forcep	102
Tuning Fork	44-48
Turvey Palatal Expander	36

Alphabetical Index

T	Pages
---	-------

Tyding Tonsil Snare	174
---------------------	-----

W	Pages
---	-------

Walter Nasal Chisel,	137
Walter Nasal Osteotome,	137
Watson-Williams Ethmoid Forcep,	109
Weber-Loch Ear Loop,	85
Webster Glabella Nasal Rasp,	149
Weil Ethmoid Forcep,	109
Weil-Blakesley Nasal Cutting Forcep,	110
West Septum Chisel,	132
Westmacott Nasal Polypus Forcep,	104
White Tonsil Grasping Forcep,	170
Whitehead Mouth Gag,	4
Wieder Tongue Depressor,	12
Wiener Nasal Rasp,	150
Wilde (Troeltsch) Ear Forcep,	51-52
Wilson Tonsil Artery Forcep,	171
Wullstein Micro Scissor,	61
Wullstein Sickle Knife,	63

Y	Pages
---	-------

Young Tongue Holding Forcep	13
-----------------------------	----

Z	Pages
---	-------

Zaufal Foreign Body Lever	90
Zoellner Application Instrument	81
Zoellner Ear Knife	79
Zoellner Ear Speculum	41
Zoellner Micro Ear Needle	78
Zoellner Respiratory	80
Zoellner Tympanum Needle	80
Zurich Nasal Scissor	106

Numerical Index

Cat. No.	Page	Cat. No.	Page	Cat. No.	Page	Cat. No.	Page
OM-011-01,	2	OM-151-14,	13	OM-463-07,	28	OT-022-01,	43
OM-015-01,	2	OM-153-15,	13	OM-465-08,	28	OT-022-02,	43
OM-015-02,	2	OM-155-17,	13	OM-465-09,	28	OT-022-03,	43
OM-015-03,	2	OM-159-19,	13	OM-465-10,	28	OT-022-04,	43
OM-018-10,	2	OM-170-17,	14	OM-501-15,	29	OT-030-01,	43
OM-030-13,	2	OM-173-17,	14	OM-501-17,	29	OT-030-02,	43
OM-032-11,	2	OM-201-01,	15	OM-503-24,	29	OT-030-03,	44
OM-035-11,	3	OM-201-02,	15	OM-507-18,	29	OT-032-00,	43-44
OM-035-14,	3	OM-201-03,	15	OM-540-15,	30	OT-032-01,	44
OM-039-12,	3	OM-251-24,	15	OM-550-20,	30	OT-032-02,	43-44
OM-041-16,	4	OM-255-24,	16	OM-554-01,	31	OT-032-03,	44
OM-041-20,	4	OM-259-21,	16	OM-554-02,	31	OT-032-04,	43-44
OM-043-18,	4	OM-261-21,	16	OM-570-01,	31	OT-040-01,	44
OM-043-20,	4	OM-265-25,	16	OM-570-02,	31	OT-042-01,	44
OM-047-09,	5	OM-271-16,	17	OM-601-21,	32	OT-042-02,	44
OM-047-11,	5	OM-271-25,	17	OM-602-21,	32	OT-042-03,	44
OM-047-13,	5	OM-271-35,	17	OM-605-21,	32	OT-042-04,	44
OM-047-15,	5	OM-271-42,	17	OM-607-19,	32	OT-042-05,	44
OM-049-11,	5	OM-271-55,	17	OM-613-10,	33	OT-042-06,	44
OM-049-13,	5	OM-271-70,	17	OM-615-20,	33	OT-050-01,	45
OM-049-15,	5	OM-271-80,	17	OM-620-06,	33	OT-052-01,	45
OM-051-13,	5	OM-273-16,	17	OM-620-08,	33	OT-052-02,	45
OM-058-00,	5	OM-273-25,	17	OM-621-06,	33	OT-052-03,	45
OM-060-00,	6	OM-273-35,	17	OM-621-08,	33	OT-052-04,	45
OM-061-00,	6	OM-273-42,	17	OM-623-08,	34	OT-052-05,	45
OM-063-01,	6	OM-273-55,	17	OM-626-04,	34	OT-052-06,	45
OM-063-02,	6	OM-273-70,	17	OM-626-06,	34	OT-052-07,	45
OM-063-03,	6	OM-273-80,	17	OM-626-08,	34	OT-054-00,	45
OM-070-00,	7	OM-275-43,	18	OM-626-10,	34	OT-054-01,	45
OM-071-00,	7	OM-275-57,	18	OM-626-12,	34	OT-054-02,	45
OM-071-01,	7	OM-275-83,	18	OM-626-14,	34	OT-054-03,	45
OM-073-01,	7	OM-279-30,	18	OM-626-16,	34	OT-054-04,	45
OM-073-02,	7	OM-279-45,	18	OM-628-04,	35	OT-054-05,	45
OM-073-03,	7	OM-279-60,	18	OM-628-06,	35	OT-054-06,	45
OM-073-04,	7	OM-281-22,	19	OM-628-08,	35	OT-054-07,	45
OM-073-05,	7	OM-283-22,	19	OM-629-06,	35	OT-054-08,	4
OM-080-00,	8	OM-287-22,	20	OM-631-06,	35	50T-060-12,	46
OM-081-00,	8	OM-289-22,	20	OM-633-48,	35	OT-064-02,	46
OM-083-01,	8	OM-291-01,	21	OM-633-56,	35	OT-064-04,	46
OM-083-02,	8	OM-291-02,	21	OM-637-12,	36	OT-064-06,	46
OM-083-03,	8	OM-301-16,	21	OM-639-14,	36	OT-070-00,	46
OM-090-00,	8	OM-305-08,	22	OM-641-17,	36	OT-072-01,	46
OM-091-00,	8	OM-315-01,	22	OM-700-17,	37	OT-072-02,	46
OM-093-01,	8	OM-315-02,	22	OM-730-19,	37	OT-072-03,	46
OM-093-02,	8	OM-319-18,	23	OM-735-04,	37	OT-100-10,	47
OM-093-03,	8	OM-331-19,	23	OM-735-06,	37	OT-100-20,	47
OM-097-01,	9	OM-333-06,	23	OM-735-08,	37	OT-100-30,	47
OM-097-02,	9	OM-333-09,	23	OM-740-08,	38	OT-102-01,	47
OM-097-03,	9	OM-335-06,	23	OM-740-12,	38	OT-102-02,	47
OM-099-01,	9	OM-335-09,	23	OM-740-16,	38	OT-102-03,	47
OM-099-02,	9	OM-337-18,	24	OM-742-08,	38	OT-102-04,	47
OM-099-03,	9	OM-339-22,	24	OM-742-12,	38	OT-102-05,	47
OM-099-04,	9	OM-341-19,	24	OM-742-16,	38	OT-102-06,	47
OM-099-05,	9	OM-361-04,	25	OM-742-20,	38	OT-102-60,	47
OM-099-06,	9	OM-361-06,	25	OM-750-03,	39	OT-104-01,	48
OM-099-07,	9	OM-361-07,	25	OM-750-04,	39	OT-104-02,	48
OM-099-08,	9	OM-361-09,	25	OM-750-06,	39	OT-104-03,	48
OM-100-00,	10	OM-361-11,	25	OT-010-01,	42	OT-104-04,	47-48
OM-100-01,	10	OM-381-21,	25	OT-010-02,	42	OT-104-05,	47-48
OM-100-48,	10	OM-383-21,	25	OT-010-03,	42	OT-110-00,	48
OM-110-01,	10	OM-401-05,	26	OT-010-04,	42	OT-112-02,	48-49
OM-110-02,	10	OM-401-07,	26	OT-012-00,	42	OT-112-03,	49
OM-110-03,	10	OM-401-09,	26	OT-012-01,	42	OT-112-04,	49
OM-121-11,	11	OM-405-01,	26	OT-012-02,	42	OT-112-05,	49
OM-123-12,	11	OM-405-02,	26	OT-012-03,	42	OT-112-06,	48-49
OM-127-13,	11	OM-415-24,	27	OT-012-04,	42	OT-122-01,	49
OM-131-15,	11	OM-421-18,	27	OT-014-00,	43	OT-122-02,	49
OM-135-14,	11	OM-451-01,	28	OT-014-01,	43	OT-122-03,	49
OM-141-15,	11	OM-453-02,	28	OT-014-02,	43	OT-122-04,	49
OM-143-15,	12	OM-455-03,	28	OT-014-03,	43	OT-122-05,	49
OM-145-14,	12	OM-457-04,	28	OT-014-04,	43	OT-122-06,	49
OM-147-01,	13	OM-459-05,	28	OT-020-01,	43	OT-126-00,	50
OM-147-02,	13	OM-461-06,	28	OT-020-02,	43	OT-128-00,	50

Numerical Index

Cat. No.	Page	Cat. No.	Page	Cat. No.	Page	Cat. No.	Page
OT-130-00,	50	OT-345-08,	64	OT-481-06,	73	OT-570-01,	85
OT-134-00,	50	OT-347-08,	64	OT-481-10,	73	OT-570-02,	85
OT-136-00,	51	OT-351-08,	65	OT-483-03,	73	OT-570-10,	85
OT-138-00,	51	OT-355-08,	65	OT-483-06,	73	OT-570-11,	85
OT-140-00,	51	OT-357-08,	65	OT-483-10,	73	OT-570-12,	85
OT-150-01,	52	OT-361-01,	65	OT-485-03,	74	OT-570-13,	85
OT-150-02,	52	OT-361-02,	65	OT-485-06,	74	OT-570-14,	85
OT-150-03,	52	OT-365-01,	65	OT-485-10,	74	OT-571-01,	85
OT-161-18,	52	OT-365-02,	65	OT-487-03,	74	OT-571-02,	85
OT-163-18,	52	OT-365-03,	65	OT-487-06,	74	OT-571-10,	85
OT-165-16,	53	OT-365-04,	65	OT-487-10,	74	OT-571-11,	85
OT-170-30,	53	OT-367-08,	65	OT-491-03,	74	OT-571-12,	85
OT-170-40,	53	OT-401-01,	66	OT-491-06,	74	OT-571-13,	85
OT-201-11,	54	OT-401-02,	66	OT-491-10,	75	OT-571-14,	85
OT-201-13,	54	OT-403-15,	66	OT-493-03,	75	OT-572-01,	86
OT-203-11,	54	OT-405-15,	66	OT-493-06,	75	OT-572-02,	86
OT-205-11,	55	OT-407-15,	66	OT-493-10,	75	OT-572-10,	86
OT-207-18,	55	OT-409-15,	67	OT-497-20,	75	OT-572-11,	86
OT-211-14,	55	OT-413-16,	67	OT-501-30,	76	OT-572-12,	86
OT-213-16,	56	OT-415-15,	67	OT-502-15,	76	OT-572-13,	86
OT-213-21,	56	OT-417-15,	67	OT-503-01,	76	OT-572-14,	86
OT-215-14,	56	OT-419-16,	68	OT-503-02,	76	OT-573-01,	86
OT-221-10,	57	OT-420-03,	68	OT-505-01,	76	OT-573-02,	86
OT-223-12,	57	OT-421-01,	68	OT-505-02,	76	OT-573-10,	86
OT-225-14,	57	OT-421-02,	68	OT-505-03,	76	OT-573-11,	86
OT-229-14,	57	OT-421-03,	68	OT-505-04,	76	OT-573-12,	86
OT-229-16,	57	OT-421-04,	68	OT-507-02,	76	OT-573-13,	86
OT-235-11,	58	OT-421-05,	68	OT-507-03,	76	OT-573-14,	86
OT-237-10,	58	OT-421-06,	68	OT-509-01,	76	OT-576-15,	87
OT-251-08,	58	OT-421-07,	68	OT-509-02,	76	OT-580-18,	87
OT-253-08,	58	OT-422-03,	69	OT-509-03,	76	OT-583-13,	87
OT-255-08,	59	OT-422-04,	69	OT-509-04,	77	OT-600-10,	87
OT-255-14,	59	OT-422-05,	69	OT-509-05,	77	OT-600-20,	87
OT-255-16,	59	OT-422-06,	69	OT-509-06,	77	OT-600-30,	87
OT-257-08,	59	OT-422-07,	69	OT-511-01,	77	OT-600-40,	87
OT-257-14,	59	OT-422-08,	69	OT-511-02,	77	OT-600-50,	87
OT-261-20,	60	OT-440-16,	69	OT-511-03,	77	OT-610-10,	88
OT-261-25,	60	OT-450-16,	70	OT-511-04,	77	OT-610-20,	88
OT-261-30,	60	OT-451-16,	70	OT-512-16,	78	OT-610-30,	88
OT-262-25,	60	OT-453-16,	70	OT-513-16,	78	OT-610-40,	88
OT-263-20,	60	OT-455-16,	70	OT-515-01,	78	OT-610-50,	88
OT-263-30,	60	OT-459-03,	70	OT-515-02,	78	OT-650-01,	88
OT-265-08,	61	OT-459-06,	70	OT-515-04,	78	OT-660-01,	88
OT-267-08,	61	OT-459-10,	70	OT-517-01,	79	OT-660-02,	88
OT-269-08,	61	OT-461-03,	70	OT-517-02,	79	OT-660-03,	88
OT-271-08,	61	OT-461-06,	70	OT-520-01,	79	OT-671-16,	89
OT-273-08,	62	OT-461-10,	70	OT-520-02,	79	OT-673-01,	89
OT-281-08,	62	OT-461-12,	70	OT-520-03,	79	OT-673-02,	89
OT-283-08,	62	OT-461-15,	70	OT-520-04,	79	OT-673-03,	89
OT-285-08,	62	OT-461-20,	70	OT-521-01,	80	OT-675-01,	89
OT-287-08,	62	OT-461-25,	70	OT-521-02,	80	OT-675-02,	89
OT-289-08,	62	OT-463-03,	71	OT-521-03,	80	OT-675-03,	89
OT-291-08,	62	OT-463-06,	71	OT-525-14,	80	OT-675-04,	89
OT-293-08,	62	OT-463-10,	71	OT-535-01,	81	OT-675-05,	89
OT-301-08,	63	OT-463-12,	71	OT-535-02,	81	OT-677-01,	90
OT-303-08,	63	OT-463-15,	71	OT-535-03,	81	OT-677-02,	90
OT-305-08,	63	OT-463-20,	71	OT-535-04,	81	OT-679-16,	90
OT-307-08,	63	OT-463-25,	71	OT-535-05,	81	OT-681-19,	90
OT-309-08,	63	OT-465-05,	71	OT-535-06,	81	OT-683-01,	91
OT-311-08,	63	OT-465-10,	71	OT-535-07,	81	OT-683-02,	91
OT-313-08,	63	OT-467-03,	71	OT-535-08,	81	OT-683-03,	91
OT-315-08,	63	OT-467-06,	71	OT-535-09,	81	OT-683-04,	91
OT-317-08,	63	OT-467-10,	71	OT-535-10,	82	OT-683-05,	91
OT-321-08,	63	OT-469-03,	72	OT-535-11,	82	OT-683-06,	91
OT-323-08,	63	OT-469-06,	72	OT-535-12,	82	OT-701-16,	92
OT-325-08,	63	OT-469-10,	72	OT-535-13,	82	OT-704-16,	92
OT-327-08,	63	OT-473-05,	72	OT-535-14,	83	OT-707-16,	92
OT-331-08,	64	OT-473-10,	72	OT-535-15,	83	OT-709-16,	92
OT-333-08,	64	OT-475-05,	72	OT-535-16,	83	OT-711-16,	93
OT-335-08,	64	OT-475-10,	72	OT-535-17,	83-84	OT-713-16,	93
OT-337-08,	64	OT-478-16,	73	OT-551-13,	84	OT-721-16,	93
OT-341-08,	64	OT-479-16,	73	OT-555-01,	84	OT-723-16,	93
OT-343-08,	64	OT-481-03,	73	OT-555-02,	84	OT-727-20,	94

Numerical Index

Cat. No.	Page						
OT-729-18,	94	RH-131-20,	110	RH-257-42,	120	RH-451-04,	131
OT-731-02,	94	RH-133-18,	110	RH-257-48,	120	RH-451-05,	131
OT-731-03,	94	RH-135-16,	110	RH-261-00,	122	RH-455-16,	131
OT-731-04,	94	RH-140-01,	111	RH-261-01,	122	RH-461-16,	131
OT-731-05,	94	RH-161-01,	111	RH-261-02,	122	RH-465-01,	132
OT-751-02,	94	RH-161-02,	111	RH-261-03,	122	RH-465-02,	132
OT-751-03,	94	RH-161-03,	111	RH-261-04,	122	RH-465-03,	132
OT-755-19,	95	RH-165-17,	111	RH-261-05,	122	RH-471-19,	132
OT-759-14,	95	RH-167-16,	111	RH-261-11,	121	RH-473-19,	132
OT-800-00,	95	RH-169-01,	112	RH-261-13,	121	RH-475-18,	132
OT-802-00,	95	RH-169-02,	112	RH-261-14,	121	RH-476-18,	133
RH-011-10,	98	RH-171-01,	112	RH-263-00,	122	RH-477-20,	133
RH-011-20,	98	RH-171-02,	112	RH-263-01,	122	RH-479-20,	133
RH-011-30,	98	RH-173-01,	112	RH-263-02,	122	RH-482-19,	133
RH-011-40,	98	RH-173-02,	112	RH-263-04,	122	RH-501-23,	134
RH-012-01,	98	RH-175-01,	113	RH-263-05,	122	RH-503-18,	134
RH-012-03,	98	RH-175-02,	113	RH-265-00,	123	RH-505-23,	134
RH-012-04,	98	RH-175-03,	113	RH-265-01,	123	RH-507-21,	134
RH-012-05,	98	RH-175-04,	113	RH-265-02,	123	RH-509-21,	135
RH-012-07,	98	RH-177-01,	113	RH-265-03,	123	RH-511-19,	135
RH-013-01,	99	RH-177-02,	113	RH-265-04,	123	RH-513-15,	135
RH-013-02,	99	RH-177-03,	113	RH-265-05,	123	RH-541-22,	135
RH-013-03,	99	RH-177-04,	113	RH-267-00,	123	RH-543-21,	136
RH-013-04,	99	RH-177-05,	113	RH-267-01,	123	RH-545-21,	136
RH-013-05,	99	RH-177-06,	113	RH-267-02,	123	RH-547-21,	136
RH-013-06,	99	RH-181-12,	114	RH-267-03,	123	RH-551-16,	136
RH-013-07,	99	RH-183-12,	114	RH-267-04,	123	RH-600-16,	137
RH-017-01,	99	RH-191-19,	114	RH-267-05,	123	RH-601-16,	137
RH-017-02,	99	RH-193-19,	114	RH-269-00,	124	RH-604-16,	137
RH-017-03,	99	RH-195-19,	114	RH-269-01,	124	RH-605-16,	137
RH-021-10,	99	RH-201-01,	115	RH-269-02,	124	RH-608-16,	137
RH-023-01,	99	RH-201-02,	115	RH-269-03,	124	RH-610-02,	138
RH-023-02,	99	RH-201-03,	115	RH-269-04,	124	RH-610-04,	138
RH-023-03,	99	RH-201-04,	115	RH-269-05,	124	RH-610-06,	138
RH-031-10,	100	RH-201-05,	115	RH-271-00,	124	RH-610-08,	138
RH-033-01,	100	RH-203-01,	115	RH-271-01,	124	RH-613-16,	138
RH-033-02,	100	RH-203-02,	115	RH-271-02,	124	RH-615-16,	138
RH-033-03,	100	RH-203-03,	115	RH-271-03,	124	RH-620-02,	139
RH-035-13,	101	RH-203-04,	115	RH-271-04,	124	RH-620-04,	139
RH-037-13,	101	RH-205-01,	115	RH-271-05,	124	RH-620-07,	139
RH-041-13,	101	RH-205-02,	115	RH-271-15,	125	RH-620-09,	139
RH-051-10,	102	RH-205-03,	115	RH-281-19,	125	RH-620-12,	139
RH-053-01,	102	RH-205-04,	115	RH-283-19,	125	RH-620-14,	139
RH-053-02,	102	RH-207-01,	116	RH-291-21,	125	RH-621-02,	139
RH-053-03,	102	RH-207-02,	116	RH-293-21,	125	RH-621-04,	139
RH-053-04,	102	RH-207-03,	116	RH-301-01,	126	RH-621-06,	139
RH-061-10,	103	RH-231-01,	116	RH-301-02,	126	RH-621-09,	139
RH-063-01,	103	RH-231-02,	116	RH-303-01,	126	RH-621-12,	139
RH-063-02,	103	RH-233-01,	116	RH-303-02,	126	RH-626-02,	140
RH-063-03,	103	RH-233-02,	116	RH-305-03,	127	RH-626-04,	140
RH-063-04,	103	RH-233-03,	116	RH-305-04,	127	RH-626-07,	140
RH-065-10,	104	RH-233-04,	116	RH-305-05,	127	RH-626-09,	140
RH-067-01,	104	RH-234-00,	117	RH-305-06,	127	RH-626-12,	140
RH-067-02,	104	RH-235-01,	117	RH-307-20,	127	RH-627-06,	140
RH-067-03,	104	RH-235-02,	117	RH-309-01,	127	RH-630-07,	141
RH-067-04,	104	RH-235-03,	117	RH-309-02,	127	RH-630-14,	140
RH-073-45,	105	RH-237-35,	118	RH-350-15,	128	RH-630-16,	140
RH-073-75,	105	RH-237-45,	118	RH-351-15,	128	RH-630-18,	140
RH-091-26,	106	RH-239-35,	118	RH-355-15,	128	RH-632-04,	141
RH-093-26,	106	RH-239-45,	118	RH-361-16,	128	RH-632-06,	141
RH-095-24,	106	RH-241-01,	119	RH-365-15,	129	RH-632-14,	141
RH-097-12,	106	RH-241-02,	119	RH-367-15,	129	RH-632-16,	141
RH-101-15,	107	RH-241-03,	119	RH-369-15,	129	RH-636-10,	141
RH-103-18,	107	RH-241-04,	119	RH-371-01,	129	RH-636-12,	141
RH-105-21,	107	RH-245-12,	119	RH-371-02,	129	RH-638-11,	141
RH-107-16,	108	RH-251-25,	120	RH-373-01,	130	RH-638-13,	141
RH-109-16,	108	RH-251-35,	120	RH-373-02,	130	RH-638-14,	141
RH-111-18,	108	RH-251-42,	120	RH-380-01,	130	RH-640-03,	142
RH-113-20,	108	RH-251-48,	120	RH-380-02,	130	RH-640-04,	142
RH-121-18,	109	RH-253-35,	120	RH-380-03,	130	RH-642-03,	142
RH-123-20,	109	RH-255-35,	120	RH-381-03,	130	RH-642-04,	142
RH-127-16,	109	RH-257-25,	120	RH-381-04,	130	RH-642-07,	142
RH-129-18,	109	RH-257-35,	120	RH-381-05,	130	RH-642-09,	142

Numerical Index

Cat. No.	Page						
RH-642-12,	142	RH-729-01,	154	RH-853-06,	166	TL-201-15,	181
RH-643-01,	142	RH-729-02,	154	RH-855-01,	167	TL-201-17,	181
RH-643-02,	142	RH-730-21,	154	RH-855-02,	167	TL-213-01,	181
RH-643-03,	142	RH-732-01,	155	RH-855-03,	167	TL-213-02,	181
RH-645-01,	143	RH-732-02,	155	RH-855-04,	167	TL-213-03,	181
RH-645-02,	143	RH-734-01,	155	RH-855-05,	167	TL-215-01,	182
RH-645-03,	143	RH-734-02,	155	RH-857-01,	167	TL-215-02,	182
RH-647-01,	143	RH-735-18,	155	RH-857-02,	167	TL-215-03,	182
RH-647-02,	143	RH-738-19,	156	RH-857-03,	167	TL-251-01,	182
RH-647-03,	143	RH-739-19,	156	RH-857-04,	167	TL-301-01,	182
RH-648-04,	144	RH-742-01,	156	RH-861-01,	168	TL-302-20,	183
RH-648-05,	144	RH-742-02,	156	RH-861-02,	168	TL-302-30,	183
RH-648-06,	144	RH-742-03,	156	RH-861-03,	168	TL-302-40,	183
RH-648-07,	144	RH-742-04,	156	RH-861-11,	168	TL-303-20,	183
RH-649-01,	143	RH-742-05,	156	RH-861-12,	168	TL-310-01,	184
RH-649-02,	143	RH-742-06,	156	RH-861-13,	168	TL-314-01,	184
RH-650-10,	144	RH-742-07,	156	RH-863-01,	169	TL-318-01,	184
RH-650-12,	144	RH-742-08,	156	RH-863-02,	169	TL-321-01,	184
RH-650-14,	144	RH-742-09,	156	RH-863-03,	169	TL-323-01,	184
RH-650-16,	144	RH-742-10,	156	RH-863-11,	169	TL-326-02,	185
RH-654-07,	145	RH-744-21,	157	RH-863-12,	169	TL-326-03,	185
RH-655-07,	145	RH-746-01,	157	RH-863-13,	169	TL-328-02,	185
RH-657-07,	145	RH-746-02,	157	RH-881-20,	170	TL-328-03,	185
RH-660-02,	145	RH-746-03,	157	RH-881-25,	170	TL-334-01,	185
RH-660-04,	145	RH-746-04,	157	RH-881-30,	170	TL-338-01,	185
RH-660-06,	145	RH-746-05,	157	RH-881-40,	170	TL-338-02,	185
RH-660-07,	145	RH-761-22,	158	RH-891-25,	170	TL-342-01,	185
RH-660-09,	145	RH-763-17,	158	RH-891-30,	170	TL-347-01,	186
RH-660-12,	145	RH-765-16,	158	RH-891-40,	170	TL-347-02,	186
RH-664-01,	146	RH-767-20,	158	TL-011-21,	172	TL-347-03,	186
RH-664-02,	146	RH-771-01,	159	TL-015-19,	172	TL-347-04,	186
RH-664-03,	146	RH-771-02,	159	TL-021-25,	172	TL-349-01,	186
RH-664-04,	146	RH-771-03,	159	TL-027-24,	172	TL-353-01,	186
RH-664-05,	146	RH-775-00,	159	TL-031-22,	173	TL-353-02,	186
RH-664-06,	146	RH-781-01,	159	TL-033-21,	173	TL-357-01,	186
RH-664-07,	146	RH-781-02,	159	TL-037-23,	173	TL-359-01,	186
RH-664-08,	146	RH-781-03,	159	TL-039-23,	173	TL-361-01,	187
RH-664-10,	146	RH-781-04,	159	TL-043-22,	174	TL-361-02,	187
RH-664-12,	146	RH-785-19,	160	TL-045-23,	174	TL-361-03,	187
RH-666-05,	147	RH-789-16,	161	TL-051-18,	174	TL-361-04,	187
RH-667-05,	147	RH-791-15,	161	TL-051-20,	174	TL-361-05,	187
RH-669-05,	147	RH-793-14,	161	TL-051-22,	174	TL-361-06,	187
RH-672-01,	147	RH-801-20,	161	TL-055-23,	175	TL-363-01,	187
RH-691-01,	148	RH-804-01,	162	TL-059-24,	175	TL-363-02,	187
RH-691-02,	148	RH-804-02,	162	TL-063-17,	175	TL-363-03,	187
RH-692-19,	149	RH-804-03,	162	TL-065-00,	175	TL-363-04,	187
RH-693-01,	149	RH-807-23,	162	TL-073-19,	176	TL-363-05,	187
RH-693-02,	149	RH-811-15,	162	TL-073-20,	176	TL-363-06,	187
RH-698-00,	149	RH-814-11,	163	TL-080-20,	176	TL-365-01,	187
RH-701-01,	150	RH-816-12,	163	TL-101-18,	176	TL-365-02,	187
RH-701-02,	150	RH-817-14,	163	TL-101-23,	176	TL-365-03,	187
RH-701-03,	150	RH-820-11,	163	TL-105-20,	177	TL-365-04,	187
RH-701-04,	150	RH-822-11,	164	TL-111-19,	177	TL-365-05,	187
RH-701-05,	150	RH-839-21,	164	TL-113-19,	177	TL-365-06,	187
RH-703-01,	150	RH-841-01,	164	TL-115-20,	177	TL-367-01,	188
RH-703-02,	150	RH-841-02,	164	TL-116-19,	178	TL-367-02,	188
RH-703-03,	150	RH-841-03,	164	TL-117-19,	178	TL-367-03,	188
RH-703-04,	150	RH-841-04,	164	TL-119-22,	178	TL-367-04,	188
RH-705-01,	151	RH-841-05,	164	TL-120-20,	178	TL-367-05,	188
RH-707-01,	151	RH-841-06,	164	TL-150-28,	179	TL-367-06,	188
RH-710-18,	151	RH-849-55,	165	TL-152-28,	179	TL-367-07,	188
RH-712-01,	152	RH-849-90,	165	TL-156-29,	179	TL-367-08,	188
RH-712-02,	152	RH-851-01,	166	TL-161-40,	179	TL-367-09,	188
RH-714-01,	152	RH-851-02,	166	TL-161-50,	179	TL-367-10,	188
RH-714-02,	152	RH-851-03,	166	TL-163-00,	179	TL-369-01,	188
RH-720-01,	153	RH-851-04,	166	TL-173-00,	180	TL-369-02,	188
RH-720-02,	153	RH-851-05,	166	TL-181-00,	180	TL-369-03,	188
RH-722-01,	153	RH-851-06,	166	TL-181-01,	180	TL-369-04,	188
RH-722-02,	153	RH-853-01,	166	TL-183-01,	180	TL-369-05,	189
RH-724-01,	153	RH-853-02,	166	TL-183-02,	180	TL-369-06,	189
RH-724-02,	153	RH-853-03,	166	TL-183-03,	180	TL-369-07,	189
RH-724-03,	153	RH-853-04,	166	TL-191-14,	181	TL-369-08,	189
RH-726-21,	154	RH-853-05,	166	TL-201-13,	181	TL-369-09,	189

Numerical Index

Cat. No.	Page	Cat. No.	Page	Cat. No.	Page	Cat. No.	Page
TL-369-10,	189	TL-406-26,	194				
TL-369-11,	189	TL-406-28,	194				
TL-369-12,	189	TL-406-30,	194				
TL-369-13,	189	TL-413-12,	195				
TL-371-01,	188	TL-413-14,	195				
TL-371-02,	188	TL-413-16,	195				
TL-371-03,	188	TL-413-18,	195				
TL-371-04,	188	TL-420-10,	195				
TL-371-05,	189	TM-101-75,	198				
TL-371-06,	189	TM-101-90,	198				
TL-371-07,	189	TM-151-14,	198				
TL-371-08,	189	TM-155-14,	199				
TL-371-09,	189	TM-161-07,	199				
TL-371-11,	189	TM-165-01,	199				
TL-371-12,	189	TM-165-02,	199				
TL-371-13,	189	TM-171-16,	200				
TL-373-01,	188	TM-173-16,	200				
TL-373-03,	188	TM-181-21,	200				
TL-373-04,	188	TM-183-13,	200				
TL-373-05,	189	TM-185-11,	200				
TL-373-06,	189	TM-189-17,	201				
TL-373-07,	189	TM-191-17,	201				
TL-373-08,	189	TM-201-15,	201				
TL-373-13,	189	TM-203-15,	201				
TL-381-30,	190	TM-250-15,	202				
TL-381-40,	190	TM-252-15,	202				
TL-381-50,	190	TM-252-20,	202				
TL-381-60,	190	TM-301-04,	202				
TL-383-30,	190	TM-301-05,	202				
TL-383-40,	190	TM-301-06,	202				
TL-383-50,	190	TM-301-07,	202				
TL-383-60,	190	TM-301-08,	202				
TL-385-30,	191	TM-301-09,	202				
TL-385-40,	191	TM-301-10,	202				
TL-385-50,	191	TM-301-11,	202				
TL-385-60,	191	TM-301-12,	202				
TL-387-30,	191	TM-301-13,	202				
TL-387-40,	191	TM-321-04,	203				
TL-387-50,	191	TM-321-05,	203				
TL-387-60,	191	TM-321-06,	203				
TL-389-30,	192	TM-321-07,	203				
TL-389-40,	192	TM-321-08,	203				
TL-389-50,	192	TM-321-09,	203				
TL-389-60,	192	TM-321-10,	203				
TL-391-30,	192	TM-321-11,	203				
TL-391-40,	192	TM-321-12,	203				
TL-391-50,	192	TM-321-13,	203				
TL-391-60,	192	TM-351-08,	203				
TL-393-30,	193	TM-351-09,	203				
TL-393-40,	193	TM-351-10,	203				
TL-393-50,	193	TM-351-11,	203				
TL-393-60,	193	TM-351-12,	203				
TL-402-08,	194	TM-351-13,	203				
TL-402-10,	194						
TL-402-12,	194						
TL-402-14,	194						
TL-402-16,	194						
TL-402-18,	194						
TL-402-20,	194						
TL-402-22,	194						
TL-402-24,	194						
TL-402-26,	194						
TL-402-28,	194						
TL-402-30,	194						
TL-406-08,	194						
TL-406-10,	194						
TL-406-12,	194						
TL-406-14,	194						
TL-406-16,	194						
TL-406-18,	194						
TL-406-20,	194						
TL-406-22,	194						
TL-406-24,	194						



SCAN FOR MORE CATALOGUES

COPYRIGHT NOTICE

The copyright, trademarks, domain names and other intellectual property rights in all materials and information in this presentation belongs to Surtex Instruments Ltd. All rights are reserved. Technical alterations are possible. This presentation may be used for no other purpose than offering, buying and selling of Surtex Instruments products. No part may be reproduced or copied in any form. In the event of misuse we retain the right to recall our catalogues and take appropriate legal action.

DISCLAIMER

This document provides general information about Surtex Instruments Ltd, its products and services, and summarises general capabilities and offerings, which we deem to be relevant to our customers. Whilst Surtex Instruments has taken appropriate steps to ensure the accuracy of the information contained in the document Surtex Instruments gives no warranty regarding the accuracy or completeness of such information.

The presentation was prepared to provide information and does not constitute a contract. Surtex reserves the right, without prior notice, to change, delete, supplement, or otherwise amend at any time the information, images and offerings contained in the presentation or other documents. Surtex shall not be liable for any loss or damage of whatever nature (whether direct, indirect, consequential or other) may arise as a result of any third party relying on the information contained in this document.

For further detailed information, prices, terms and conditions, customer should contact Surtex Instruments at sales@surtex-instruments.co.uk.

© 2024 SURTEX INSTRUMENTS. ALL RIGHTS RESERVED. PRINTED IN UNITED KINGDOM. EC-01 Rev 1.1



Unit 202 Kingspark Business Centre
152-178 Kingston Road, New Malden
Surrey KT3 3ST United Kingdom
Phone: +44 (0) 208 123 3513

E-mail: sales@surtex-instruments.co.uk
Website: www.surtex-instruments.com



A member of the ABHI
British HealthTech
Innovation