

SURTEX
INSTRUMENTS

Crafted to Excellence





SURTEX
INSTRUMENTS

Crafted to Excellence

OUR PRODUCT BRANDS



A new dimension in cutting and dissecting performance. Our class leading integral self-sharpening technology standard on all Surtex scissors enabling precise cutting with minimum force even after pro-longed use.



These Surtex premium scissors feature finely ground micro serrations along one blade of the scissor preventing tissue slippage. The ultra sharp edge of the other blade cuts through the tissue precisely and atraumatically.



Instruments with this symbol benefit from tungsten carbide (TC) inserts expertly bonded on the cutting edge of the instrument. This material is harder than stainless steel and enables the instruments to last up to 6 times longer than standard instruments.



A true craftsmanship masterpiece and surgeons favourite. These are the best scissors in our range and combine the features of SuperEdge™ & UltraCut™ into a perfect pair of scissors that is precise, durable and reliable.



The only choice when the perfect grip matters. Utilising Creep Feed technology individual teeth are grinded into hard wear resistant tungsten carbide (TC). The result is super sharp teeth with well-defined point and inserts with impeccable gripping properties and instruments that last up to 6 times longer.



These instruments have jaws, which are coated with diamond dust resulting in improved gripping enabling superior manipulation of sutures, needles and tissues.



Our specially designed high strength fibre grip handles resulting in lightweight instruments. These ergonomic handles provide the surgeon with excellent feel, better grip, reduced shock, and less hand fatigue.



Surtex ENT diagnostic and anaesthesia instruments featuring direct bulb illumination of xenon grade bulbs.



This range features distal fibre optic to carry light from xenon bulb allowing cold illumination on operating site in our ENT diagnostics and anaesthesia instruments.



Premium grade ENT diagnostics and anaesthesia instruments combining features of distal fibre optic technology and ultra bright LED illumination ensuring the best in clear vision of operating site.



SURTEX INSTRUMENTS

Crafted to Excellence

● INTRODUCTION TO THE CATALOGUE

Thank you for your interest in Surtex ophthalmic instruments. We remain committed to providing the highest quality products and customer service.

In this special sets catalogue you will find a comprehensive range of ophthalmic surgery sets, which have been developed following detailed consultation with leading ophthalmic surgeons. In addition to a large range of general surgical instruments for use in ophthalmic surgery, we also manufacture specific instruments for ophthalmology. These can be located in Surtex Ophthalmic Catalogue & Surtex General Instruments Catalogue.

We appreciate your patronage and look forward to a successful relationship with you.

● ABOUT SURTEX INSTRUMENTS

Surtex Instruments is a privately owned medical technology company dedicated to innovation and excellence in crafting of surgical instruments. The company started as a specialist in reconstructive surgery instruments and based on its success and innovative drive gradually expanded to cover a vast range of surgical instruments.

We now manufacture over 8,000 surgical instruments and offer innovative solutions in trauma, gynecology, cardiovascular, microsurgery and dental surgery. All our products are available to view online where you can download our entire range of catalogues or request a printed version.

● SURTEX QUALITY STANDARDS

In European Union medical devices are subject to a series of uniform laws and the EC directive places rigid demands on surgical instruments. Surtex Instruments is proud to claim that we have always surpassed the legal requirements to ensure that each instrument reaching the surgeons hand is indeed a testament to surgical craftsmanship.

We are **BS EN ISO 13485:2012/ISO 9001:2008 (SGS UK)** certified and **CE** approved to ensure compliance with the standards defined in line with medical devices directive **93/42/EEC**.

Surtex Instruments are manufactured to conform to the following standards (where applicable)

- BS 5194-2**
- BS-5194-3**
- BS-5194-4**
- BS EN ISO 13402**

Our instruments are manufactured from steel conforming to **BSEN ISO 7153-1** or **BS 5194-1** (where applicable) and steel of European origin is used to form the blanks required for surgical instrument manufacturing.

WHY IS SURTEX FOCUSED ON INSTRUMENTS QUALITY?

Surgeons are highly talented and skilled individuals. They operate under significant time constraints and are required to make split second decisions while performing complex operations. At Surtex we understand that in such a scenario the surgeon requires tools that are not just fit for the purpose but are reliable, dependable and crafted by a company that has a track record for relentless pursuit for excellence in the making of surgical instruments.



● AFFORDABLE QUALITY

Surgical instruments represent a significant proportion of hospital's total assets. Hospitals and health centers are often required to minimize costs of surgical instruments without impacting quality or patient's health. Decision makers tend to respond via group purchasing contracts & multiple quotes to shop around for the lowest price. However this approach has significant drawbacks as numerous instrument companies have focused towards competing on price alone with a disregard towards instruments quality while claiming their instruments to be guaranteed and free of defect.

The outcome for the buyer being a perceived notion of a quality instrument, which is actually little more than a knock-off and delivers only poor value for money.

At **Surtex Instruments** our success has been based on helping customers navigate successfully through the instrument procurement process by identifying real value for money through fair comparison. We believe that quality should be affordable and you should not pay too much for over priced brands and neither should you pay too much for a low quality instrument.

● LIFE TIME WARRANTY

We encourage our customers to calculate lifetime costs prior to making purchase decisions and investing in surgical instruments.

Surtex Instruments are guaranteed for life* against manufacturing defects of material and workmanship. This means we will repair or replace free of charge any surgical instrument that does not meet its functional requirements when used for its intended purpose and maintained properly.




*Conditions apply.




SURTEX
INSTRUMENTS

Crafted to Excellence



OP



OPHTHALMIC SURGERY INSTRUMENTS

Ophthalmology involves the care of eye conditions in patients of all ages from premature neonates right through to the elderly. The range of conditions encountered in ophthalmology is wide - trauma, cataracts, diabetic eye disease, congenital and genetic eye problems.

Surtex realises that ophthalmology is an extremely specialised field, which places the highest demands on its surgeons therefore Surtex Ophthalmic instruments are meticulously designed, lightweight, non-corrosion, rust-proof titanium and stainless steel to ensure perfect ergonomics for the surgeons hand.

Our offering includes a broad range of instruments for anterior and posterior eye segment surgeries. At your disposal are meticulously crafted eye surgery forceps, scissors, cannulas, kives, diamond knives, phaco choppers, needle holders available as stand alone instruments or in the form of procedure sets for Glaucoma, Enucleation and Cataract etc.

Augenheilkunde Instrumente

Oftalmológica cirugía Instrumentos

Ophthalmique chirurgie Instruments

Oftalmica chirurgia Strumenti



Eye Speculum
Lidsperrer
Blefarostatos
Blépharostats
Blefarostati


Barraquer Wire Speculum
Standard Wire

Catalog #	Blade Size
OP-010-10	10 mm
OP-010-15	15 mm



Barraquer Wire Speculum
Heavy Wire

Catalog #	Blade Size
OP-012-10	15 mm



Barraquer Wire Speculum
Solid Blades

Catalog #	Blade Size
OP-020-10	10 mm
OP-020-15	15 mm




Feaster Wire Speculum
Solid Blades

Catalog #	Blade Size
OP-022-20	20 mm




Kartz-Barraquer Wire Speculum
With Open Blades, Standard Wire

Catalog #	Blade Size
OP-030-13	13 mm
OP-030-15	15 mm



OP-030-13




OP-030-15



Eye Speculum
Lidsperrer
Blefarostatos
Blépharostats
Blefarostati

Kartz-Barraquer Wire Speculum
With Open Blades, Heavy Wire

Catalog #	Blade Size
OP-032-15	15 mm




Frankel Infant Lid Speculum
Adjustable

Catalog #	Blade Size
OP-040-10	10 mm




Alphonson Infant Lid Speculum
With Finger Grip

Catalog #	Size
OP-042-00	N/A



Lid Speculum (For Infants)

Catalog #	Blade Size
OP-046-00	4 mm



Barraquer Infant Lid Speculum

Catalog #	Blade Size
OP-050-09	9 mm





Eye Speculum
Lidsperrer
Blefarostatos
Blépharostats
Blefarostati

Cook Eye Speculum

Catalog #	Blade Size
OP-060-08	8 mm
OP-060-11	11 mm
OP-060-15	15 mm



Murdoch Eye Speculum

Catalog #	Blade Size
OP-064-12	12 mm
OP-064-15	15 mm




Melinger Eye Speculum

Catalog #	Blade Size
OP-070-15	15 mm




Williams Eye Speculum
With Locking Screw

Catalog #	Blade Size
OP-074-11	11 mm
OP-074-15	15 mm



Lancaster Eye Speculum
Adjustable

Catalog #	Blade Size
OP-076-16	16 mm





Eye Speculum
Lidsperrer
Blefarostatos
Blépharostats
Blefarostati


Lester-Burch Eye Speculum
Adjustable

Catalog #	Size
OP-078-00	Standard



Tanna-Lieberman Eye Speculum
Adjustable

Catalog #	Blade Size
OP-081-15	15 mm



Lieberman Eye Speculum
Adjustable

Catalog #	Blade Size
OP-083-10	10 mm
OP-083-15	15 mm



Lieberman Eye Speculum
Adjustable

Catalog #	Blade Size
OP-085-15	15 mm



Lieberman Eye Speculum
Adjustable, Solid Blades

Catalog #	Blade Size
OP-087-10	10 mm






Eye Speculum, Scleral Fixation Ring
Lidsperrer, Scleral Fixation Ring
Blefarostatos, Anillo de fijación escleral
Blépharostats, Anneau de fixation scléral
Blefarostati, Sclerale Fixation Anello

Lieberman Eye Speculum
Adjustable




Catalog #	Blade Size
OP-087-15	15 mm

Castroviejo Eye Speculum
Adjustable Mechanism




Catalog #	Blade Size
OP-091-14	14 mm
OP-091-16	16 mm

Flieringa Scleral Fixation Ring
Set of 9 Ref: OP-100-14 to OP-100-22



Catalog #	Size
OP-100-00	Standard






Scleral Fixation Ring, Depressors
Scleral Fixation Ring, Depressors
Anillo de fijación escleral, Depresores
Anneau de fixation scléral, Dépresseurs
Sclerale Fixation Anello, Depressore


Flieringa Scleral Fixation Ring

Catalog #	Blade Size
OP-100-14	14 mm
OP-100-15	15 mm
OP-100-16	16 mm
OP-100-17	17 mm
OP-100-18	18 mm
OP-100-19	19 mm
OP-100-20	20 mm
OP-100-21	21 mm
OP-100-22	22 mm




McNeill-Goldman Scleral Fixation Ring & Blepharostat
Small, I.D

Catalog #	Blade Size
OP-111-15	15 mm




McNeill-Goldman Scleral Fixation Ring & Blepharostat
Medium, I.D

Catalog #	Blade Size
OP-111-16	16 mm



McNeill-Goldman Scleral Fixation Ring & Blepharostat
Large, I.D

Catalog #	Blade Size
OP-111-17	17 mm





Depressors, Eyelid Retractors
 Depressors, Augenlid Wundhaken
 Depresores, Retractor de párpados
 Dépresseurs, Écarteur de paupières
 Depressore, Palpebra Divaricatore

Gass Scleral Marker

Catalog #	Tip Diameter
OP-151-01	Ø 1.5 mm



Schepens Scleral Depressor
Small Thimble

Catalog #	Tip Width
OP-153-01	7 mm



Schepens Scleral Depressor
Medium Thimble

Catalog #	Tip Width
OP-153-02	7 mm




Schepens Scleral Depressor
Large Thimble

Catalog #	Tip Width
OP-153-03	7 mm



Schoket Scleral Depressor
Double End With Clip

Catalog #	Size
OP-160-14	14 cm (5 1/2")





Eyelid Retractors
Augenlid Wundhaken
Retractor de párpados
Écarteur de paupières
Palpebra Divaricatore

Schepens Orbital Retractor



Catalog #	Size
OP-171-15	15 cm (6")

Jaffe Wire Lid Retractor
Pair



Catalog #	Blade Size
OP-173-13	13 mm
OP-173-15	15 mm

OP-173-13, 13 mm

OP-173-15, 15 mm

Desmarres Lid Retractor
Thin Solid Blades



OP-175-11

OP-175-15

OP-175-13

OP-175-17

Catalog #	Size, Blade Size, Profile
OP-175-11	13 cm (5"), 11 mm, Size 0
OP-175-13	13 cm (5"), 13 mm, Size 1
OP-175-15	13 cm (5"), 15 mm, Size 2
OP-175-17	13 cm (5"), 17 mm, Size 3



Scleral Shield

Set of 2 Ref: OP-180-01 to OP-180-02



Catalog #	Size
OP-180-00	Standard

Scleral Shield
Pair



Catalog #	Blade Size, Profile
OP-180-01	20 x 28 mm, Small
OP-180-02	22 x 30 mm, Large



Eyelid Retractors
Augenlid Wundhaken
Retractor de párpados
Écarteur de paupières
Palpebra Divaricatore

Jaeger Lid Plate
Stainless Steel



Catalog #	Size
OP-192-11	11 cm (4 1/4")

Wright Fascia Needle
Oval Hole



Catalog #	Size, Width, Length
OP-197-14	14 cm (5 1/2"), 1 x 6 mm

Walton Conjunctiva Retractor
Thin Solid Blade



Catalog #	Size, Blade Width
OP-201-12	12.5 (5"), 6 mm

Helveston Tissue Retractor
Thin Curved Blade



Catalog #	Size, Blade Width
OP-211-07	12.5 (5"), 7 mm
OP-211-09	12.5 (5"), 9 mm
OP-211-11	6.5 (2 1/2"), 11 mm

Stevenson Lacrimal Sac Retractor
3 x 3 Blunt Prongs, Adjustable



Catalog #	Size
OP-220-00	Standard

INCHES
MILLIMETERS



Keratome
Keratome
Queratomos
Kératotomes
Cheratomo

Stevenson Lacrimal Sac Retractor


Solid Blades With Serrated Edge, Adjustable



Catalog #	Blade Size
OP-222-00	13 mm

Agricola Lacrimal Sac Retractor


3 x 3 Blunt Prongs



Catalog #	Size
OP-226-00	Standard

Jaeger Keratome

Straight




Catalog #	Size, Profile
OP-240-01	13 cm (5"), Fig. 1
OP-240-02	13 cm (5"), Fig. 2
OP-240-03	13 cm (5"), Fig. 3



Jaeger Keratome

Angled



Catalog #	Size, Profile
OP-241-01	13 cm (5"), Fig. 1
OP-241-02	13 cm (5"), Fig. 2
OP-241-03	13 cm (5"), Fig. 3



Took Corneal Knife

Straight, Blade With Curved Cutting Edge



Catalog #	Size
OP-250-11	11 cm (4 1/4")



Corneal Knives
Hornhaut-Messer
Cuchillos corneales
Couteaux cornée
Coltelli corneali

Corneal Dissector
Angled, Straight Blade



Catalog #	Size
OP-253-12	12 cm (4 3/4")

Corneal Dissector
Angled, Curved Blade



Catalog #	Size
OP-255-12	12 cm (4 3/4")

Gill Corneal Knife
Curved Cutting Edge, Straight



Catalog #	Size
OP-260-12	12.5 cm (5")

Gill Corneal Knife
Curved Cutting Edge, Light Curved



Catalog #	Size
OP-261-12	12.5 cm (5")

Fukasaku Hockey Knife
Sharp Front Edge, Semi Sharp Bottom Edge



Catalog #	Size, Front Edge, Bottom Edge
OP-265-12	12 cm (4 3/4"), 3 mm, 6 mm

INCHES
MILLIMETERS



Gouges and Spuds
Hohlmeißel und meißel
Gubias para Cinceles
Gouges et Ciseaux-burins
Sgorbie e Scalpelli

Golf Club Foreign Body Spud
Angled, Lightly Rounded Tip



Catalog #	Size
OP-271-12	12 cm (4 3/4")

Francis Foreign Body Spud
Angled, Flag Shaped Tip




Catalog #	Size
OP-273-12	12 cm (4 3/4")

Ellis Foreign Body Spud
Lightly Curved and Rounded Tip



Catalog #	Size
OP-275-11	11 cm (4 1/4")

Spud & Guage
In Reversible Handle



Catalog #	Size
OP-280-18	18 cm (7")

Wills Eye Hospital Foreign Body Spud & Burr
Delicately Curved



Catalog #	Size
OP-283-11	11 cm (4 1/4")



Diamond Knives
Diamantmesser
Bisturì de diamante
Couteaux de diamant
Coltelli di diamanti

LRI Diamond Step Knife
Pre-Calibrated



Catalog #	Size
OP-300-12	12 cm (5")

LRI Diamond Micrometer Knife



Catalog #	Size
OP-302-12	12 cm (4 3/4")

AK Diamond Knife



Catalog #	Size
OP-305-11	11 cm (4 1/4")

Diamond Micrometer Knife
45 ° Cutting Edge



Catalog #	Size, Width
OP-311-01	13 cm (5"), 1.0 mm

Diamond Micrometer Knife
Double Vertical and 45 ° Cutting Edge



Catalog #	Size, Width
OP-311-02	13 cm (5"), 1.0 mm



Diamond Micrometer Knife
Triple Vertical and 15 ° Cutting Edge




Catalog #	Size, Width
OP-311-03	13 cm (5"), 1.0 mm

Diamond Knife
30° Cutting Edge




Catalog #	Size, Width
OP-316-01	12.5 cm (5"), 1.0 mm

Diamond Knife
45° Cutting Edge



Catalog #	Size, Width
OP-316-02	12.5 cm (5"), 1.0 mm

Diamond Knife
Tri-Facet Cutting Edge



Catalog #	Size, Width
OP-316-03	12.5 cm (5"), 1.0 mm

Diamond Knife
Spear Cutting Edge




Catalog #	Size, Width
OP-316-04	12.5 cm (5"), 1.0 mm




Diamond Knives
Diamantmesser
Bisturi de diamante
Couteaux de diamant
Coltelli di diamanti

Diamond Knife
10-Facet Cutting Edge



Catalog #	Size, Width
OP-316-05	12.5 cm (5"), 1.0 mm

Diamond Knife
Single 45° Cutting Edge




Catalog #	Size, Width
OP-320-01	12 cm (4 3/4"), 1.0 mm

Diamond Knife
Tri-Facet Cutting Edge



Catalog #	Size, Width
OP-320-02	12 cm (4 3/4"), 1.0 mm

Diamond Knife
10-Facet Cutting Edge



Catalog #	Size, Width
OP-320-03	12 cm (4 3/4"), 1.0 mm

Crescent Diamond Knife



Catalog #	Size, Width
OP-331-25	11.5 cm (4 1/2"), 2.5 mm



Diamond Knives, Trephines
Diamantmesser, Trephinen
Bisturí de diamante, Trépanos
Couteaux de diamant, Trépanos
Cottelli di diamanti, Trepani

Diamond Phaco Knife
Straight, Keratome Shaped




Catalog #	Size, Width
OP-340-28	13 cm (5"), 2.8 mm
OP-340-30	13 cm (5"), 3.0 mm
OP-340-32	13 cm (5"), 3.2 mm

Diamond Phaco Knife
Angled, Keratome Shaped




Catalog #	Size, Width
OP-341-25	13 cm (5"), 2.8 mm
OP-341-30	13 cm (5"), 3.0 mm
OP-341-32	13 cm (5"), 3.2 mm

Diamond Phaco Knife
Keratome Shaped, 4 Cutting Edges



Catalog #	Size, Width
OP-351-25	13 cm (5"), 2.5 mm
OP-351-28	13 cm (5"), 2.8 mm
OP-351-30	13 cm (5"), 3.0 mm
OP-351-32	13 cm (5"), 3.2 mm

Diamond Phaco Knife
Trapezoid Cutting Edge



Catalog #	Size, Width
OP-355-27	12 cm (4 3/4"), 2.7 mm to 3.2 mm



Trephines
Trephinen
Trépanos
Trépan
Trepani

Elliot Scleral Trephine Blade

Catalog #	Size, Diameter
OP-370-10	3.5 cm (1 1/2"), Ø 1.0 mm
OP-370-15	3.5 cm (1 1/2"), Ø 1.5 mm
OP-370-20	3.5 cm (1 1/2"), Ø 2.0 mm



Elliot Handle

Catalog #	Size
OP-370-90	4.5 cm (1 3/4")



Camellin Lasek Trephine

Catalog #	Size, Diameter
OP-381-08	12 cm (4 3/4"), Ø 8 mm
OP-381-09	12 cm (4 3/4"), Ø 9 mm
OP-381-10	12 cm (4 3/4"), Ø 10 mm



Camellin Lasek Trephine & Alcohol Retaining Well

Catalog #	Size, Diameter
OP-385-08	12.5 cm (5"), Ø 8 mm, 8.5 mm
OP-385-09	12.5 cm (5"), Ø 9 mm, 9.5 mm
OP-385-10	12.5 cm (5"), Ø 10 mm, 10.5 mm





Trephines
Trephinen
Trépanos
Trépons
Trepani


Gass Corneoscleral Punch
Rotatable

Catalog #	Diameter, Deep Bite
OP-391-00	Ø 1.5 mm, 0.75 mm




Luntz-Dodick Trabeclectomy Punch Complete
With Puch Ref: OP-455-01

Catalog #	Diameter, Deep Bite
OP-395-00	Ø 1.0 mm, 0.5 mm




Luntz-Dodick Trabeclectomy Punch Only

Catalog #	Diameter, Deep Bite
OP-395-01	Ø 1.0 mm, 0.5 mm



Kelly Descemet's Membrane Punch

Catalog #	Diameter, Deep Bite
OP-397-00	Ø 1.0 mm, 0.5 mm



Christensen Punctum Punch
For Opening Stenotic Punctums

Catalog #	Size, Width
OP-399-00	10 cm (4"), 1.0 mm





Trephines, Cannulas
Trepinen, Kanülen
Trépanos, Cánulas
Trépan, Canules
Trepani, Cannule

Corneal Trephine Blade

Catalog #	Diameter
OP-400-01	Ø 6.0 mm
OP-400-02	Ø 6.5 mm
OP-400-03	Ø 6.75 mm
OP-400-04	Ø 7.0 mm
OP-400-05	Ø 7.25 mm
OP-400-06	Ø 7.5 mm
OP-400-07	Ø 7.75 mm
OP-400-08	Ø 8.0 mm
OP-400-09	Ø 8.25 mm
OP-400-10	Ø 8.5 mm
OP-400-11	Ø 8.75 mm
OP-400-12	Ø 9.0 mm
OP-400-13	Ø 9.5 mm
OP-400-14	Ø 14.0 mm



Corneal Trephine Handle

Catalog #	Diameter
OP-401-01	Ø 6.0 mm
OP-401-02	Ø 6.5 mm
OP-401-03	Ø 6.75 mm
OP-401-04	Ø 7.0 mm
OP-401-05	Ø 7.25 mm
OP-401-06	Ø 7.5 mm
OP-401-07	Ø 7.75 mm
OP-401-08	Ø 8.0 mm
OP-401-09	Ø 8.25 mm
OP-401-10	Ø 8.5 mm
OP-401-11	Ø 8.75 mm
OP-401-12	Ø 9.0 mm
OP-401-13	Ø 9.5 mm





Lacrimal Cannula
Straight



Catalog #	Gauge
OP-410-23	23 G

Lacrimal Cannula
Curved



Catalog #	Gauge
OP-411-23	23 G

McIntyre Lacrimal Cannula
Straight



Catalog #	Gauge
OP-412-23	23 G

McIntyre Lacrimal Cannula
Curved



Catalog #	Gauge
OP-413-23	23 G

West Lacrimal Cannula
Blunt Tip



Catalog #	Gauge
OP-415-23	23 G



Cannulas
Kanülen
Cánulas
Canules
Cannule

Air Injection Cannula
Angled at 5 mm

Catalog #	Gauge
OP-421-27	27 G
OP-421-30	30 G



Air Injection Cannula
Angled at 7 mm

Catalog #	Gauge
OP-423-27	27 G
OP-423-30	30 G



Randolph Cyclodialysis Cannula
Flat Tip

Catalog #	Gauge
OP-427-19	19 G
OP-427-21	21 G



Viscoelastic Aspirating Cannula
Angled at 10 mm

Catalog #	Gauge
OP-429-19	19 G
OP-429-22	22 G



Troutman Olive Tip Cannula
Delicate

Catalog #	Gauge
OP-433-25	25 G





Rainin Air Injection Cannula
Flat Tip, Angled at 7 mm

Catalog #	Gauge
OP-437-27	27 G
OP-437-30	30 G



Lasek Alcohol Cannula
Delicate Olive Shaped Tip

Catalog #	Gauge
OP-439-25	25 G



Stevens Sub Tenon's Anesthesia Cannula
Gently Curved, Flat Tip

Catalog #	Gauge
OP-441-19	19 G



Arkinson Retrobulbar Needle

Catalog #	Gauge
OP-443-23	23 G



Corbin Sub Tenon's Anesthesia Cannula
Smooth Rounded Tip

Catalog #	Gauge
OP-445-23	23 G





Cannulas
Kanülen
Cánulas
Canules
Cannule


Fukasaku Tenon's Anesthesia Cannula
Gently Curved, Smooth Tip

Catalog #	Gauge
OP-447-24	24 G



Simcoe (Reversed) Irrigation Aspirating Cannula

Catalog #	Gauge
OP-451-23	23 G



Simcoe (Reversed) Irrigation Aspirating Cannula

Catalog #	Gauge
OP-453-23	23 G



Simcoe Irrigation Aspirating Cannula
Right

Catalog #	Gauge
OP-455-23	23 G



Simcoe Irrigation Aspirating Cannula
Left


Catalog #	Gauge
OP-457-23	23 G





Gills Irrigation Aspirating Cannula

Catalog #	Gauge
OP-459-23	23 G




Kelman Cystotome

Catalog #	Gauge
OP-462-22	23 G




Kartz Cystotome

Catalog #	Gauge
OP-463-22	22 G



McIntyre Cystotome

Catalog #	Gauge
OP-465-23	23 G



Bishop-Harmon Anterior Chamber Irrigator

Catalog #	Gauge
OP-471-19	19 G





Cannulas
Kanülen
Cánulas
Canules
Cannule


Hydrodissection Cannula
Angled at 8 mm, Flat Tip

Catalog #	Gauge
OP-475-25	25 G



Seeley Hydrodissection Cannula
Angled at 12 mm, Flat Tip

Catalog #	Gauge
OP-477-26	25 G



Knolle Anterior Chamber Irrigation Cannula
Angled at 2 mm

Catalog #	Gauge
OP-481-23	23 G
OP-481-27	27 G
OP-481-30	30 G



Sauter Hydrodissection Cannula
Verrically Flattened Tip

Catalog #	Gauge
OP-485-27	27 G



Bracken Irrigation Cannula
Smooth Beveled Tip

Catalog #	Gauge
OP-487-19	19 G





Irrigation Cannula
Thin Wall

Catalog #	Gauge
OP-489-23	23 G




Alio Lasik Irrigation Cannula
With One Front and Two Side Ports

Catalog #	Gauge
OP-491-23	23 G



Irrigation Cannula

Catalog #	Gauge
OP-493-21	21 G



Spatulated Lasik Cannula
Flattened Spatulated Tip, End Opening

Catalog #	Gauge
OP-495-26	26 G



Fukasaku Lasik Cannula
Low Profile Tip, End Opening

Catalog #	Gauge
OP-501-23	23 G





Cannulas
Kanülen
Cánulas
Canules
Cannule

Carter Lasik Cannula
Flattened Tip, End Opening

Catalog #	Gauge
OP-503-25	25 G




Banaji Lasik Cannula
With 6 Ports

Catalog #	Gauge
OP-507-25	25 G



Buratto Lasik Cannula
With 3 Ports , Slightly Flattened

Catalog #	Gauge
OP-509-23	23 G
OP-509-25	25 G



Brierley Lasik Cannula
With 6 Ports , Spatulated Flattened

Catalog #	Gauge
OP-511-25	25 G



McIntyre Anterior Chamber Cannula
Straight


Catalog #	Gauge
OP-516-23	23 G
OP-516-26	26 G
OP-516-30	30 G






McIntyre Anterior Chamber Cannula
Angled

Catalog #	Gauge
OP-517-23	23 G
OP-517-26	26 G
OP-517-30	30 G



Gimbel Cortex Aspirating Cannula
End Opening

Catalog #	Gauge
OP-519-30	30 G




Simcoe Cortex Extractor
Smooth Tip, 13 mm Long

Catalog #	Gauge
OP-521-23	23 G




Simcoe Cortex Extractor
Smooth Tip, 21 mm Long

Catalog #	Gauge
OP-523-23	23 G



Simcoe Cortex Extractor
Smooth Closed Tip, 21 mm Long

Catalog #	Gauge
OP-525-21	21 G





Cannulas
Kanülen
Cánulas
Canules
Cannule

Gills Aspirating Cannula

Catalog #	Gauge, Tip Size
OP-529-05	25 G, 5 mm
OP-529-07	25 G, 7 mm
OP-529-08	25 G, 8 mm
OP-529-10	25 G, 10 mm



Fukasaku HydroChop Cannula

Vertically Flattened Tip, Dual Side Ports

Catalog #	Gauge
OP-531-23	23 G




Drews Irrigation Vectis

Catalog #	Gauge
OP-533-23	23 G



Sheets Irrigation Vectis

Catalog #	Gauge
OP-535-20	20 G



Knolle-Pearce Irrigation Vectis

Catalog #	Gauge
OP-537-23	23 G






Cannulas, Chalazion Curette
Kanülen, Chalazionslöffel
Cánulas, Cucharillas para Chalazión
Canules, Cuillers à chalazion
Cannule, Cucchiai per calazio


Simcoe Irrigation Lens Nucleus Loop
Right

Catalog #	Gauge
OP-539-25	25 G



Simcoe Irrigation Lens Nucleus Loop
Left

Catalog #	Gauge
OP-541-25	25 G



Drof Infusion Cannula
3 mm Long Tip

Catalog #	Gauge 1, Gauge 2
OP-543-22	18 G, 22 G



Heslin Infusion Cannula
Textured Tip

Catalog #	Gauge
OP-545-23	23 G
OP-545-25	25 G



Infusion Cannula
Self Retaining Tip

Catalog #	Gauge
OP-547-22	22 G





Ophthalmology Instruments
Ophthalmologie Instrumente
Ofthalmologia Instrumentos
Ophthalmologie Instruments
Ofthalmologia Strumenti

Meyerhoefer Chalazion Curette



Catalog #	Size, Diameter, Profile
OP-550-00	12.5 cm (5"), Ø 1.75 mm, Size 0
OP-550-01	12.5 cm (5"), Ø 1.2 mm, Size 1
OP-550-02	12.5 cm (5"), Ø 2.25 mm, Size 2
OP-550-03	12.5 cm (5"), Ø 2.5 mm, Size 3
OP-550-04	12.5 cm (5"), Ø 3.5 mm, Size 4



Nightingale Posterior Capsule
Polishing Curette, With Semi Sharp Edge



Catalog #	Size, Diameter
OP-571-12	12.5 cm (5"), Ø 1.75 mm

Shepard Capsule Polishing Curette
Open Ring Semi Sharp Edges on Top and Bottom




Catalog #	Size, Diameter
OP-573-12	12.5 cm (5"), Ø 1.75 mm

Blumenthal Conjunctiva Dissector
With Blunt Disc Shaped Tip



Catalog #	Size
OP-587-13	13 cm (5")





Burattor Lasik Flap Protector



Catalog #	Size
OP-591-12	12.5 cm (5")

Burattor Lasik Flap Protector
Blades Angled 15°



Catalog #	Size
OP-594-11	11.5 cm (4 1/2")

Mannis-Buratto Lasik Flap Protector
Blades Angled 45°



Catalog #	Size
OP-597-11	11.5 cm (4 1/2")

Rowen Combination Lasik Flap
Protector and Spatula



Catalog #	Size
OP-599-13	13 cm (5")

Sloane Lasik Micro Hoe
Semi-Circular Tip With Beveled Edges



Catalog #	Size
OP-611-12	12 cm (4 3/4")



Ophthalmology Instruments
Ophthalmologie Instrumente
Oftalmologia Instrumentos
Ophthalmologie Instruments
Oftalmologia Strumenti

Sloane Epi Peeler
Semi-Sharp Edges



Catalog #	Size
OP-615-11	11.5 cm (4 1/2")

Sloane Lasik Micro Hoe & Epi Peeler



Catalog #	Size
OP-619-13	13.5 cm (5 1/4")

Sloane Lasik Epi-Dissector
Spatulated Tip With Semi-Sharp Edges



Catalog #	Size
OP-623-11	11.5 cm (4 1/2")

Sloane Lasik Flap Reposer
Gently Curved Spatulated Tip



Catalog #	Size
OP-627-12	12 cm (4 3/4")

Culler Iris Spatula



Catalog #	Size, Width
OP-630-01	13 cm (5"), 1 mm
OP-630-02	13 cm (5"), 2 mm



Wheeler Double-End Spatula
Lightly Curved



Catalog #	Size
OP-632-13	13 cm (5")

Barraquer Iris Spatula
Round, Smooth



Catalog #	Size
OP-635-01	12 cm (4 3/4"), 0.25 mm
OP-635-02	12 cm (4 3/4"), 0.5 mm

Gimbel Nucleus Spatula
Round, Smooth, Angled 30°



Catalog #	Size, Tip Length, Tip Diameter
OP-637-12	12.5 cm (5"), 10 mm, Ø 0.25 mm

Escaf Nucleus Manipulator
Angled 45°



Catalog #	Size, Tip Length
OP-639-12	12 cm (4 3/4"), 13 mm

Knolle Nucleus Spatula
Malleable Flat, Stainless Steel



Catalog #	Size, Width
OP-640-12	12.5 cm (5"), 0.5 mm



Ophthalmology Instruments
Ophthalmologie Instrumente
Ofthalmologia Instrumentos
Ophthalmologie Instruments
Ofthalmologia Strumenti

**Koch Nucleus Spatula
Smooth**



Catalog #	Size
OP-643-01	11.5 cm (4 1/2")

**Koch Nucleus Spatula
With Notches**



Catalog #	Size
OP-643-02	11.5 cm (4 1/2")

Akahoshi Nucleus Sustainer



Catalog #	Size
OP-645-11	11.5 cm (4 1/2")

**Escaf Nucleus Sustainer
Smooth Ball Tip, Texture interior Surface**



Catalog #	Size
OP-647-12	12 cm (4 3/4")

Lieberman Nucleus Manipulator



Catalog #	Size, Profile
OP-649-01	12 cm (4 3/4"), Right Hand
OP-649-02	12 cm (4 3/4"), Left Hand



Lieberman Nucleus Manipulator and Spatula



Catalog #	Size, Profile
OP-651-01	14 cm (5 1/2"), Right Hand
OP-651-02	14 cm (5 1/2"), Left Hand

Fukasaku-Lieberman Nucleus Manipulator and Hook



Catalog #	Size, Profile
OP-653-01	13.5 cm (5 1/4"), Right Hand
OP-653-02	13.5 cm (5 1/4"), Left Hand

Masket Phaco Spatula
Round Spatula With Disc Shaped Tip



Catalog #	Size, Spatula Length
OP-657-13	13 cm (5"), 12 mm

Chang Phaco Chopper
Right Hand



Catalog #	Size, Tip Length, Profile
OP-659-01	13 cm (5"), 0.9 mm, Right Hand
OP-659-02	13 cm (5"), 0.9 mm, Left Hand

Minardi Phaco Chopper



Catalog #	Size, Tip Length
OP-661-11	11.5 cm (4 1/2"), 0.7 mm



Ophthalmology Instruments
Ophthalmologie Instrumente
Oftalmologia Instrumentos
Ophthalmologie Instruments
Oftalmologia Strumenti

Hwang Phaco Chopper



Catalog #	Size
OP-663-13	13.5 cm (5 1/4")

Agarwal Karate Phaco Chopper



Catalog #	Size
OP-665-11	11.5 cm (4 1/2")

Chang Micro Finger Chopper



Catalog #	Size
OP-667-11	11.5 cm (4 1/2")

Brierley Nucleus Splitter



Catalog #	Size
OP-669-12	12 cm (4 3/4")

Dodick-Kammann Nucleus Chopper
Straight



Catalog #	Size, Tip Length
OP-670-15	12 cm (4 3/4"), 1.5 mm
OP-670-20	12 cm (4 3/4"), 2.0 mm



Dodick-Kammann Nucleus Chopper
Angled

Catalog #	Size, Tip Length
OP-671-15	12 cm (4 3/4"), 1.5 mm
OP-671-20	12 cm (4 3/4"), 2.0 mm



Ecsaf Nucleus Divdider
Wedge Shaped

Catalog #	Size, Tip Length
OP-673-12	12 cm (4 3/4"), 1.25 mm



Neuhann Nucleus Divdider
Strongly Curved

Catalog #	Size,
OP-675-12	12 cm (4 3/4")



Nagahara Nucleus Spatula
With Flat Interior Edge

Catalog #	Size,
OP-677-11	11 cm (4 1/4")



Fakasaku Snapper Hook
Wedge Shaped

Catalog #	Size,
OP-679-11	11.5 cm (4 1/2")





Ophthalmology Instruments
Ophthalmologie Instrumente
Oftalmologia Instrumentos
Ophthalmologie Instruments
Oftalmologia Strumenti

Fakasaku Snapper Hook & Small Pupil
Wedge Shaped



Catalog #	Size
OP-681-11	11.5 cm (4 1/2")

Rosen Phaco Splitter
Wedge Shaped



Catalog #	Size
OP-683-11	11.5 cm (4 1/2")

Rosen Phaco Splitter
Angled 60°, Wedge Shaped



Catalog #	Size
OP-685-11	11.5 cm (4 1/2")

Rowen Phaco Spatula



Catalog #	Size
OP-687-13	13 cm (5")

Bechert Nucleus Rotator
Blunt Forked Tip, Angled at 7 mm



Catalog #	Size
OP-689-01	11.5 cm (4 1/2")



Jeff-Bechert Nucleus Rotator
Blunt Forked Tip, Angled at 9 mm



Catalog #	Size
OP-689-02	11.5 cm (4 1/2")

Clayman Nucleus Rotator
With V-Notched Tip



Catalog #	Size
OP-691-11	12 cm (4 3/4")

Drysdale Nucleus Manipulator
Paddle Shaped Tip



Catalog #	Size
OP-693-11	11 cm (4 1/4")

Pisacano Nucleus Rotator
Horizontal Flat Disc Shaped Tip



Catalog #	Size
OP-695-11	11.5 cm (4 1/2")

Pisacano Nucleus Rotator
Vertical Flat Disc Shaped Tip



Catalog #	Size
OP-697-11	11.5 cm (4 1/2")



Ophthalmology Instruments
Ophthalmologie Instrumente
Oftalmologia Instrumentos
Ophthalmologie Instruments
Oftalmologia Strumenti

Siepsner Endocapsular Controller
Trumpet Shaped Tip



Catalog #	Size
OP-699-11	11.5 cm (4 1/2")

Anis Nucleus Manipulator
Ball Tip, Straight



Catalog #	Size, Tip Diameter
OP-700-12	12 cm (4 3/4"), Ø 0.5 mm

Anis Nucleus Manipulator
Ball Tip, Angled



Catalog #	Size, Tip Diameter
OP-701-11	11.5 cm (4 1/2"), Ø 0.5 mm

Tennant Nucleus Rotator
Ball Shaped Tip



Catalog #	Size, Tip Diameter
OP-703-12	12 cm (4 3/4"), Ø 0.6 mm

Lippman Nucleus Spatula
Shaft With Blunt Forked Tip



Catalog #	Size, Tip Diameter
OP-705-11	11.5 cm (4 1/2"), Ø 0.5 mm



Girard Synechiaie Spatula
Semi-Sharp Tip



Catalog #	Size, Tip Diameter
OP-707-11	11.5 cm (4 1/2"), Ø 0.5 mm

Obstbaum Synechiaie Spatula
Semi-Sharp Tip



Catalog #	Size, Tip Diameter
OP-709-12	12 cm (4 3/4"), Ø 0.25 mm

Moran Lasik Enhancement Spatula



Catalog #	Size
OP-711-13	13.5 cm (5 1/4")

Hersh Flap Spatula



Catalog #	Size
OP-713-14	14 cm (5 1/4")

Binder Flap Elevator



Catalog #	Size
OP-715-13	13 cm (5")



Ophthalmology Instruments
Ophthalmologie Instrumente
Ofthalmologia Instrumentos
Ophthalmologie Instruments
Ofthalmologia Strumenti

Fukasaku Lasik Spatula



Catalog #	Size
OP-717-12	12 cm (4 3/4")

Guimaraes Lasik Spatula
Gently Curved



Catalog #	Size
OP-719-12	12 cm (4 3/4")

Katzen Lasik Flap Unzipper



Catalog #	Size
OP-720-13	13.5 cm (5 1/4")

Maddox Lasik Spatula



Catalog #	Size
OP-723-12	12.5 cm (5")

Lasik Flap Manipulator
Delicate Olive Tip



Catalog #	Size
OP-725-11	11.5 cm (4 1/2")



Jaff Lens Spatula
Flat



Catalog #	Size, Width
OP-727-12	12 cm (4 3/4"), 0.5 mm

Lester IOL Manipulator
Straight



Catalog #	Size, Tip Diameter
OP-728-12	12 cm (4 3/4"), Ø 0.25 mm

Lester IOL Manipulator
Angled



Catalog #	Size, Tip Diameter
OP-729-12	12 cm (4 3/4"), Ø 0.25 mm

Lester IOL Manipulator
Vaulted Shank



Catalog #	Size, Tip Diameter
OP-731-11	11.5 cm (4 1/2"), Ø 0.25 mm

Hunkeler IOL Manipulator
Hemi-Ball Tip



Catalog #	Size, Tip Diameter
OP-733-12	12 cm (4 3/4"), Ø 0.3 mm



Ophthalmology Instruments
Ophthalmologie Instrumente
Oftalmologia Instrumentos
Ophthalmologie Instruments
Oftalmologia Strumenti

Alio-Kelman IOL Hook



Catalog #	Size, Tip Diameter
OP-735-12	12 cm (4 3/4"), Ø 0.2 mm

Aker Lens Pusher

Smooth, Fish Mouth Shaped Tip, With Guard



Catalog #	Size
OP-738-12	12.5 cm (5")

Bechert Lens Pusher

Vertical, Y Shaped Tip, Angled



Catalog #	Size
OP-739-11	11.5 cm (4 1/2")

Elsching Cyclodialysis Spatula

Flat



Catalog #	Size
OP-741-13	13 cm (5")

Castroviejo Synechiaie Spatula

Semi Sharp Tips



Catalog #	Size, Width
OP-743-13	13 cm (5"), 0.5 mm



Castroviejo Cyclodialysis Spatula

Catalog #	Size, Width
OP-745-25	13 cm (5"), 0.25 mm
OP-745-50	13 cm (5"), 0.5 mm



Lewis Lens Loop

Catalog #	Size, Profile
OP-750-01	13 cm (5"), Small
OP-750-02	14 cm (5 1/2"), Large



New Orleans Lens Loop

Catalog #	Size
OP-752-13	13 cm (5")



Wilder Lens Loop
Serrated

Catalog #	Size
OP-754-13	13.5 cm (5 1/4")



Kirby Lens Loop and Muscle Hook

Catalog #	Size
OP-757-15	15 cm (6")





Ophthalmology Instruments
Ophthalmologie Instrumente
Ofthalmologia Instrumentos
Ophthalmologie Instruments
Ofthalmologia Strumenti

Morrison Lens Loop and Probe



Catalog #	Size
OP-758-14	14.5 cm (5 3/4")

Kansas Lens Loop
Serrated



Catalog #	Size, Width
OP-760-13	13.5 cm (5 1/4"), 3 mm

Keener-Arlt Lens Loop
Serrated



Catalog #	Size, Width
OP-761-13	13.5 cm (5 1/4"), 3 mm

Alfonso Nucleus Trisector
For Use With Ref: OP-775-02



Catalog #	Size, Width
OP-765-10	10.5 cm (4 1/4"), 1.5 mm

Bucher-Kansas Nucleus Trisector
For Use With Ref: OP-775-03



Catalog #	Size, Width
OP-767-10	10.5 cm (4 1/4"), 2.5 mm



Kansas Nucleus Trisector
For Use With Ref: OP-773-12



Catalog #	Size, Width
OP-769-11	11 cm (4 1/4"), 4 mm

Kansas Nucleus Bisector



Catalog #	Size
OP-771-11	11.5 cm (4 1/2")

Kansas Nucleus Vectis



Catalog #	Size
OP-773-12	12 cm (4 3/4")

Alfonso Cutting Platform



Catalog #	Size, Width
OP-775-02	12.5 cm (5"), 2 mm
OP-775-03	12.5 cm (5"), 3 mm

Giannetti Nucleus Capture Spoon



Catalog #	Size
OP-779-12	12.5 cm (5")



Ophthalmology Instruments
Ophthalmologie Instrumente
Oftalmologia Instrumentos
Ophthalmologie Instruments
Oftalmologia Strumenti

Giannetti Nucleus Capture Spatula



Catalog #	Size
OP-781-12	12.5 cm (5")

Giannetti Nucleus Capture Hook



Catalog #	Size
OP-783-12	12.5 cm (5")

Culler Lens Spoon
Sterling Silver



Catalog #	Size
OP-790-13	13 cm (5")

Paton Spatula & Spoon
For Corneal Button Transfer



Catalog #	Size
OP-793-14	14 cm (5 1/2")

Rosenwasser Lamella Donor Shovel



Catalog #	Size, Width, Profile
OP-795-04	11.5 cm (4 1/2"), 4 mm, With Retaining Rim on 2 Sides
OP-795-05	11.5 cm (4 1/2"), 5 mm, With Retaining Rim on 3 Sides





Bunge Evisceration Spoon



Catalog #	Size
OP-800-01	13.5 cm (5 1/4"), Small
OP-800-02	14 cm (5 1/2"), Large

Wells Enucleation Spoon



Catalog #	Size
OP-802-15	15 cm (6")

Graether Collar Button
Straight



Catalog #	Size
OP-804-12	12.5 cm (5")



Ophthalmology Instruments
Ophthalmologie Instrumente
Ofthalmologia Instrumentos
Ophthalmologie Instruments
Ofthalmologia Strumenti

Graether Collar Button
Angled



Catalog #	Size
OP-805-12	12.5 cm (5")

Z Hook Manipulator
Z Shaped Blunt Tip



Catalog #	Size, Tip Diameter
OP-811-11	11.5 cm (4 1/2"), Ø 0.15 mm

Reversed Sinskey Hook
Blunt Tip



Catalog #	Size, Tip Diameter
OP-813-11	11.5 cm (4 1/2"), Ø 0.2 mm

Bonn Micro Iris Hook
Blunt Tip



Catalog #	Size
OP-815-12	12.5 cm (5")

Shephard Iris Hook
Blunt Tip



Catalog #	Size
OP-817-12	12.5 cm (5")



Jaff Lens Manipulating Hook
Blunt Tip, Straight



Catalog #	Size, Tip Diameter
OP-819-01	11.5 cm (4 1/2"), Ø 0.15 mm

Clayman-Jaffe Lens Manipulating Hook
Blunt Tip, Angled



Catalog #	Size, Tip Diameter
OP-819-02	11.5 cm (4 1/2"), Ø 0.15 mm

Sinsky Lens Manipulating Hook
Blunt Tip, Straight



Catalog #	Size, Tip Diameter
OP-821-01	11.5 cm (4 1/2"), Ø 0.25 mm

Harris-Sinsky Lens Manipulating Hook
Blunt Tip, Angled



Catalog #	Size, Tip Diameter
OP-821-02	11.5 cm (4 1/2"), Ø 0.25 mm

Lewicky Lens Manipulating Hook
Blunt Tip, Angled



Catalog #	Size, Tip Diameter
OP-823-12	12 cm (4 3/4"), Ø 0.5 mm



Ophthalmology Instruments
Ophthalmologie Instrumente
Oftalmologia Instrumentos
Ophthalmologie Instruments
Oftalmologia Strumenti

Hirschman Iris Hook
Blunt, Angled



Catalog #	Size
OP-825-11	11 cm (4 1/4")

Sinskey Lens Manipulating Hook



Catalog #	Size, Tip Diameter, Profile
OP-826-01	11.5 cm (4 1/2"), Ø 0.25 mm, Straight
OP-826-02	11.5 cm (4 1/2"), Ø 0.25 mm, Angled

Sinskey II Lens Manipulating Hook
Straight



Catalog #	Size, Tip Diameter, Profile
OP-826-03	11.5 cm (4 1/2"), Ø 0.2 mm, Straight
OP-826-04	11.5 cm (4 1/2"), Ø 0.2 mm, Angled

Kuglen Iris Hook and Lens Manipulator
Push Pull Model



Catalog #	Size, Profile
OP-827-01	11.5 cm (4 1/2"), Straight
OP-827-02	11.5 cm (4 1/2"), Angled

Kuglen Iris Hook and Lens Manipulator
Clover Leaf Model



Catalog #	Size
OP-829-01	11.5 cm (4 1/2")



Rentsch Boat Hook
Straight With Guard



Catalog #	Size
OP-831-12	12 cm (4 3/4")

Osher Y Hook
Straight With Guard



Catalog #	Size
OP-833-12	12 cm (4 3/4")

Fenzl Lens Manipulating Hook
Straight With Guard



Catalog #	Size, Tip Diameter
OP-835-01	12 cm (4 3/4"), Ø 0.15 mm

Fenzl Lens Manipulating Hook
Angled



Catalog #	Size, Tip Diameter
OP-835-02	12 cm (4 3/4"), Ø 0.15 mm

Maltzman-Fenzl Lens Manipulating Hook
Straight With Guard



Catalog #	Size, Tip Diameter
OP-837-01	12 cm (4 3/4"), Ø 0.15 mm



Ophthalmology Instruments
Ophthalmologie Instrumente
Ofthalmologia Instrumentos
Ophthalmologie Instruments
Ofthalmologia Strumenti

Maltzman-Fenzl Lens Manipulating Hook
Angled



Catalog #	Size, Tip Diameter
OP-837-02	12 cm (4 3/4"), Ø 0.15 mm

Jaffe-Maltzman Lens Manipulating Hook
Angled



Catalog #	Size, Tip Diameter
OP-839-02	12 cm (4 3/4"), Ø 0.15 mm

Giannetti Push Pull Manipulator



Catalog #	Size
OP-841-01	12.5 cm (5")

Deitz ICL Manipulator
Blunt Forked Tip



Catalog #	Size
OP-843-02	12 cm (4 3/4")

Fixation Hook
Sharp



Catalog #	Size, Profile
OP-845-01	12.5 cm (5"), Small
OP-845-02	13 cm (5"), Large




Twist Hook





Catalog #	Size, Profile
OP-849-01	11 cm (4 1/4"), Right
OP-849-02	11 cm (4 1/4"), Left

Silverman Fixation Hook



Catalog #	Size
OP-851-09	9.5 cm (3 3/4")

Lu-Mendez LRI Guide and Fixation Ring




Catalog #	Size, I.D
OP-861-11	11 cm (4 1/4"), 11 mm

Gimbel Stabilization Ring
With Swivel Handle



Catalog #	Size, Tip Diameter
OP-863-11	11.5 cm (4 1/2"), 13 mm

Fine-Thronton Swivel Fixation Ring
With Swivel Handle



Catalog #	Size, Tip Diameter
OP-865-13	12 cm (4 3/4"), 13 mm
OP-865-16	12 cm (4 3/4"), 16 mm



Ophthalmology Instruments
Ophthalmologie Instrumente
Ofthalmologia Instrumentos
Ophthalmologie Instruments
Ofthalmologia Strumenti

Thronton Globe Fixation Ring



Catalog #	Size, Tip Diameter
OP-867-10	10 cm (4"), 16 mm

Hofmann-Thronton Globe Fixation Ring
With Swivel Handle



Catalog #	Size, Tip Diameter
OP-869-11	11 cm (4 1/4"), 16 mm

Stevens Tenotomy Hook
Blunt, Small



Catalog #	Size
OP-871-12	12.5 cm (5")

Helveston Teaser Hook
Extra Delicate, Blunt



Catalog #	Size
OP-873-12	12.5 cm (5")

Gass Retinal Detachment Hook
With Oval Hole



Catalog #	Size
OP-875-13	13 cm (5")



Helveston Finder Hook



Catalog #	Size, Profile
OP-877-01	12.5 cm (5")?, Large
OP-877-02	12 cm (4 3/4"), Small

Helveston Muscle Hook



Catalog #	Size, Tip Length, Profile
OP-879-01	13 cm (5"), 8 mm, Fig. 1
OP-879-02	13 cm (5"), 10 mm, Fig. 2
OP-879-03	13 cm (5"), 12 mm, Fig. 3

Jameson Muscle Hook



Catalog #	Size, Profile
OP-881-01	13 cm (5"), Small
OP-881-02	13.5 cm (5 1/4"), Large

Green Muscle Hook



Catalog #	Size
OP-883-13	13.5 cm (5 1/4")



Ophthalmology Instruments
Ophthalmologie Instrumente
Oftalmologia Instrumentos
Ophthalmologie Instruments
Oftalmologia Strumenti

Kellen Capsulorrhexis Marker

6 mm Diameter Marking Ring With 2.5 mm Extensions



Catalog #	Size
OP-891-11	11.5 cm (4 1/2")

Maloney Intra-Operative Keratometer

For Qualitative Measurement of Astigmatism



Catalog #	Size
OP-900-00	11.5 cm (4 1/2")

Mendez Degree Gauge



Catalog #	Size
OP-901-11	11.5 cm (4 1/2")

Greene Visual Axis Marker

Concave Tip



Catalog #	Size, Tip Diameter
OP-905-11	11.5 cm (4 1/2"), 0.5 mm



McDonald Optic Zone Marker
Double End With Cross Hairs



Catalog #	Size, Tip Diameter 1, Tip Diameter 2
OP-907-11	11 cm (4 1/4"), 6.5 mm, 7.0 mm

Bores Optic Zone Marker

Catalog #	Size, Tip Diameter
OP-911-01	10.5 cm (4 1/4"), 3.0 mm
OP-911-02	10.5 cm (4 1/4"), 3.25 mm
OP-911-03	10.5 cm (4 1/4"), 3.5 mm
OP-911-04	10.5 cm (4 1/4"), 3.75 mm
OP-911-05	10.5 cm (4 1/4"), 4.0 mm
OP-911-06	10.5 cm (4 1/4"), 4.25 mm
OP-911-07	10.5 cm (4 1/4"), 4.5 mm
OP-911-08	10.5 cm (4 1/4"), 4.75 mm
OP-911-09	10.5 cm (4 1/4"), 5.0 mm
OP-911-10	10.5 cm (4 1/4"), 6.0 mm
OP-911-11	10.5 cm (4 1/4"), 7.0 mm
OP-911-12	10.5 cm (4 1/4"), 8.0 mm
OP-911-13	10.5 cm (4 1/4"), 9.0 mm





Ophthalmology Instruments
Ophthalmologie Instrumente
Ofthalmologia Instrumentos
Ophthalmologie Instruments
Ofthalmologia Strumenti

Hoffer Optic Zone Marker
With Cross Hairs

Catalog #	Size, Tip Diameter
OP-913-01	10.5 cm (4 1/4"), 3.0 mm
OP-913-02	10.5 cm (4 1/4"), 3.25 mm
OP-913-03	10.5 cm (4 1/4"), 3.5 mm
OP-913-04	10.5 cm (4 1/4"), 3.75 mm
OP-913-05	10.5 cm (4 1/4"), 4.0 mm
OP-913-06	10.5 cm (4 1/4"), 4.25 mm
OP-913-07	10.5 cm (4 1/4"), 4.5 mm
OP-913-08	10.5 cm (4 1/4"), 4.75 mm
OP-913-09	10.5 cm (4 1/4"), 5.0 mm
OP-913-10	10.5 cm (4 1/4"), 5.5 mm
OP-913-11	10.5 cm (4 1/4"), 6.0 mm
OP-913-12	10.5 cm (4 1/4"), 7.0 mm
OP-913-13	10.5 cm (4 1/4"), 8.0 mm
OP-913-14	10.5 cm (4 1/4"), 9.0 mm



Anis Radial Marker

Catalog #	Size, Center Opening, O.D, Profile
OP-915-08	10 cm (4"), 4.0 mm, Ø 12.5 mm, 8 Blades
OP-915-12	10 cm (4"), 4.0 mm, Ø 12.5 mm, 12 Blades
OP-915-16	10 cm (4"), 4.0 mm, Ø 12.5 mm, 16 Blades



Anis Suture Placement Marker
8 Radial Blades

Catalog #	Size, Center Opening
OP-917-10	10 cm (4"), 8 mm





Fechtner Trabeclectomy Marker



Catalog #	Size
OP-919-09	9 cm (3 1/2")


Helveston Scleral Ruler

15 mm Long in Notches in 5 mm Increments




Catalog #	Size
OP-930-12	12 cm (4 3/4")

Stahl Corneal Caliper




OP-932-01



OP-932-02

Catalog #	Size
OP-932-01	Standard
OP-932-02	Standard

Harms Trabeculotomy Probe



Catalog #	Size, Profile
OP-941-01	5 cm (2"), Right
OP-941-02	5 cm (2"), Left

Infant Lacrimal Dilator

Very Fine Blunt Tip



Catalog #	Size
OP-944-08	8 cm (3")



Ophthalmology Instruments
Ophthalmologie Instrumente
Ofthalmologia Instrumentos
Ophthalmologie Instruments
Ofthalmologia Strumenti

Muldoon Lacrimal Dilator



Catalog #	Size
OP-946-10	10 cm (4")

Wilder Lacrimal Dilator



Catalog #	Size, Profile
OP-948-01	10.5 cm (4 1/4"), Fine Taper, Size 1
OP-948-02	10.5 cm (4 1/4"), Medium Taper, Size 2
OP-948-03	10.5 cm (4 1/4"), Heavy Taper, Size 3

Castroviejo Lacrimal Dilator

Double Ended Fine and Medium Taper



Catalog #	Size
OP-950-14	14 cm (5 1/2")



COPYRIGHT NOTICE

The copyright, trademarks, domain names and other intellectual property rights in all materials and information in this presentation belongs to Surtex Instruments Ltd. All rights are reserved. Technical alterations are possible. This presentation may be used for no other purpose than offering, buying and selling of Surtex Instruments products. No part may be reproduced or copied in any form. In the event of misuse we retain the right to recall our catalogues and take appropriate legal action.

DISCLAIMER

This documents provides general information about Surtex Instruments Ltd, its products and services, and summarises general capabilities and offerings, which we deem to be relevant to our customers. Whilst Surtex Instruments has taken appropriate steps to ensure the accuracy of the information contained in the document Surtex Instruments gives no warranty regarding the accuracy or completeness of such information.

The presentation was prepared to provide information and does not constitute a contract. Surtex reserves the right, without prior notice, to change, delete, supplement, or otherwise amend at any time the information, images and offerings contained in the presentation or other documents. Surtex shall not be liable for any loss or damage of whatever nature (whether direct, indirect, consequential or other) may arise as a result of any third party relying on the information contained in this document.

For further detailed information, prices, terms and conditions, customer should contact Surtex Instruments at sales@surtex-instruments.co.uk.

© 2020 SURTEX INSTRUMENTS. ALL RIGHTS RESERVED. PRINTED IN UNITED KINGDOM. BT1-01 Rev 1.0



Unit 202 Kingspark Business Centre
152-178 Kingston Road, New Malden
Surrey KT3 3ST United Kingdom
Phone: +44 (0) 208 123 3513

E-mail: sales@surtex-instruments.co.uk
Website: www.surtex-instruments.co.uk



A member of the ABHI



Association of
British HealthTech
Industries

Health Tech **For Life**